

HallmarkSurgicalTM
When you care.

**PROVING
INDISPENSABLE
TO THOSE WHO CARE**

CATALOGUE EDITION 1

Our Story

When you care.

The Healthcare Profession, as caring for life, is arguably the most responsible and most respected professions known to man. Life is precious and perilous. To fall short when a life, (or lives), are at stake is unthinkable - on both your part and ours.

Hallmark Surgical - the name stands for excellence - has deliberately set the bar high. The instruments and devices we provide are devoted to preserving that life. But, more than that, Hallmark Surgical upholds the skill and experience of people by providing you with trusted knowledge and support, at all levels of health care.

Because, like you, every day we both protect and serve.

What's in it for you?



NZ OWNED &
OPERATED



TRUSTED
KNOWLEDGE



INDISPENSABLE
SOLUTIONS



PROVEN PREMIUM
SERVICE



INNOVATIVE
DEVICES



GLOBAL
SOURCING



SUPPLY CHAIN
EXCELLENCE



DELIVERY
GUARANTEE

Proving indispensable to those who care.

HallmarkSurgical™
When you care.

BIOPSY DEVICES AND ACCESSORIES



5

BIOPSY DEVICES AND ACCESSORIES

ELECTROSURGICAL



11

ELECTROSURGICAL

SURGICAL ACCESSORIES



41

SURGICAL ACCESSORIES

SURGICAL SUCTION



49

SURGICAL SUCTION

HAND SURGERY



61

HAND SURGERY

COLORECTAL & ENDOSCOPY DEVICES



67

COLORECTAL & ENDOSCOPY DEVICES

HOLLOWARE



83

HOLLOWARE

SOLO™ SINGLE USE INSTRUMENTS



93

SOLO™ SINGLE USE INSTRUMENTS

WOMENS HEALTH DEVICES



103

WOMENS HEALTH DEVICES


Biopsy Devices & Accessories



SHARPS

Punch Biopsies	6-7
BiopSafe	8-9

SOLO™ Punch Biopsy Range



Punch Biopsies

Full range of colour coded sizes for easy reference.

SOLO™ PUNCH BIOPSY

- 2mm - 12mm diameter range
- Precision sharpened blade
- Colour coded by size
- Single use, latex free, sterile packed
- UOM: Each, 25/Dispenser Box

Order Code	Product
B802XS	2mm
B803S	3mm
B804SM	4mm
B805M	5mm
B806ML	6mm
B808XL	8mm
B8123XL	12mm



SOLO™ PUNCH BIOPSY WITH EJECTOR

- 2mm, 3mm and 4mm diameter range with push button ejector to easily remove biopsy
- Precision sharpened blade
- Colour coded by size
- Single use, latex free, sterile packed
- UOM: Each, 25/Dispenser Box

Order Code	Product
B822XS	2mm
B823S	3mm
B824SM	4mm



BiopSafe

Innovative Prefilled Histology Containers

The secret is in the lid

With BiopSafe, there are no difficult formalin bags and wall holders that often leads to formalin waste. The system consists of a small, easy to use container (vial) with formalin capsuled in the lid. When the biopsy is placed at the bottom of the container, you screw the lid on and gently apply pressure with your thumb, allowing the formalin to flow out and cover the biopsy. Everything is done inside the container ensuring the clinician is protected from both the liquid and vapour of formalin.

Hermetically Sealed System

As well as causing irritation to mucous membranes and allergy, formalin is a carcinogenic substance. Therefore it is very important to secure a safe working environment when handling formalin. With BiopSafe, the handling is done in a hermetically sealed system providing the best protection for staff. The staff are never exposed to formalin, not even during the process of screwing the lid on or off.

The new road to safer biopsy handling

One of the most common problems in connection with biopsy handling is the risk of being exposed to formalin, either through touch or inhalation, a risk doctors, laboratory technicians and nurses are exposed to every day. With BiopSafe, the problem is finally solved.


Ready in seconds

Besides from being safe, the method is also fast. BiopSafe is always at hand and the process of placing the biopsy, adding formalin and making it ready for transportation to the laboratory is done in a matter of seconds. BiopSafe saves valuable time!



Order Code	Product
B8320LE	20ml Container
B8360RY	60ml Container

Electrosurgical



ELECTROSURGICAL

Electrisure™ Premium Disposable Bipolar Forceps	12-14
ZAP™ Disposable Bipolar Forceps	15-19
Stingray Reusable Bipolar Forceps	20-25
Bipolar Cables	26
Disposable Laparoscopic Electrodes	27-30
Diathermy Pencils	31
Monopolar Electrodes	32-33
Patient Return Pad Electrodes	34
Monopolar Cables	35
Electrosurgical Accessories	36-37
ALSA Diathermy Generators & Accessories	38-40

Electrisure™

Premium Non-Stick Disposable Bipolar Forceps

Silver alloy tips for exceptional non-stick performance



Perfect tip quality for every procedure

Angled cable connector reduces cable drag and weight on the forcep

Large selection of forcep styles to suit any surgical procedure

Electrisure™ Premium Disposable Bipolar Forceps

Silver-plated non-stick tips

ELECTRISURE™ MCPHERSON BIPOLAR FORCEPS

- NSI™ Silver alloy non-stick, low smoke tips
- Durable baked nylon insulation
- 3m cable with molded generator plug
- Angled cable connection to Forcep for easier movement
- Sterile packed, latex free, single use
- UOM: Each or 10/Dispenser Box

Order Code	Product
E801105MF	8.9cm x 0.5mm tips
E801110FG	8.9cm x 1.0mm tips



ELECTRISURE™ ADSON BIPOLAR FORCEPS

- NSI™ Silver alloy non-stick, low smoke tips
- Durable baked nylon insulation
- 3m cable with molded generator plug
- Angled cable connection to Forcep for easier movement
- Sterile packed, latex free, single use
- UOM: Each or 10/Dispenser Box

Order Code	Product
E802105AS	12.7cm x 0.5mm tips
E802110DA	12.7cm x 1.0mm tips



Electrisure™ Premium Disposable Bipolar Forceps

Silver-plated non-stick tips

ELECTRISURE™ BAYONET BIPOLAR FORCEPS



- NSI™ Silver alloy non-stick, low smoke tips
- Durable baked nylon insulation
- 3m cable with molded generator plug
- Angled cable connection to Forcep for easier movement
- Sterile packed, latex free, single use
- UOM: Each or 10/Dispenser Box

Order Code	Product
E804105YB	15.2cm x 0.5mm tips
E804110AN	15.2cm x 1.0mm tips
E804115TE	15.2cm x 1.5mm tips
E804205EG	20.3cm x 0.5mm tips
E804210HT	20.3cm x 1.0mm tips
E804215JW	20.3cm x 1.5mm tips

ELECTRISURE™ STRAIGHT BIPOLAR FORCEPS



- NSI™ Silver alloy non-stick, low smoke tips
- Durable baked nylon insulation
- 3m cable with molded generator plug
- Angled cable connection to Forcep for easier movement
- Sterile packed, latex free, single use
- UOM: Each or 10/Dispenser Box


Order Code	Product
E803110RG	15.9cm x 1.0mm tips

Zap™

Your ideal everyday bipolar forcep

Finely machined
stainless steel tips

Ideal for short procedures and
high usage applications



Uncompromised tip quality
for every procedure
guaranteeing consistent
outcomes

Angled cable connector
reduces cable drag and
weight on the forcep

Large selection of forcep
styles to suit any surgical
procedure

HallmarkSurgical™
When you care.

ZAP™ Disposable Bipolar Forceps

Perfect tip condition for every procedure, every time

ZAP™ MCPHERSON BIPOLAR FORCEPS




- Stainless Steel tips
- Durable baked nylon insulation
- 3m cable with molded generator plug
- Angled cable connection to Forcep for easier movement
- Sterile packed, latex free, single use
- UOM: Each or 50/Dispenser Box



Order Code	Product
E8200MC00	10.2 cm x Needle tips
E8201MC05	10.2cm x 0.5mm tips
E8202MC10	10.2cm x 1.0mm tips
E8203MA00	10.2cm x Needle tips Angled
E8204MA05	10.2cm x 0.5mm tips Angled

ZAP™ JEWELLER BIPOLAR FORCEPS



- Stainless Steel tips
- Durable baked nylon insulation
- 3m cable with molded generator plug
- Angled cable connection to Forcep for easier movement
- Sterile packed, latex free, single use
- UOM: Each or 50/Dispenser Box

Order Code	Product
E8210J00	11.4cm x Needle tips
E8211J05	11.4cm x 0.5mm tips
E8212J10	11.4cm x 1.0mm tips
E8213JA00	11.4cm x Needle tips Angled
E8214JA05	11.4cm x 0.5mm tips Angled

ZAP™ Disposable Bipolar Forceps

Perfect tip condition for every procedure, every time

ZAP™ ADSON BIPOLAR FORCEPS

- Stainless Steel tips
- Durable baked nylon insulation
- 3m cable with molded generator plug
- Angled cable connection to Forcep for easier movement
- Sterile packed, latex free, single use
- UOM: Each or 50/Dispenser Box

Order Code	Product
E8215A05	12.1cm x 0.5mm tips
E8217A10	12.1cm x 1.0mm tips



ZAP™ STRAIGHT BIPOLAR FORCEPS

- Stainless Steel tips
- Durable baked nylon insulation
- 3m cable with molded generator plug
- Angled cable connection to Forcep for easier movement
- Sterile packed, latex free, single use
- UOM: Each or 50/Dispenser Box

Order Code	Product
E8219L05	15.2cm x 0.5mm tips
E8220L10	15.2 cm x 1.0mm tips
E8221L15	15.2cm x 1.5mm tips
E8222L20	15.2cm x 2.0mm tips
E8240L05	18.0cm x 0.5mm tips
E8241L10	18.0cm x 1.0mm tips
E8242L15	18.0cm x 1.5mm tips
E8243L20	18.0cm x 2.0mm tips



ZAP™ Disposable Bipolar Forceps

Perfect tip condition for every procedure, every time

ZAP™ BAYONET BIPOLAR FORCEPS

- Stainless Steel tips
- Durable baked nylon insulation
- 3m cable with molded generator plug
- Angled cable connection to Forcep for easier movement
- Sterile packed, latex free, single use
- UOM: Each or 50/Dispenser Box



Order Code	Product
E8225H05	15.2cm x 0.5mm tips
E8226H10	15.2cm x 1.0mm tips
E8227H15	15.2cm x 1.5mm tips
E8228H20	15.2cm x 2.0mm tips
E8230H05	20.3cm x 0.5mm tips
E8231H10	20.3cm x 1.0mm tips
E8232H15	20.3cm x 1.5mm tips
E8233H20	20.3cm x 2.0mm tips

ZAP™ MCPHERSON BIPOLAR FORCEPS - NO CABLE

- Stainless Steel tips
- Durable baked nylon insulation
- No cable - US 2 pin plug
- Sterile packed, latex free, single use
- UOM: Each or 50/Dispenser Box



Order Code	Product
E831105WN	8.9cm x 0.5mm tips

ZAP™ Disposable Bipolar Forceps

Perfect tip condition for every procedure, every time

ZAP™ ADSON BIPOLAR FORCEPS - NO CABLE

- Stainless Steel tips
- Durable baked nylon insulation
- No cable - US 2 pin plug
- Sterile packed, latex free, single use
- UOM: Each or 50/Dispenser Box

Order Code	Product
E832105BQ	12.7cm x 0.5mm tips
E832110ST	12.7cm x 1.0mm tips



ZAP™ DISPOSABLE BIPOLAR FORCEPS

ZAP™ STRAIGHT BIPOLAR FORCEPS - NO CABLE

- Stainless Steel tips
- Durable baked nylon insulation
- No cable - US 2 pin plug
- Sterile packed, latex free, single use
- UOM: Each or 50/Dispenser Box


Order Code	Product
E833110PM	17.8cm x 1.0mm tips



Stingray Bayonet Bipolar Forcep, Silver Alloy Tips

Minimal viewing obstruction with fin alignment for precision tip control and the elimination of scissoring action

BAYONET 20.3CM KEY FORCEP - REUSABLE




- US 2-Pin Style Connection
- Keyfin design available to assure the alignment of tines and prevent scissoring action
- Non-Stick Silver Alloy Tips
- Ergonomic Handles
- Available with angled tips at 20, 30 and 45 degrees
- Made in USA, Reusable

Order Code	Product
EUS03-230NS	0.4mm Tip, Non-Stick
EUS03-233NS	0.7mm Tip, Non-Stick
EUS03-235NS	1.0mm Tip, Non-Stick
EUS03-237NS	1.2mm Tip, Non-Stick
EUS03-238NS	1.5mm Tip, Non-Stick

Order Code	Product
EUS03-230FG	0.4mm Tip, Non-Stick with Fin Guide
EUS03-233FG	0.7mm Tip, Non-Stick with Fin Guide
EUS03-235FG	1.0mm Tip, Non-Stick with Fin Guide
EUS03-237FG	1.2mm Tip, Non-Stick with Fin Guide
EUS03-238FG	1.5mm Tip, Non-Stick with Fin Guide

BAYONET KEY FORCEP 22.9CM - REUSABLE



- US 2-Pin Style Connection
- Keyfin design available to assure the alignment of tines and prevent scissoring action
- Non-Stick Silver Alloy Tips
- Ergonomic Handles
- Available with angled tips at 20, 30 and 45 degrees
- Made in USA

Order Code	Product
EUS03-240NS	0.4mm Tip, Non-Stick
EUS03-243NS	0.7mm Tip, Non-Stick
EUS03-245NS	1.0mm Tip, Non-Stick
EUS03-247NS	1.2mm Tip, Non-Stick
EUS03-248NS	1.5mm Tip, Non-Stick

Order Code	Product
EUS03-240FG	0.4mm Tip, Non-Stick with Fin Guide
EUS03-243FG	0.7mm Tip, Non-Stick with Fin Guide
EUS03-245FG	1.0mm Tip, Non-Stick with Fin Guide
EUS03-247FG	1.2mm Tip, Non-Stick with Fin Guide
EUS03-248FG	1.5mm Tip, Non-Stick with Fin Guide

See page 26 for Bipolar Cables.

Stingray Reusable Bipolar Forceps

Premium Bipolar Forceps, made in USA

POTTS SMITH 17.8CM STRAIGHT BIPOLAR FORCEP - REUSABLE

- US 2-Pin Style Connection
- Available in Stainless Steel or Non-Stick Silver Alloy Tips
- Coated with premium nylon insulation
- Made in USA




Order Code	Product
EUS03-300	0.5mm Tip, Stainless Steel
EUS03-302	0.7mm Tip, Stainless Steel
EUS03-303	1.0mm Tip, Stainless Steel
EUS03-304	1.5mm Tip, Stainless Steel

Order Code	Product
EUS03-300NS	0.5mm Tip, Non-Stick Silver Alloy
EUS03-302NS	0.7mm Tip, Non-Stick Silver Alloy
EUS03-303NS	1.0mm Tip, Non-Stick Silver Alloy
EUS03-304NS	1.5mm Tip, Non-Stick Silver Alloy

SEMKIN 15.2CM STRAIGHT BIPOLAR FORCEP - REUSABLE

- US 2-Pin Style Connection
- Available in Stainless Steel or Non-Stick Silver Alloy Tips
- Coated with premium nylon insulation
- Available with a stopper to prevent scissoring action - Add 'S' to the code to order with this feature
- Made in USA



Order Code	Product
EUS03-368	0.25mm Tip, Stainless Steel
EUS03-370	0.5mm Tip, Stainless Steel
EUS03-371	0.7mm Tip, Stainless Steel
EUS03-372	1.2mm Tip, Stainless Steel
EUS03-374	1.5mm Tip, Stainless Steel

Order Code	Product
EUS03-368NS	0.25mm Tip, Non-Stick Silver Alloy
EUS03-370NS	0.5mm Tip, Non-Stick Silver Alloy
EUS03-371NS	0.7mm Tip, Non-Stick Silver Alloy
EUS03-372NS	1.2mm Tip, Non-Stick Silver Alloy
EUS03-374NS	1.5mm Tip, Non-Stick Silver Alloy

See page 26 for Bipolar Cables.

Stingray Reusable Bipolar Forceps

Premium Bipolar Forceps, made in USA

ADSON 12.7CM BIPOLAR FORCEPS - REUSABLE




- US 2-Pin Style Connection
- Available in Stainless Steel or Non-Stick Silver Alloy Tips
- Coated with premium nylon insulation
- Made in USA

Order Code	Product
EUS03-376	0.4mm Tip, Stainless Steel
EUS03-377	0.7mm Tip, Stainless Steel
EUS03-380	1.0mm Tip, Stainless Steel
EUS03-381	1.2mm Tip, Stainless Steel
EUS03-382	1.5mm Tip, Stainless Steel
EUS03-084	2.0mm Tip, Stainless Steel

Order Code	Product
EUS03-376NS	0.4mm Tip, Non-Stick Silver Alloy
EUS03-377NS	0.7mm Tip, Non-Stick Silver Alloy
EUS03-380NS	1.0mm Tip, Non-Stick Silver Alloy
EUS03-381NS	1.2mm Tip, Non-Stick Silver Alloy
EUS03-382NS	1.5mm Tip, Non-Stick Silver Alloy
EUS03-084NS	2.0mm Tip, Non-Stick Silver Alloy

BAYONET 15.5CM BIPOLAR FORCEP - REUSABLE



- US 2-Pin Style Connection
- Available in Stainless Steel or Non-Stick Silver Alloy Tips
- Coated with premium nylon insulation
- Available with irrigation - Add 'R' to the code to order this feature
- Made in USA

Order Code	Product
EUS03-403	0.5mm Tip, Stainless Steel
EUS03-404	0.7mm Tip, Stainless Steel
EUS03-405	1.0mm Tip, Stainless Steel
EUS03-406	1.2mm Tip, Stainless Steel

Order Code	Product
EUS03-403NS	0.5mm Tip, Non-Stick Silver Alloy
EUS03-404NS	0.7mm Tip, Non-Stick Silver Alloy
EUS03-405NS	1.0mm Tip, Non-Stick Silver Alloy
EUS03-406NS	1.2mm Tip, Non-Stick Silver Alloy


See page 26 for Bipolar Cables.

Stingray Reusable Bipolar Forceps

Premium Bipolar Forceps, made in USA

BAYONET 19.1CM BIPOLAR FORCEP - REUSABLE

- US 2-Pin Style Connection
- Available in Stainless Steel or Non-Stick Silver Alloy Tips
- Coated with premium nylon insulation
- Available with irrigation - Add 'R' to the code to order this feature
- Made in USA



Order Code	Product
EUS03-410	0.25mm Tip, Stainless Steel
EUS03-412	0.5mm Tip, Stainless Steel
EUS03-414	0.7mm Tip, Stainless Steel
EUS03-416	1.0mm Tip, Stainless Steel
EUS03-418	1.2mm Tip, Stainless Steel
EUS03-420	1.5mm Tip, Stainless Steel
EUS03-421	2.0mm Tip, Stainless Steel

Order Code	Product
EUS03-410NS	0.25mm Tip, Non-Stick Silver Alloy
EUS03-412NS	0.5mm Tip, Non-Stick Silver Alloy
EUS03-414NS	0.7mm Tip, Non-Stick Silver Alloy
EUS03-416NS	1.0mm Tip, Non-Stick Silver Alloy
EUS03-418NS	1.2mm Tip, Non-Stick Silver Alloy
EUS03-420NS	1.5mm Tip, Non-Stick Silver Alloy
EUS03-421NS	2.0mm Tip, Non-Stick Silver Alloy

BAYONET 21CM BIPOLAR FORCEP - REUSABLE

- US 2-Pin Style Connection
- Available in Stainless Steel or Non-Stick Silver Alloy Tips
- Coated with premium nylon insulation
- Available with irrigation - Add 'R' to the code to order this feature
- Made in USA



Order Code	Product
EUS03-423	0.25mm Tip, Stainless Steel
EUS03-425	0.5mm Tip, Stainless Steel
EUS03-427	0.7mm Tip, Stainless Steel
EUS03-429	1.0mm Tip, Stainless Steel
EUS03-431	1.2mm Tip, Stainless Steel
EUS03-433	1.5mm Tip, Stainless Steel
EUS03-435	2.0mm Tip, Stainless Steel


Order Code	Product
EUS03-423NS	0.25mm Tip, Non-Stick Silver Alloy
EUS03-425NS	0.5mm Tip, Non-Stick Silver Alloy
EUS03-427NS	0.7mm Tip, Non-Stick Silver Alloy
EUS03-429NS	1.0mm Tip, Non-Stick Silver Alloy
EUS03-431NS	1.2mm Tip, Non-Stick Silver Alloy
EUS03-433NS	1.5mm Tip, Non-Stick Silver Alloy
EUS03-435NS	2.0mm Tip, Non-Stick Silver Alloy

See page 26 for Bipolar Cables.

Stingray Reusable Bipolar Forceps

Premium Bipolar Forceps, made in USA


MCPHERSON 9.4CM STRAIGHT FORCEP - REUSABLE



- US 2-Pin Style Connection
- Available in Stainless Steel and Non-Stick Silver Alloy Tips
- Coated with premium nylon insulation
- Premium lightweight design
- Made in USA

Order Code	Product
EUS03-600	0.4mm Tip, Stainless Steel
EUS03-600NS	0.4mm Tip, Non-Stick Silver Alloy

MCPHERSON 9.4CM ANGLED FORCEP - REUSABLE



- US 2-Pin Style Connection
- Available in Stainless Steel and Non-Stick Silver Alloy Tips
- Coated with premium nylon insulation
- Premium lightweight design
- Made in USA

Order Code	Product
EUS03-601	0.4mm Tip, Stainless Steel
EUS03-601NS	0.4mm Tip, Non-Stick Silver Alloy

JEWELLER 11.9CM STRAIGHT FORCEPS - REUSABLE



- US 2-Pin Style Connection
- Coated with premium nylon insulation
- Premium lightweight design
- Made in USA

Order Code	Product
EUS03-628	0.25mm Tip, Stainless Steel
EUS03-630	0.4mm Tip, Stainless Steel
EUS03-632	0.7mm Tip, Stainless Steel
EUS03-628NS	0.25mm Tip, Non-Stick Silver Alloy
EUS03-630NS	0.4mm Tip, Non-Stick Silver Alloy
EUS03-632NS	0.7mm Tip, Non-Stick Silver Alloy
EUS03-630T	0.4mm, Titanium

See page 26 for Bipolar Cables.


Stingray Eco-Series Reusable Bipolar Forceps

Premium Bipolar Forceps, made in USA

ECO-SERIES 19.1CM BAYONET BIPOLAR FORCEP - REUSABLE

- US 2-Pin Style Connection
- Non-stick, silver alloy tips
- Pin Guide for ultimate tip alignment
- Coated with premium nylon insulation
- Lightweight design
- Made in USA


Order Code	Product
EUS03-842EF	0.5mm Tip, Non-Stick Silver Alloy
EUS03-845EF	1.0mm Tip, Non-Stick Silver Alloy
EUS03-848EF	1.5mm Tip, Non-Stick Silver Alloy



ECO-SERIES 12.7CM STRAIGHT BIPOLAR FORCEP - REUSABLE

- US 2-Pin Style Connection
- Designed for cost-conscious departments
- Great alternative to Adson or Straight Forceps
- Coated with premium nylon insulation
- Lightweight Design
- Made in USA

Order Code	Product
EUS03-380EF	1.0mm Tip



Bipolar Cables

ZAP™ FLYING LEAD BIPOLAR CABLE - SINGLE USE



- Standard US 2 pin Forceps plug, angled and shrouded
- 3m pliable sheathed cable
- Single use, latex free, sterile easipeel packed
- UOM: Each or 50/Dispenser Box

Order Code	Product
E2021BM3	Single Use cable, flying leads

ZAP™ MOLDED PLUG BIPOLAR CABLE - SINGLE USE



- Standard US 2 pin Forceps plug, angled and shrouded
- 3m pliable sheathed cable
- Single use, latex free, sterile easipeel packed
- UOM: Each or 50/Dispenser Box

Order Code	Product
E2022XP3	Single Use cable, molded plug

ELECTRISURE™ FLYING LEADS BIPOLAR CABLE - REUSABLE



- Flying leads and fixed connector options
- Standard US 2 pin, slimline angled, US plug
- UK flat pins forceps plug type available
- 3m pliable silicone sheathed cable
- Reusable, latex free, sterilize at 134°C
- UOM: Each

Order Code	Product
E52443UK	Silicone cable, UK flat pins forceps plug
E52443US	Silicone cable, US 2-Pin Plug

ELECTRISURE™ MOLDED PLUG BIPOLAR CABLE - REUSABLE



- Silicone 3m, 2 Pin Plug, Molded Plug

Order Code	Product
EV52243US	Silicone cable, molded plug

Introducing the SOLO™

Single Use Laparoscopic Electrode Range!

Extensive Tip Designs

L-Hook, J-Hook, Spatula, Ball, Needle, Cone, Knife and integrated suction designs available to allow your entire requirement to be transferred to Single Use



Versatile

A variety of universal cable connections, that can be used with you current cables, or through separate ports to current laparoscopic instruments you are using




Quality Performance

Solid 4.5mm diameter 'One-Piece' stainless steel shaft for exceptional heat transfer and performance

Tough Protection

Hard wearing 0.25mm baked epoxy insulation providing prevention of scratches from trocars and assured protection throughout the surgical field.



Disposable Laparoscopic Electrodes

SOLO™ Laparoscopic Electrodes, for pristine tip condition with consistent outcomes

SOLO™ L-HOOK LAPAROSCOPIC ELECTRODE



- Moulded handle with contours for positive and comfortable grip
- High quality baked epoxy insulated shaft
- Solid one-piece stainless steel shaft for exceptional heat transfer and performance
- Latex free, sterile packed

Order Code	Product
E5050AB	36cm, 3m Cable, 4mm Generator Plug
E5051DF	36cm, 3m Cable, 8mm Generator Plug
E5052BH	36cm, no cable

Order Code	Product
E5055FW	45cm, 3m Cable, 4mm Generator Plug
E5056PL	45cm, 3m Cable, 8mm Generator Plug
E5057KK	45cm, no cable

SOLO™ J-HOOK LAPAROSCOPIC ELECTRODE



- Moulded handle with contours for positive and comfortable grip
- High quality baked epoxy insulated shaft
- Solid one-piece stainless steel shaft for exceptional heat transfer and performance
- Latex free, sterile packed

Order Code	Product
E5060EJ	36cm, 3m Cable, 4mm Generator Plug
E5061HB	36cm, 3m Cable, 8mm Generator Plug
E5062DC	36cm, no cable

Order Code	Product
E5065TR	45cm, 3m Cable, 4mm Generator Plug
E5066GV	45cm, 3m Cable, 8mm Generator Plug
E5067FD	45cm, no cable

SOLO™ SPATULA LAPAROSCOPIC ELECTRODE



- Moulded handle with contours for positive and comfortable grip
- High quality baked epoxy insulated shaft
- Solid one-piece stainless steel shaft for exceptional heat transfer and performance
- Latex free, sterile packed

Order Code	Product
E5070HV	36cm, 3m Cable, 4mm Generator Plug
E5071TM	36cm, 3m Cable, 8mm Generator Plug
E5072AE	36cm, no cable

Order Code	Product
E5075KL	45cm, 3m Cable, 4mm Generator Plug
E5076QX	45cm, 3m Cable, 8mm Generator Plug
E5077EC	45cm, no cable

Disposable Laparoscopic Electrodes

SOLO™ Laparoscopic Electrodes, for pristine tip condition with consistent outcomes

SOLO™ BALL LAPAROSCOPIC ELECTRODE

- Moulded handle with contours for positive and comfortable grip
- High quality baked epoxy insulated shaft
- Solid one-piece stainless steel shaft for exceptional heat transfer and performance
- Latex free, sterile packed



Order Code	Product
E5080PB	36cm, 3m Cable, 4mm Generator Plug
E5081JJ	36cm, 3m Cable, 8mm Generator Plug
E5082UY	36cm, no cable

Order Code	Product
E5085UC	45cm, 3m Cable, 4mm Generator Plug
E5086WK	45cm, 3m Cable, 8mm Generator Plug
E5087EU	45cm, no cable

SOLO™ NEEDLE LAPAROSCOPIC ELECTRODE

- Moulded handle with contours for positive and comfortable grip
- High quality baked epoxy insulated shaft
- Solid one-piece stainless steel shaft for exceptional heat transfer and performance
- Latex free, sterile packed



Order Code	Product
E5090EM	36cm, 3m Cable, 4mm Generator Plug
E5091CT	36cm, 3m Cable, 8mm Generator Plug
E5092CH	36cm, no cable

Order Code	Product
E5095GQ	45cm, 3m Cable, 4mm Generator Plug
E5096GN	45cm, 3m Cable, 8mm Generator Plug
E5097GT	45cm, no cable

SOLO™ CONE LAPAROSCOPIC ELECTRODE

- Moulded handle with contours for positive and comfortable grip
- High quality baked epoxy insulated shaft
- Solid one-piece stainless steel shaft for exceptional heat transfer and performance
- Latex free, sterile packed



Order Code	Product
E5100NY	36cm, 3m Cable, 4mm Generator Plug
E5101XA	36cm, 3m Cable, 8mm Generator Plug
E5102BB	36cm, no cable

Order Code	Product
E5105ZL	45cm, 3m Cable, 4mm Generator Plug
E5106UV	45cm, 3m Cable, 8mm Generator Plug
E5107DT	45cm, no cable

Disposable Laparoscopic Electrodes

SOLO™ Laparoscopic Electrodes, for pristine tip condition with consistent outcomes

SOLO™ KNIFE LAPAROSCOPIC ELECTRODE



- Moulded handle with contours for positive and comfortable grip
- High quality baked epoxy insulated shaft
- Solid one-piece stainless steel shaft for exceptional heat transfer and performance
- Latex free, sterile packed

Order Code	Product
E5110XR	36cm, 3m Cable, 4mm Generator Plug
E5111RF	36cm, 3m Cable, 8mm Generator Plug
E5112FQ	36cm, no cable

Order Code	Product
E5115RM	45cm, 3m Cable, 4mm Generator Plug
E5116TY	45cm, 3m Cable, 8mm Generator Plug
E5117GP	45cm, no cable

SOLO™ HAND-SWITCHING LAPAROSCOPIC ELECTRODES



- Moulded handle with contours for positive and comfortable grip
- High quality baked epoxy insulated shaft
- Solid one-piece stainless steel shaft for exceptional heat transfer and performance
- Cut and Coag push button control
- Remove unnecessary cables from your operating theatre

Order Code	Product
E5053LW	L-Hook, 33cm
E5058EB	L-Hook, 45cm
E5063KG	J-Hook, 33cm
E5068RP	J-Hook, 45cm

Order Code	Product
E5073TN	Spatula, 33cm
E5078ED	Spatula, 45cm
E5093VG	Needle, 33cm
E5098RK	Needle, 45cm

Diathermy Pencils

ELECTRISURE™ DIATHERMY PENCIL - SINGLE USE

- 3m cable with universal 3 pin generator plug
- Supplied with Teflon coated non-stick blade electrode
- Single use, latex free, sterile easipeel packed
- UOM: Each or 40/Dispenser Box



Order Code	Product
E102PP	Diathermy Pencil, hand switching

ELECTRISURE™ DIATHERMY PENCIL - REUSABLE

- 3m cable with universal 3 pin generator plug
- Blade Electrode
- Hand switching
- Reusable
- Sterilize at 134°C
- UOM: Each



Order Code	Product
ES104141HS	Diathermy Pencil, hand switching

ELECTRISURE™ DIATHERMY PENCIL - REUSABLE

- 3m cable
- 4mm Banana Plug
- Blade Electrode
- Foot Control
- Reusable
- Sterilize at 134°C
- UOM: Each



Order Code	Product
ES104814FC	Diathermy Pencil, foot switching

Monopolar Electrodes

ZAP™ BLADE ELECTRODES - SINGLE USE



- Universal Hex connection for monopolar pencils
- One-piece stainless steel core for premium heat transfer
- Teflon Coated for non-stick performance
- Sterile easipeel packed
- UOM: Each or 25/Dispenser box

Order Code	Product
E6021070HNT	Standard 7cm
E6023070DVT	Standard 7cm Fine
E6025070NMT	Modified 7cm

ZAP™ NEEDLE ELECTRODES - SINGLE USE



- Universal Hex connection for monopolar pencils
- One-piece stainless steel core for premium heat transfer
- Teflon Coated for non-stick performance
- Sterile easipeel packed
- UOM: Each or 25/Dispenser box

Order Code	Product
E6031070HHT	Standard 7cm
E6032070MYT	Modified 7cm

ZAP™ BALL ELECTRODES - SINGLE USE



- Universal Hex connection for monopolar pencils
- One-piece stainless steel core for premium heat transfer
- Teflon Coated for non-stick performance
- Sterile easipeel packed
- UOM: Each or 25/Dispenser box

Order Code	Product
E6045055DZ	5.5cm x 5mm diameter
E6045130DAT	13cm x 5mm diameter

Monopolar Electrodes

ZAP™ ROUND LOOP ELECTRODES - SINGLE USE

- Universal Hex connection for monopolar pencils
- 5.5cm Length
- Sterile easipeel packed
- UOM: Each or 25/Dispenser box

Order Code	Product
E60605506KQ	6mm diameter
E60605508TL	8mm diameter
E60605510MP	10mm diameter



ELECTRISURE™ MONOPOLAR ELECTRODES - REUSABLE

- Universal Hex connection for monopolar pencils
- One-piece stainless steel core for premium heat transfer
- Reusable
- Autoclavable to 134°C
- UOM: Each

Order Code	Product
E6021070HNR	Blade 7cm
E6031070HHR	Needle 7cm
E6045055DZR	Ball 7cm x 5mm



ZAP™ LLETZ LOOP ELECTRODES - SINGLE USE

- 13cm Length, Universal Hex connection for monopolar pencils
- High grade tungsten wire
- Sterile easipeel packed
- UOM: Each or 25/Dispenser box

Order Code	Product
E60511010PT	10x10mm
E60511510WA	15x10mm
E60511512EU	15x12mm
E60511515BV	15x15mm
E60512008NR	20x8mm
E60512010BB	20x10mm
E60512012ZB	20x12mm
E60512015BG	20x15mm
E60512020PR	20x20mm



Patient Return Pad Electrodes

PATIENT ELECTRODE RETURN PAD - SINGLE USE

- Self-adhesive with conduction gel
- Adult 129cm² area
- 3m pliable cable with Valleylab type generator plug
- Single use, latex free, sterile Aluminium foil packed
- UOM: Each or 25/Dispenser Box



Order Code	Product
E4022HH	Adult
E4023QD	Paediatric
E4024VX	Neonatal

PATIENT ELECTRODE RETURN PAD, WITHOUT CABLE - SINGLE USE

- Self-adhesive with hydro gel and foam support
- Adult 132cm² area
- No cable, connects to reusable patient return cables
- UOM: Each or 25/Dispenser Box



Order Code	Product
E4016TK	Adult
E4017BZ	Paediatric
E4018PH	Neonatal

Monopolar Cables

ZAP™ MONOPOLAR CABLES - SINGLE USE

- 4mm or 8mm pin generator plug connector
- 4mm female Forcep plug
- 3m pliable sheathed cable
- Single use, latex free, sterile easipeel packed
- UOM: Each or 10/Dispenser Box

Order Code	Product
E2200143	4mm generation pin, 4mm female banana plug
E2200243	8mm generator plug, 4mm female banana plug



ELECTRISURE™ MONOPOLAR CABLES - REUSABLE

- 8mm generator plug connector
- 3m pliable Silicone sheathed cable
- Reusable, latex free, sterilize at 134°C
- UOM: Each

Order Code	Product
E598206ST	5/32" screw thread plug
E598340BP	4mm female banana plug



Electrosurgical Accessories

DIATHERMY QUIVERS



- Screw off base for easy cleaning
- Stainless Steel holding ring
- High grade medical plastic - will not split, crack or powderise
- Belled top design for easy instrument placement
- Reusable, autoclavable at 134°C
- UOM: Each

Order Code	Product
EDQS15	150mm
EDQS20	200mm
EDQS25	250mm
EDQS30	300mm
EDQS35	350mm
EDQS40	400mm
EDQSC	Quiver end cap only
EDQSR	Quiver ring only

ELECTROCAUTERIES - SINGLE USE



- Fast and effective precision cauterization
- Designed for wide variety of applications
- Battery powered and cordless
- Available in varying temperatures and tip sizes
- Single use, sterile packed
- UOM: Each or 10/Box

Order Code	Product
EF7233	1200°C Semi-Extended Fine 50mm Tip
EF7234	1200°C Extended Fine 75mm Tip
EF7244	1200°C Fine 28mm Tip
EF7255	600°C Fine 28mm Tip
EF7266	1200°C Large 28mm Tip
EF7288	800°C Thick 28mm Tip

Electrosurgical Accessories

FIAB NERVE STIMULATOR LOCATOR

- Designed to identify exposed motor nerves during microsurgery operations
- Special applications in plastic surgery and head and neck surgery
- Supplied with an active electrode and subcutaneous reference needle
- Alternating and pulsating current guarantees safe and effective stimulation
- Single Use, Sterile Packed
- UOM: Each or 10/Box.

Order Code

EF1745



ALSA Diathermy Generators

Multi-functional Diathermy Generators for minor surgical procedures.

ALSA 80W ELECTROSURGICAL GENERATOR



- 80 watt output
- Footswitch activation
- 9 monopolar functions; Bipolar function
- Digital control panel
- Touchpad switches for easy cleaning
- Supplied with Monopolar Kit ESB700/A
- Made in Italy

Order Code	Product
ESU50MPC	Multi Function 80W

ALSA MONOPOLAR ACCESSORY KIT



- 1 x footswitch pencil with 3m cable
- 1 x patient return plate with 3m cable
- 1 x silicone patient strap
- 6 x electrodes including blades, needles and balls
- All items are reusable and autoclaveable.

Order Code	Product
ESB700/A	Suits Generator ESU50-MPC

ALSA 140W ELECTROSURGICAL GENERATOR



- 140 watt output
- Handswitch and Footswitch activation
- 9 monopolar functions
- Bipolar function
- Digital control panel
- Touchpad switches for easy cleaning
- Supplied with Monopolar Kit ESB730/A
- Made in Italy

Order Code	Product
ESU140MPC	Multi Function 140W

ALSA Diathermy Generator Accessories

ALSA MONOPOLAR ACCESSORY KIT

- 1 x handswitch pencil with 3m cable
- 1 x patient return plate with 3m cable
- 1 x silicone patient strap
- 6 x electrodes including blades, needles and balls
- All items are reusable and autoclaveable.



Order Code	Product
ESB730/A	Suits Generator ESU140-MPC

BIPOLAR ACCESSORIES - REUSABLE

- See pages 12-25 for a complete range of Bipolar Forcep options

Order Code	Product
EV52243US	Bipolar cable, 3m, Silicone
EUS03-380EF	Stingray Adson bipolar forcep; 12cm x 1.0mm tips



ELECTRISURE™ MONOPOLAR ELECTRODES - REUSABLE

- Universal Hex connection for monopolar pencils
- One-piece stainless steel core for premium heat transfer
- Reusable
- Autoclavable to 134°C
- UOM: Each

Order Code	Product
E6021070HNR	Blade 7cm
E6031070HHR	Needle 7cm
E6045055DZR	Ball 7cm x 5mm



ALSA Diathermy Generator Accessories

MONOPOLAR PENCILS - REUSABLE



Order Code	Product
ES104141HS	Hand switching pencil
ES801364	Foot switching pencil

PATIENT RETURN PAD ELECTRODES



Order Code	Product
ES801274	Stainless Steel patient return pad electrode with 3m cable, reusable
ES449195	Patient return cable only, 3m, reusable
ES712353	Patient return plate only, stainless steel, reusable
ES800842	Silicone fixing belt, reusable
E801147	Patient return pad, self-adhesive, without cable, single use

ALSA FOOTSWITCH



Order Code	Product
ES900761	Alsa footswitch for ESU50MPC machine

Surgical Accessories



SURGICAL ACCESSORIES

Visceral Retainer	42
Magnetic Mats	42
Light Handle Cover	43
Lubricating Gel	43
Jumbo Swabs	43
Gowns & Aprons	44-45
Masks & Face-Shields	46
Gloves	47
Hand Sanitizer	47

Surgical Accessories

THE STURGEON FISH VISCERAL RETAINER



- Mechanical barrier to reduce the risk of needle perforation injury to the bowel.
- Impregnated with barium sulphate to render it radiopaque for additional safety.
- Surgical grade silicone
- Autoclaveable, latex free, reusable
- UOM: Each

Order Code	Product
R3020MVR	Standard, 300 x 130 x 1.5mm
R3021FVR	Large, 350 x 156 x 2.0mm

SOLO™ MAGNETIC INSTRUMENT MAT



- Raised magnets with non-slip foam backing
- Single use, latex free, sterile easipeel packed
- UOM: Each, 10/Dispenser Box

Order Code	Product
G3013EP	Small, 25cm x 40cm
G3014KS	Large, 50cm x 40cm

SOLO™ MAGNETIC INSTRUMENT MAT WITH DEMAGNETIZED ZONE



- With central demagnetised zone for holding fine instruments and sharps
- Raised magnets with non-slip foam backing
- Single use, latex free, sterile easipeel packed
- UOM: Each, 10/Dispenser Box

Order Code	Product
G3011MM	Small, 25 x 40cm, Demag Zone

Surgical Accessories

SOLO™ LIGHT HANDLE COVER

- Universal Light Handle cover suiting most light handles
- Large Foam Hand Guard for easy fitting and removal without touching the light handle
- Individually sterile easipeel packed
- UOM: Each, 50/Dispenser Box, 400/Outer

Order Code

G3001LC



AGEL LUBRICATING GEL - 500ML

- Water soluble, greaseless, transparent, non-staining, non-irritating
- High viscosity for excellent adherence
- Latex free, clean packed
- UOM: Each

Order Code

G100100



SURGIGEL PLUS LUBRICATING GEL SACHETS

- Sterile, latex free 2.7g Satchet
- Water soluble, greaseless, transparent, non-staining, non-irritating
- High viscosity for excellent adherence
- UOM: 144/Dispenser Box

Order Code

G100120



SOLO™ JUMBO SWAB COTTON APPLICATOR

- Single use, latex free, clean packed
- UOM: 20/Pack

Order Code

G054200JU



SOLO™ Isolation Gown

Velcro neck closure for easy fitting and removal



Long waist ties for secure fitting and complete protection



Durable PP Material with splashproof PE coating for excellent protection



Comfortable knitted cuffs designed for assistance when donning and doffing

One size fits all.
Length: 114cm;
Width: 153cm;
Arm Length: 56cm

HallmarkSurgical™
When you care.

Surgical Accessories

SOLO™ ISOLATION GOWNS

- PP + PE Isolation Gown, 35gsm with excellent water repellent features
- Knitted Cuffs for enhanced donning and doffing
- Waist Ties and Velcro Neck Closure for increase security, user comfort and convenience
- UOM: 50/Carton

Order Code

G100131GN



SOLO™ PE THUMBLOOP GOWNS

- Durable Polyethylene construction
- Convenient Thumbloop design for enhanced donning and doffing
- Waterproof
- Waist Ties and neck loop for great protection level
- UOM: 200/Carton

Order Code

G100132CP



SOLO™ DISPOSABLE APRONS

- Durable polyethylene (55mu)
- 106 x 69cm size
- 25 Aprons/cardboard header for easy mounting and tear off dispensing
- Latex Free, Clean-packed
- UOM: 25/Header Card, 1000/Carton

Order Code

G100130PA



Surgical Accessories

SOLO™ 3-PLY EARLOOP MASKS



- 3-Ply General Use Mask
- Ultrasonic Welded Construction for durable long-term use
- Strong comfortable earloops
- Nose-wire seal designed to suit all users
- UOM: 50/Carton

Order Code

G2104BM

SOLO™ LEVEL 3 SURGICAL MASK



- ASTM Level 3 Surgical Mask
- Ultrasonic Welded Construction for durable long-term use
- Bacterial Filtration Efficiency (BFE): >99%
- Particle Filtration Efficiency (PFE): >99% @ 0.1 micron
- High Fluid Resistance: 160mmHG
- Bacterial Filtration Efficiency (BFE): >99%
- Breathability - Delta P: <3.0mmH2O/cm²
- Flammability: Class 1
- UOM: 40/Carton

Order Code Product

G2107HD Earloop

G2108TM Tie-Back

SOLO™ ANTIFOG FACE-SHIELD



- Anti-Fog coating on the inside of the shield guarantees high visibility during use
- Soft open PU foam gives comfortable contact against your forehead
- 33cm wide x 19cm Long shield
- Made in NZ
- UOM: 60/Carton

Order Code

G266JB

Surgical Accessories

BASTION LATEX MEDICAL EXAMINATION GLOVES, LIGHTLY POWDERED

- Natural biodegradable rubber
- Greatest material durability and flexibility
- Lightly powdered to aid fitting
- Superior tactile sensitivity, comfort, and fit
- UOM: 100/Dispenser Box

Order Code	Product
G100151	Small
G100152	Medium
G100153	Large
G100154	Extra Large



BASTION LATEX MEDICAL EXAMINATION GLOVES, POWDER FREE

- Natural biodegradable rubber
- Greatest material durability and flexibility
- Superior tactile sensitivity, comfort, and fit
- UOM: 100/Dispenser Box

Order Code	Product
G100161	Small
G100162	Medium
G100163	Large
G100164	Extra Large



INSTANT HAND SANITIZER - 500ML

- 70% Ethyl Alcohol with Biocide Formulation
- Contains Aloe extract for moisturizing your skin
- Effective in killing up to 99.9% of common germs and bacteria
- Manufactured in Singapore for hospital and professional use
- UOM: 12/Carton

Order Code
G100102



Surgical Suction



SURGICAL SUCTION

Suction Liners	50-51
Frazier Suctions	52-53
Fukushima Teardrop Suctions	54-55
Specialty Suctions	56-60

SOLO™ Suction Liner System

World Class innovative suction management system for every level of healthcare including OR, Endoscopy, ICU, Outpatients, Clinics and Wards

Fluid activated shut-off valve for safety during procedures



Integrated suction connection on the liners providing lowest set-up time among all suction liners

Highly flexible, durable and compact suction liners reducing storage space and environmental waste

Robust, strong canister's available in 1L, 1.5L, 2L and 3L sizes, with measurements readable from 270°

'Series' installation enables automatic switch over to multiple suction canisters for high fluid applications

Available to be installed in Series, Rail Mounted, Wall mounted or trolley mounted for your perfect set-up

Suction Liners

Compact, fast and durable suction liner system

SUCTION CANISTERS

- Durable, long-lasting reusable canister
- Wide measuring scale visible from 270 degrees for easy reading
- Built-in brackets on canister provides easy mounting options
- Available in 4 sizes
- Large range of connection options available for your ideal suction system setup
- Reusable UOM: Each, 12 or 25/Carton

Order Code	Product
H7100JK	1000ml
H7101WL	1500ml
H7102RH	2000ml
H7103DV	3000ml



SUCTION LINERS


- Super durability suction liner providing exceptional safety.
- Compact soft design meaning less storage space in your department
- Easy single handed removal of suction liners
- Fast setup - do not need to line up suction port on canister allowing for immediate insertion and connection
- Fastest liner in the market to reach 80kPa pressure
- Available in 4 sizes
- Single Use, Clean-packed
- UOM: Each, 70, 90 or 100/Carton

Order Code	Product
H7110DW	1000ml
H7111VP	1500ml
H7112DG	2000ml
H7113UA	3000ml



New Zealand's Premium Frazier Suction

Available with polished, atraumatic or fenestrated tips depending on your surgical need



Rigid stainless steel lumen which doubles as a probe during surgery



Precise control through the ergonomically shaped handle and suction port



Non-blocking handle design - The aperture of the channel through the handle is significantly larger than the cannula, ensuring no blockages occur



Manufactured in the UK

HallmarkSurgical™
When you care.

Surgical Suction Cannula

Eliminate risk with single use Suctions but with no compromise in quality, control and performance of Suction.

FRAZIER SUCTION - POLISHED TIP

- Premium quality suction with strong Stainless Steel lumen which also acts as a probe.
- Surgical grade ABS moulded handle, non-blocking design
- Precise suction control port
- Handle ergonomically shaped for firm grip and control
- Made in UK
- Single use, latex free, sterile easipeel packed
- UOM: Each, 25/Dispenser Box


Order Code	Product
H3106F	6fg x 9cm polished tip
H311613BE	6fg x 13cm polished tip
H320909F	9fg x 9cm polished tip
H320913L	9fg x 13cm polished tip
H33012F	12fg x 9cm polished tip
H3301213TH	12fg x 13cm polished tip



FRAZIER SUCTION - ATRAUMATIC TIP

- Rounded cannula tip providing atraumatic suction and probing at the surgical site.
- Premium quality suction with strong Stainless Steel lumen which also acts as a probe.
- Surgical grade ABS moulded handle, non-blocking design
- Precise suction control port
- Handle ergonomically shaped for firm grip and control
- Made in UK
- Single use, latex free, sterile easipeel packed
- UOM: Each, 25/Dispenser Box


Order Code	Product
H3116FA	6fg x 9cm atraumatic tip
H3219FA	9fg x 9cm atraumatic tip
H321913KZ	9fg x 13cm atraumatic tip
H33112FA	12fg x 9cm atraumatic tip
H3311213NB	12fg x 13cm atraumatic tip



SOLO™ Fukushima Teardrop Suctions

Premium quality suction with rigid stainless steel lumen which can be used as a probe during surgery

Surgical grade ABS moulded handle, ergonomically shaped for firm control and grip



Available in 6fg, 9fg and 12fg
- Colour-coded for easy size identification

Precise Teardrop suction control port providing unparalleled suction control for delicate Spinal and Neuro procedures

Designed and
Manufactured in Italy

Cleaning Stylet included

HallmarkSurgical™
When you care.

Surgical Suction Cannula

Eliminate risk with single use Suctions but with no compromise in quality, control and performance of Suction.

FRAZIER SUCTION - FENESTRATED TIP

- Premium quality suction with strong Stainless Steel lumen which also acts as a probe.
- Surgical grade ABS moulded handle, non-blocking design
- Precise suction control port
- Handle ergonomically shaped for firm grip and control
- Made in UK
- Single use, latex free, sterile easipeel packed
- UOM: Each, 25/Dispenser Box

Order Code	Product
H3120609KK	6fg x 9cm, fenestrated tip
H3230909RQ	9fg x 9cm, fenestrated tip



SOLO™ FUKUSHIMA SUCTION - TEARDROP CONTROL PORT

- Teardrop shaped suction control port for precise suction adjustment and control
- Especially suited for Spinal and Neuro Surgery
- Cleaning stylet included
- Ergonomic handle design
- Single use, latex free, sterile easipeel packed
- UOM: Each, 25/Dispenser Box

Order Code	Product
H401006HG	6fg x 16cm
H401009RL	9fg x 16cm
H401012NF	12fg x 16cm



Surgical Suction Cannula

Eliminate risk with single use Suctions but with no compromise in quality, control and performance of Suction.

SINUS FESS SUCTION



- Premium quality suction with strong Stainless-Steel lumen
- Curved tips for maximum visibility and access
- Precise suction control port on some models
- ABS handle ergonomically shaped for firm grip and control
- Made in UK
- Single use, latex free, sterile easipeel packed
- UOM: 25/Dispenser Box

Order Code	Product
H210970SA	9fg x 70 deg. Rounded tip, suction control
H210972TY	9fg x 70 deg. Ball tip, suction control
H210990PT	9fg x 90 deg. Rounded tip, suction control

Order Code	Product
H211270CL	12fg x 70 deg. Rounded tip, suction control
H211272NW	12fg x 70 deg. Ball end, suction control
H211270EN	12fg x 70 deg. Rounded tip without suction control

LONG SUCTION



- Premium quality suction with strong Stainless-Steel lumen
- Surgical grade ABS moulded handle, non-blocking design
- Precise suction control port
- Handle ergonomically shaped for firm grip and control
- Polished tip
- Made in UK
- Single use, latex free, sterile easipeel packed
- UOM: 25/Dispenser Box

Order Code	Product
H3350620MH	6fg x 20cm
H3350627TS	6fg x 27cm
H3350645SL	6fg x 45cm
H3350920NN	9fg x 20cm
H3350927WE	9fg x 27cm
H3350935HV	9fg x 35cm

Order Code	Product
H3350945RT	9fg x 45cm
H3350955YF	9fg x 55cm
H3351220MQ	12fg x 20cm
H3351240YT	12fg x 40cm
H3351255DV	12fg x 55cm

Surgical Suction Cannula

Eliminate risk with single use Suctions but with no compromise in quality, control and performance of Suction.

SIGMOIDOSCOPE SUCTION

- Premium quality suction with strong Stainless Steel lumen
- Non-blocking handle
- Atraumatic tip
- Single use, latex free, sterile easipeel packed
- Made in UK
- UOM: 25/Dispenser Box

Order Code

H3361244SG



SCHUKNECHT SUCTION

- Premium quality suction with strong Stainless-Steel lumen
- Precise suction control port
- ABS handle ergonomically shaped for firm grip and control
- Made in UK
- Single use, latex free, sterile easipeel packed
- UOM: 25/Dispenser Box

Order Code

Product

H30218SG 18g x 6cm

H30220NA 20g x 6cm



Surgical Suction Cannula

Eliminate risk with single use Suctions but with no compromise in quality, control and performance of Suction.

LASER SUCTION



- Laser suction combines two instruments in one
- A 0.9mm fibre optic laser port bonded to the suction tube allows simultaneous aspiration with cutting and coagulation
- For use with 600 micron core laser fibre with 850 micron sheath
- Premium quality suction with strong Stainless Steel lumen
- Surgical grade ABS moulded handle, non-blocking design
- Precise suction control port
- Handle ergonomically shaped for firm grip and control
- Made in UK
- Single use, latex free, sterile easipeel packed
- UOM: 25/Dispenser Box

Order Code	Product
H34100609GW	6fg x 9cm, extended laser channel
H34102FT	6fg x 9cm, flush finish

NERVE ROOT SUCTION



- Premium quality suction with strong Stainless-Steel lumen
- Surgical grade ABS moulded handle, non-blocking design
- Precise suction control port
- Handle ergonomically shaped for firm grip and control
- Curved lip at the tip allows simultaneous retraction and aspiration
- The tube is flattened for easy access to the surgical site
- Made in UK.
- Single use, latex free, sterile easipeel packed
- UOM: 25/Dispenser Box

Order Code	Product
H33706121NR	12fg x 12cm, 60 deg
H33706122RT	12fg x 12cm, 60 deg, reverse tip

Surgical Suction Cannula

Eliminate risk with single use Suctions but with no compromise in quality, control and performance of Suction.

NASAL ASPIRATING DISSECTOR

- Premium quality suction with strong Stainless-Steel lumen
- Surgical grade ABS moulded handle, non-blocking design
- Precise suction control port
- Handle ergonomically shaped for firm grip and control
- Flattened and curved tip allows for ideal access and elevation in the nasal cavity
- Small fenestration on both sides of the tip allows for simultaneous suction
- Made in UK.
- Single use, latex free, sterile easipeel packed
- UOM: 25/Dispenser Box

Order Code

H215157NA



DANDY DRAINAGE CANNULA

- Premium quality Stainless Steel cannula
- Flexible silicone tube connector
- Etched at approximate 1cm markings, as a guide only
- This is not a calibrated instrument
- Single use, latex free, sterile easipeel packed
- UOM: 25/Dispenser Box

Order Code

H401236CD



Surgical Suction Cannula

Eliminate risk with single use Suctions but with no compromise in quality, control and performance of Suction.

TIBBS CANNULA SET



- Premium curved Stainless Steel cannula with atraumatic tips
- For use in vascular surgery, for irrigation of blood vessels
- Unique olive and cone hybrid design, fits a wide range of different sized veins
- Standard luer connection with non-blocking design
- Made in UK
- Single use, sterile packed, latex free
- UOM: 25/Dispenser Box

Order Code

H60113TB

ASPIRATION NEEDLE



- Sharp bevelled tip for safe and controlled penetration when aspirating fluids from cysts and organs
- Luer lock fitting
- Premium quality Stainless Steel cannula
- Made in UK
- Single use, latex free, sterile easipeel packed
- UOM: 25/Dispenser Box

Order Code Product

H5000332VM 3mm OD x 32cm

H5000532US 5mm OD x 32cm

Hand Surgery



HAND SURGERY

Hand Immobilisers	62-64
Finger Traps	65

Introducing FlexiHand

Silicone-encased Lead Hands

Soft touch silicone surface for gentle patient contact

Lead core provides weight and malleability of the traditional Lead Hand

Child and Adult sizes available

Heaviest Hand at 1435g - no extra weighting or Ballast is required

Very flexible for accurate and firm positioning and returning the Hand to the flat position can be done with ease

Reusable

Lower cost per use compared to conventional Lead Hands

High density surgical Silicone protects yourself, the patient and other Instruments from Lead

Resistant to all commonly used cleaning chemicals



SAFER. SMARTER. STRONGER.

HallmarkSurgical™
When you care.

Hand Immobilisers

New Zealand's most comprehensive lead-free hand immobiliser range

SILICONE FLEXI-HAND T™

- 2.2mm Lead core providing weight and malleability of traditional lead hands
- 8mm seamless surgical silicone encasing for extended life-time use and protection from lead
- Up to 30% heavier than traditional lead hands
- Resistant to all commonly used cleaning chemicals for surgical instruments
- Autoclaveable, maximum of 136°C
- Latex free, lead free, reusable

Order Code	Product
N802SL	Adult, 37 x 27cm, 1435g
N803BK	Child, 32 x 23cm, 945g



SILICONE FLEXI-HAND X™

- Malleable Aluminium core encased in Surgical Silicone
- Anti-Rollover wings to prevent, arm and hand rotation during surgery
- Radiolucent for use with fluoroscopy and X-Ray
- Latex free, lead free, reusable
- Autoclavable 134°C

Order Code	Product
N805XR	41cm x 33cm, 640g



Hand Immobilisers

New Zealand's most comprehensive lead-free hand immobiliser range

DURAHAND



- Pliable Aluminium with closed cell epoxy coating
- X-Ray compatible, Latex free, reusable
- Autoclavable 134°C

Order Code	Product
NDH225	Large, 24 x 22cm, 83g
NDH350	X-Large, 33 x 27cm, 132g

DURAHAND - SINGLE USE



- Malleable Aluminium core with Hi-Grip Foam atraumatic covering
- Lead free, Latex free, X-Ray compatible
- Single use. Sterile easipeel packed
- Available in Large and X-Large

Order Code	Product
N760	Large
N761	X-Large

Finger Traps

Nylon and Stainless Steel Reusable Finger Traps.

NYLON FINGER TRAP

- Nylon construction for optimum patient comfort
- Gentle on skin but firm grip.
- Colour coded by size. Washable, reusable.
- Available in 5 Sizes: Extra Small, Small, Medium, Large, Extra Large UOM Each.

Order Code	Size
NK100XS	X-Small
NK200S	Small
NK300M	Medium
NK400L	Large
NK500XL	X-Large



DOUBLE NYLON FINGER TRAP

- Nylon construction for optimum patient comfort
- Gentle on skin but firm grip.
- Colour coded by size. Washable, reusable.
- Available in 2 Sizes: Medium, Large UOM Each.

Order Code	Size
NKD300M	Medium
NKD400L	Large



STAINLESS STEEL FINGER TRAP

- Maximum strength and grip
- Stranded Stainless Steel wire.
- Reusable, sterilisable at 134°C.
- Available in 3 Sizes: Small, Medium, Large UOM Each.

Order Code	Product
NK3280	Small
NK3285	Medium
NK3290	Large



Colorectal & Endoscopy Devices



COLORECTAL DEVICES

SOLO™ Proctoscopes and Anoscopes	68-69
Haemorrhoid Banding	70-72
Sigmoidoscope Systems	73-75
Lighting Accessories	76-79
Endoscopy Accessories	80-81
Bowel Washout Systems	82

SOLO™ Proctoscopes and Anoscopes

Premium European manufacture

SOLO™ SELF-LIGHTING PROCTOSCOPES



- Crystal clear plastic guarantees excellent light transmission
- Inbuilt LED light source for bright white illumination
- Light is precision focused for optimum site visibility
- Light system compatible
- Exceptionally smooth transition from obturator to proctoscope for ease of insertion
- Single use, latex free, clean easipeel packed
- UOM: Each, 30/Dispenser Box

Order Code	Product
P202S	Paediatric 12mm ID/16mm OD
P203M	Standard 18mm ID/22mm OD
P204L	Procedure 22mm ID/26mm OD

SOLO™ Proctoscopes and Anoscopes

Premium European manufacture

SOLO™ ANOSCOPES - BEVELLED EDGE

- Exceptionally smooth transition from obturator to proctoscope for ease of insertion
- Opaque white PEHD material directs light forward for superior illumination
- Most similar design to common reusable Anoscopes
- Single use, latex free, clean easipeel packed
- UOM: Each, 30/Dispenser Box



Order Code	Product
P2322M	Examination/Standard 22mm ID
P2324L	Procedure/Large 24mm ID

SOLO™ ANOSCOPES - STRAIGHT EDGE

- Surgical quality, with Obturator
- Opaque white PEHD material directs light forward for superior illumination
- Most similar design to common reusable Anoscopes
- Single use, latex free, clean easipeel packed
- UOM: Each, 30/Dispenser Box



Order Code	Product
P2326FL	Procedure/Large 24mm ID

Haemoband Pre-loaded Ligator

Pre-loaded with 4 bands
reducing procedure times
by 40-50%!

Slimline bander shaft allows greater
visibility for the surgeon



Individually moulded bands,
with exceptional tensile
strength guaranteeing
strong secure fitting and
eliminating the need to
double-band

Wide internal tip diameter
captures large amount of
tissue for very effective
band placement

Fast positive trigger action
activates suction, fires the
band and re-loads all in
one seamless motion

Available with Suction Tubing
or Proctoscopes

Integrated suction
control through
the trigger for
total control to
the surgeon

Haemorrhoid Banding

Pre-Loaded for fast efficient banding

HAEMOBAND PLUS PRE-LOADED SUCTION BANDER

- Fast, positive trigger activation
- Multi-action trigger that controls suction, fires and reloads the bands
- Pre-loaded with 4 bands for fast procedure time
- Connects to standard suction system
- SOLO™ PVC Tubing available as optional extra - 6mm ID x 9mm OD x 2.0M Length, fitted with universal female connectors
- Single use, latex free, clean easipeel packed
- UOM: Each, 8/Dispenser Box



Order Code

P305T

HAEMOBAND PLUS AND SUCTION TUBING KIT

- Haemoband Plus Pre-Loaded Suction Bander
- Suction tubing 6mm ID x 9mm OD x 2.0M Length, fitted with universal female connectors
- Single use, clean packed, latex free
- UOM: Each, 8/Dispenser Box



Order Code

P305TGR

HAEMOBAND PLUS AND PROCTOSCOPE KIT

- Haemoband Plus Pre-Loaded Suction Bander
- SOLO™ Large Procedural Proctoscope with LED light piece
- Single use, clean packed, latex free
- UOM: Each, 8/Dispenser Box



Order Code

P307TP

Haemorrhoid Banding

Pre-Loaded for fast efficient banding

HAEMOBAND PLUS, SUCTION TUBING AND PROCTOSCOPE KIT



- Haemoband Plus Pre-Loaded Suction Bander
- SOLO™ Large Procedural Proctoscope with LED light piece
- Suction tubing 6mm ID x 9mm OD x 2.0M Length, fitted with universal female connectors
- Single use, clean packed, latex free
- UOM: Each, 8/Dispenser Box

Order Code

P307TPGR

ASKIR 30 SUCTION ASPIRATOR



- Airflow 40L/min, 600mm/Hg
- Continuously variable suction control
- Recommended suction strength for Haemorrhoid banding
- Autoclavable 1L canister and silicone tubing
- Anti-bacterial filter
- Quiet operation
- Made in Italy, 2 year warranty; UOM: Each

Order Code

PA30

CA-MI SUCTION FILTER



- Suction Filter for PA30 Askir 30 Ca-mi Suction Aspirator
- Latex free, clean packed
- UOM: 20/Dispenser Box

Order Code

PSP0046

Self Lighting Sigmoidoscope System


6 integrated LED's
provide brightest
illumination of any
sigmoidoscope

Integrated lighthead
with magnifying lens

Exceptional distal
tip illumination for
optimum visibility

Removable powercell for
eco-friendly disposal

Double bagged disposable
bellows protects from over
insufflation of the bowel



Sigmoidoscope Systems

The widest range of Sigmoidoscopy kits available

ROCKET SELF LIGHTING SIGMOIDOSCOPE SYSTEM



- 6 powerful LED lights provide excellent distal illumination
- Brightest illumination of any Sigmoidoscope available, for optimum visibility
- Removable power cell, ON/OFF switch
- Single use, latex free, clean easipeel packed
- UOM: Each, 4/Dispenser Box

Order Code

P2211SS

SOLO™ SIGMOIDOSCOPE AND BELLOW SET



- Crystal clear plastic guarantees excellent light transmission
- Connects to all common light source systems
- Does not require the Insufflation Filter
- Single use, latex free, clean easipeel packed
- *Lighthead not included
- UOM: Each, 20/Dispenser Box

Order Code

P2102BX

SOLO™ SIGMOIDOSCOPE



- Crystal clear plastic guarantees excellent light transmission
- Compatible with standard reusable Light Heads and light source systems
- Single use, latex free, clean packed
- UOM: Each, 25/Dispenser Box

Order Code

P211SG

Sigmoidoscope Systems

The widest range of Sigmoidoscopy kits available

SOLO™ SIGMOIDOSCOPE FILTER

- For use with reusable Bellows
- Anti-bacterial, Luer lock connection
- Single use, latex free
- UOM: Each, 25/Dispenser Box

Order Code

P210105



REUSABLE BELLOWS

- Surgical grade silicone; Latex free, autoclavable
- UOM: Each

Order Code

P300020



SIGMOIDOSCOPE LIGHT HEAD


- To connect SOLO™ Sigmoidoscope to light source units
- Single piece design with hinged cap
- Stainless steel construction
- Latex free, reusable, autoclavable
- UOM: Each

Order Code

PBC702



Cordless Illuminator by Hallmark Surgical



Lighting Accessories

The exact lighting system to meet your need

SOLO™ SIGMOIDOSCOPE LIGHT HEAD

- Resposable – can be soaked, wiped down, sanitised but not autoclaved, Suitable for 20 uses
- Fits to SOLO™ Sigmoidoscope P211SG and P2102BX Kit
- Connects to cable illuminator
- Non-sterile, latex free
- UOM: Each

Order Code

P21001LH



LED CORDLESS ILLUMINATOR

- For SOLO™ Proctoscopes, Anoscopes and Sigmoidoscopes
- Long lasting bright white LED bulb (50,000 hours)
- Stainless Steel construction
- Autoclavable
- 5 year warranty, made in Canada
- UOM: Each

Order Code

P222LD



CABLE ILLUMINATOR

- 6 Volts x 6 Watts mains power supply with bright, white halogen bulb for excellent colour rendition and illumination
- Ideal for GP and general examination use
- Connects to SOLO™ Proctoscopes, Anoscopes and Sigmoidoscopes via lighthouse or adaptor
- UOM: Each

Order Code

P302900



Lighting Accessories

The exact lighting system to meet your need

SOLO™ PROCTOSCOPE ADAPTOR



- Resposable – can be soaked, wiped down, sanitised but not autoclaved
- Non-sterile, latex free
- Connects SOLO™ Proctoscopes to Cable Illuminator
- UOM: Each

Order Code	Product
P2001AD	Straight Adaptor

SOLO™ PROCTOSCOPE ADAPTOR CURVED



- Resposable – can be soaked, wiped down, sanitised but not autoclaved
- Non-sterile, latex free
- Connects SOLO™ Proctoscopes to Cable Illuminator and LED Cordless Illuminator
- UOM: Each

Order Code	Product
P333CV	Curved Proctoscope Adaptor

PROCTOSCOPE ADAPTOR - STAINLESS STEEL



- Connects SOLO™ Proctoscopes to Cold Light source units
- Stainless steel construction
- Latex free, reusable, autoclavable
- UOM: Each

Order Code	Product
PGA302	Straight adaptor, screw thread fitting

Lighting Accessories

The exact lighting system to meet your need

PROCTOSCOPE ADAPTOR - STAINLESS STEEL

- Connects SOLO™ Proctoscopes to LED Cordless Illuminator & Cable Illuminator
- Stainless steel construction
- Latex free, reusable, autoclavable
- UOM: Each

Order Code	Product
PGA302WA	Straight adaptor, square fitting



SOLO™ ANOSCOPE LIGHT SOURCE ADAPTOR

- Connects SOLO™ Anoscopes to LED Cordless Illuminator and Cable Illuminator
- Stainless Steel & Fibre Optic construction
- Reusable, autoclavable
- UOM:Each

Order Code	Product
PGA303S	Short Handle
PGA305F	Long Handle



SOLO™ LED ILLUMINATOR FOR PROCTOSCOPES

- Suits SOLO™ Proctoscope range
- High light output, latest generation LED
- Pull tab for light switch
- Single use, latex free, clean packed
- UOM: Each, 30/Dispenser Box

Order Code
P2002SL



SOLO™ Polyp Trap

In-flow Polyp Trap which enables uninterrupted suctioning during the procedure.

Clever magnifying glass and transparent construction allows for quick visual verification of specimen

Silicone polyp trap baskets with perforations to separate endoscopic fluid from polyps

Eco-friendly packaging. Packed in recyclable cardboard only

Measurement guide on the polyp trap baskets to assist with polyp sizing and recording.



Endoscopy Accessories

Clever, patented endoscopy accessories

SOLO™ POLYP TRAP

- Enables uninterrupted suctioning
- Separates Polyps from endoscopic effluent
- Magnifying lens and transparent construction allows quick visual verification of specimen
- Supplied with 2 removable specimen baskets for fast and uninterrupted procedure
- Scale on inside of specimen strainers shows size of specimen collected
- Single use, latex free, clean packed
- UOM: Each, 24/Box

Order Code

P312PT



SOLO™ BIOPSY VALVE

- Patented functional design
- Compatible with Olympus and Fujinon endoscopes
- Reduces cross contamination
- Leakage and splash back minimization - patented design
- Punched lid prevents splash back and ensures easy entry for accessories
- Single use, latex free, clean or sterile easipeel packed
- UOM: 100/Pack

Order Code

Product

P321VC

Clean-packed

P322BJ

Sterile-packed



Bowel Washout Systems

Innovative, enclosed bowel washout systems for odour free and fast procedures

SOLO™ BOWELBATHE™ ON-TABLE BOWEL WASHOUT KIT



- All-in-one ready to use kit for intra-operative colonic lavage
- Fast, easy, odourless, sealed collection system
- 10-15 minutes' average washout time
- No second incision required at appendix
- No leakage or soiling
- Enables on-table colonoscopy following lavage
- Single use, latex free, sterile easipeel packed
- UOM: Each

Order Code	Product
P501RW	35mm hose
P502WR	46mm hose

SOLO™ BOWELBATHE™ RECTAL WASHOUT KIT



- All-in-one ready to use kit, for washing out rectum and distal colon of faecal matter and exfoliated cancer cells
- Soft atraumatic proctoscope tip instead of obturator
- Fast, easy, odourless and leak-free sealed collection system
- Large bore waste tubing
- 3ltr sealed waste container
- Single use, latex free, sterile easipeel packed
- UOM: Each

Order Code	Product
P505PR	Rectal Washout Kit

Holloware



HOLLOWARE

Disposable Holloware	84-85
Reusable Quivers	86
Autoclaveable Kidney Dishes	87
Autoclaveable Procedure and Dressing Trays	88
Autoclaveable Perforated Trays	89-90
Autoclaveable Medical Bowls	91-92

Disposable Holloware

SOLO™ PULP TRAY, 750ML



- Made from sugar-cane pulp
- Natural colour meaning no toxic dyes or chemicals used in production
- 29cm x 16cm x 4cm deep
- Excellent for general use in healthcare facilities
- Non-sterile
- UOM 500/carton

Order Code

Q120BP

SOLO™ PROCEDURE TRAY, YELLOW, 500ML



- 100% pharmaceutical grade polypropylene
- Graduated with millilitres
- Strong, sturdy design
- 27cm x 13.5cm x 5cm deep
- Sterile packed
- UOM 30/carton

Order Code

Q200PT

Disposable Holloware

SOLO™ BOWL SET #1 - 6 PIECE

- 100% pharmaceutical grade polypropylene
- Graduated with millilitres
- Strong, sturdy design
- Contains: 1 x 1000ml Blue Kidney Dish, 2 x 700ml Yellow Kidney Dish, 3 x 250ml Blue Bowls
- Sterile packed
- UOM 4/inner, 24/outer



Order Code
Q900SC

SOLO™ DIATHERMY QUIVER, 20CM

- 100% pharmaceutical grade polypropylene
- Ring on rim of quiver for securing to trolleys or tables
- Strong, sturdy design
- 20cm long, 7cm diameter
- UOM 50/carton or 150/carton



Order Code	Product
Q500ST	Sterile
Q500VR	Non-Sterile

Reusable Holloware


ELECTRISURE™ DIATHERMY QUIVERS



- Screw off base for easy cleaning
- Stainless Steel holding ring
- High grade medical plastic - will not split, crack or powderise
- Belled top design for easy instrument placement
- Reusable, autoclavable at 134°C
- UOM: Each


Order Code	Product
EDQS15	150mm
EDQS20	200mm
EDQS25	250mm
EDQS30	300m
EDQS35	350mm
EDQS40	400mm
EDQSC	Quiver end cap only
EDQSR	Quiver ring only

Autoclaveable Kidney Dishes




KIDNEY DISH - 160MM

Order Code	Product
QKD16GR	Green
QKD16BL	Blue
QKD16YE	Yellow




KIDNEY DISH - 220MM

Order Code	Product
QKD22GR	Green
QKD22BL	Blue
QKD22YE	Yellow



KIDNEY DISH - 255MM

Order Code	Product
QKD25GR	Green
QKD25BL	Blue
QKD25YE	Yellow



KIDNEY DISH - 300MM

Order Code	Product
QKD30GR	Green
QKD30BL	Blue
QKD30YE	Yellow

Autoclaveable Procedure and Dressing Trays



INJECTION TRAY - 200 X 75 X 30MM

Order Code	Product
QNJ200GR	Green
QNJ200YE	Yellow



INJECTION TRAY - 270 X 100 X 40MM

Order Code	Product
QNJ270GR	Green
QNJ270YE	Yellow



3 COMPARTMENT TRAY - 160 X 140MM

Order Code	Product
QRD1614GR	Green



3 COMPARTMENT TRAY - 205 X 165MM

Order Code	Product
QRD2016GR	Green
QRD2016BL	Blue



3 COMPARTMENT TRAY - 270 X 150MM

Order Code	Product
QRD2715GR	Green



4 COMPARTMENT TRAY - 270 X 205MM

Order Code	Product
QRD2720G	Green

Autoclaveable Perforated Trays



PERFORATED TRAYS - 180 X 150 X 30MM

Order Code	Product
QRP1815GR	Green




PERFORATED TRAYS - 200 X 100 X 40MM

Order Code	Product
QRP2010GR	Green



PERFORATED TRAYS - 270 X 100 X 40MM

Order Code	Product
QRP2710GR	Green
QRP2710YE	Yellow



PERFORATED TRAYS - 270 X 150 X 30MM

Order Code	Product
QRP2715GR	Green
QRP2715BL	Blue



PERFORATED TRAYS - 300 X 270 X 60MM

Order Code	Product
QRP3027GR	Green
QRP3027BL	Blue



PERFORATED TRAYS - 400 X 135 X 60MM

Order Code	Product
QRP4013GR	Green

Autoclaveable Perforated Trays



PERFORATED TRAYS - 400 X 270 X 60MM

Order Code	Product
QRP4027GR	Green
QRP4027BL	Blue



PERFORATED TRAYS - 540 X 300 X 60MM

Order Code	Product
QRP5430GR	Green
QRP5430BL	Blue



PERFORATED TRAYS - 660 X 300 X 60MM

Order Code	Product
QRP6630GR	Green

Autoclaveable Medical Bowls




BOWL, 60MM

Order Code	Product
QWL60GR	Green




BOWL, 80MM

Order Code	Product
QWL80GR	Green




BOWL, 100MM

Order Code	Product
QWL10GR	Green
QWL10BL	Blue




BOWL, 140MM

Order Code	Product
QWL14GR	Green
QWL14BL	Blue



BOWL, 185MM


Order Code	Product
QWL18GR	Green



BOWL, 210MM


Order Code	Product
QWL21GR	Green

Autoclaveable Medical Bowls



BOWL, 240MM

Order Code	Product
QWL24GR	Green
QWL24BL	Blue



BOWL, 305MM

Order Code	Product
QWL30GR	Green
QWL30BL	Blue



BOWL, 345MM

Order Code	Product
QWL34GR	Green

SOLO™ Single Use Instruments



SOLO™ Single Use Instruments



BLADE HANDLE #3

Order Code	Product
RSU120103	12cm



IRIS SCISSORS

Order Code	Product
RSU130014	11cm Straight
RSU130015	11cm Curved



SURGICAL SCISSORS

Order Code	Product
RSU130302	13cm Sharp/Blunt, Straight
RSU130332	13cm Sharp/Sharp, Straight



MAYO SCISSORS

Order Code	Product
RSU130400	14cm Straight



MAYO STILLE SCISSORS

Order Code	Product
RSU130434	17cm Straight



METZENBAUM BABY SCISSORS

Order Code	Product
RSU130598	14.5cm Straight

SOLO™ Single Use Instruments



METZENBAUM SCISSORS

Order Code	Product
RSU130606	23cm Straight
RSU130607	23cm Curved



SIMS SCISSORS

Order Code	Product
RSU132101	20cm Curved
RSU132102	23cm Straight
RSU132100	20cm Straight



USA UMBILICAL SCISSORS

Order Code	Product
RSU132220	10.5cm



IRIS DRESSING FORCEPS

Order Code	Product
RSU140010	10cm Straight



ADSON DRESSING FORCEPS

Order Code	Product
RSU140030	12cm Micro
RSU140040	12cm Standard

SOLO™ Single Use Instruments



MCINDOE DRESSING FORCEPS

Order Code	Product
RSU140070	16cm



DRESSING FORCEPS

Order Code	Product
RSU140103	14.5cm



FORCEPS ADSON TISSUE

Order Code	Product
RSU140530	12cm, 1x2 Teeth Micro
RSU140540	12cm, 1x2 Teeth Standard



GILLIES TISSUE FORCEPS

Order Code	Product
RSU140570	15cm 1x2 Teeth



TISSUE FORCEPS

Order Code	Product
RSU140603	14.5cm 1x2 Teeth



FEILCHENFIELD SPLINTER FORCEPS

Order Code	Product
RSU141534	11cm

SOLO™ Single Use Instruments



JEWELLER FORCEP #1

Order Code	Product
RSU141600	11.5cm



MOSQUITO MICRO ARTERY FORCEPS

Order Code	Product
RSU150101	10cm, Curved



MOSQUITO ARTERY FORCEPS

Order Code	Product
RSU150120	12.5cm Straight
RSU150121	12.5cm Curved



KELLY ARTERY FORCEP

Order Code	Product
RSU150200	14cm Straight



SPENCER WELLS ARTERY FORCEPS

Order Code	Product
RSU150412	14cm Straight
RSU150413	14cm Curved
RSU150416	18cm Straight
RSU150418	20cm Straight



ROCHESTER PEAN ARTERY FORCEPS

Order Code	Product
RSU150460	14cm Straight
RSU150462	16cm Straight
RSU150464	18cm Straight

SOLO™ Single Use Instruments



BALL & SOCKET TOWEL CLAMP

Order Code	Product
RSU160135	14cm



PRESBYTERIAN TUBING CLAMP

Order Code	Product
RSU160202	20cm



FOERSTER SPONGE FORCEPS

Order Code	Product
RSU160416	25cm Straight



BONNEY POLYP FORCEPS

Order Code	Product
RSU160414	24cm, Straight



POZZI TENACULUM FORCEPS

Order Code	Product
RSU550944	25cm Straight



GILLIES SKIN HOOK

Order Code	Product
RSU170041	#2, 18cm

SOLO™ Single Use Instruments



TWO PRONG RETRACTOR

Order Code	Product
RSU170412	16cm, Blunt



SENN MILLER DOUBLE ENDED RETRACTOR

Order Code	Product
RSU170580	15cm



DESMARRES RETRACTOR

Order Code	Product
RSU300522	13cm, 12mm



WEITLANER RETRACTOR

Order Code	Product
RSU175303	3x4 teeth Blunt, 13cm



DERF NEEDLE HOLDER

Order Code	Product
RSU181010	12.5cm



CRILE WOOD NEEDLE HOLDER

Order Code	Product
RSU181070	15cm

SOLO™ Single Use Instruments



MAYO HEGAR NEEDLE HOLDER

Order Code	Product
RSU181078	15cm
RSU181080	16cm
RSU181082	18cm
RSU181084	20cm



DOUBLE ENDED PROBE

Order Code	Product
RSU190122	14cm x 2.0mm



TROELTSCH ENT FORCEP

Order Code	Product
RSU320642	12cm



HARTMANN EAR FORCEP

Order Code	Product
RSU320700	8cm



THUDICHUM NASAL SPECULUM

Order Code	Product
RSU330083	Small sterile
RSU330083PK10	Small, Pack of 10, Non-Sterile
RSU330085	Medium sterile
RSU330085PK10	Medium, Pack of 10, Non-Sterile
RSU330087	Large sterile



BARRAQUER EYE SPECULUM

Order Code	Product
RSU300004	4cm

SOLO™ Single Use Instruments



BONE CURETTE

Order Code	Product
RSU160135	Double-Ended 14cm



MAGILL CATHETER FORCEPS

Order Code	Product
RSU110300	25cm Adult
RSU110304	17cm Adult

MATERNITY PACK

Order Code	Product
RSU903MT	Contains: <ul style="list-style-type: none"> • 1 x Surgical Scissors, 13cm, Sharp/Blunt • 1 x Gillies Tissue Forceps, 15cm • 1 x Rochester Pean Artery Forcep, 22cm Straight • 1 x Mayo Hegar Needleholder, 22cm



DELIVERY PACK

Order Code	Product
RSU902DP	Contains: <ul style="list-style-type: none"> • 2 x Spencer Wells Forcep, 18cm Straight • 1 x Umbilical Scissors, 10.5cm

IUD REMOVAL PACK

Order Code	Product
RSU906ID	Contains: <ul style="list-style-type: none"> • 1 x Pozzi Uterine Tenaculum, 25cm Straight • 1 x Foerster Sponge Forcep, 24cm Straight • 1 x Metzenbaum Scissor, 23cm Straight

SOLO™ Single Use Instruments



SUTURE REMOVAL SET

Order Code	Product
RSU901SF	Contains: <ul style="list-style-type: none"> • 1 x Iris Scissor, 11cm Straight • 1 x Adson Dressing Forcep, 12cm



SUTURE SET

Order Code	Product
RSU904WL	Contains: <ul style="list-style-type: none"> • 1 x Iris Scissor, 11cm Straight • 1 x Adson Tissue Forcep, 12cm, 1x2 Teeth • 1 x Mayo Hegar Needleholder, 14cm



SUTURE REMOVAL SET #2

Order Code	Product
RSU905RP	Contains: <ul style="list-style-type: none"> • 1 x Iris Scissor, 11cm Straight • 1 x Adson Tissue Forcep, 12cm, 1x2 Teeth



FINE SUTURE PACK

Order Code	Product
RSU907XM	Contains: <ul style="list-style-type: none"> • 1 x Iris Scissor, 11cm Straight • 1 x Adson-Micro Tissue Forcep, 12cm, 1x2 Teeth • 1 x Halsey Needleholder, 13cm

Womens Health Devices



WOMENS HEALTH DEVICES

Disposable Speculum	104-107
Lletz Procedure Accessories	108
Single Use Gynaecology Instruments	109-111
Pessaries	112-119
General Consumables	120-121
Uterine Manipulators	122-125

SOLO™ Self-Lighting Speculum

Strong crystal clear medical polymer construction for excellent lighting transmission



Versatile Graves style speculum available in 3 sizes.



Vertical and pivoting movement for optimum site visualisation

Self-contained LED illumination provides excellent visibility



Cordless lighting, no ancillary charger required



Disposable Speculum

Wide range of innovative speculum for every procedure

SOLO™ GRAVES SELF ILLUMINATING SPECULUM

- Inbuilt bright white LED
- Blades both pivot and lift vertically
- Strong, crystal clear medical polymer construction
- Self-contained LED Illumination for cordless lighting
- No battery, no charger, no cable
- Complete freedom of movement!
- Single use, latex free, clean easipeel packed
- UOM Each or 25/dispenser-box

Order Code	Product
W5910LS	Small
W5911LM	Medium
W5914LL	Large



SOLO™ EASYSPEC™ SILENT SPECULUM

- Silent opening & closing mechanism, very strong design
- Fast locking key mechanism, colour coded for size
- Stepless locking positions
- FrictionLESS™ surfaced blades for smooth insertion
- Made in Sweden
- Single use, latex free, clean packed
- UOM Each or 40/dispenser-box


Order Code	Product
W110XS	Extra Small
W111S	Small
W112M	Medium
W113L	Large
W114XL	Extra Large



Smoke-X™ Disposable Diathermy Speculum

Special diathermy plastic
providing safe insulation,
with no chance of failure

Full length smoke tube under
the top blade with multiple
ports for effective extraction
and clear view.



■ Excellent visibility
through wide
opening design

Pearl White colour for
lighting reflection and
accurate colour rendition

Disposable smoke-extracting
speculum, your complete
solution to costly reusable
insulated speculum failures

Disposable Speculum

Wide range of innovative speculum for every procedure

SIMS SPECULUM

- Surgical quality speculum
- 80mm and 90mm speculum lengths
- Strong, crystal clear medical polymer construction
- Single use, latex free, sterile easipeel packed
- UOM Each or 15/dispenser-box

Order Code

W132MS



SMOKE-X™ DIATHERMY SPECULUM - SINGLE USE

- Surgical quality, special diathermic plastic
- Effective smoke extraction with full length multi-port tube
- Pearl white colour for lighting efficacy
- Single use, latex free, sterile easipeel packed
- UOM Each or 20/dispenser-box

Order Code

Product

W153MM	Medium
W154ML	Medium Long
W155BR	Broad



SOLO™ POLYSPEC SPECULUM FOR HYSTEROSCOPY

- The Polyspec is ergonomically designed with eight different opening widths and is often used for larger examinations.
- Stable construction
- Silent locking device
- Large access opening
- Large treatment and examination space
- UOM Each or 10/dispenser-box

Order Code

Product

W162WP	Polyspec
--------	----------



Lletz Procedure Accessories

ELECTRISURE™ DIATHERMY PENCIL



- 3M cable with universal 3 pin generator plug
- Supplied with Teflon coated non-stick blade electrode
- Single use, latex free, sterile easipeel packed
- UOM Each or 40/dispenser-box

Order Code

E102PP

SOLO™ LLETZ LOOP ELECTRODES



- High grade tungsten wire
- 13cm length with universal hex connection
- Single use, latex free, sterile easipeel packed
- UOM Each or 10/dispenser-box

Order Code	Product
E60511010PT	10x10mm
E60511510WA	15x10mm
E60511512EU	15x12mm
E60511515BV	15x15mm
E60512008NR	20x8mm
E60512010BB	20x10mm
E60512012ZB	20x12mm
E60512015BG	20x15mm
E60512020PR	20x20mm

SOLO™ BALL ELECTRODE



- 13cm length x 5mm dia. Teflon coated
- Single use, latex free, sterile packed
- UOM Each or 10/dispenser-box

Order Code

E6045130DAT

Single Use Gynaecology Instruments

Innovative disposable devices perfect for your clinic

IUD THREAD RETRIEVER

- 245mm length
- Well defined hooks for fast effective use
- Single use, latex free, sterile easipeel packed
- UOM Each or 25/dispenser box

Order Code

W2111UD



RAMPLEY FORCEPS

- 250mm length
- Reinforced nylon construction, ratchet locking handles
- Alignment shaft for strong grip
- Single use, latex free, sterile easipeel packed
- UOM Each or 30/dispenser-box

Order Code

W21141RM



UTERINE SOUND

- 290mm length with reinforced ergonomic handle
- Sliding depth indicator, rounded atraumatic tip
- Single use, latex free, sterile easipeel packed
- UOM Each or 25/dispenser-box

Order Code

W2112ES



Single Use Gynaecology Instruments

Innovative disposable devices perfect for your clinic

UTERINE SOUND



- 265mm length, flexible medical HD Polyethylene construction
- Rounded atraumatic tip and notched 1cm depth indicators
- Single use, latex free, sterile easipeel packed
- UOM Each or 25/dispenser-box

Order Code

W2115PN

UTERINE SOUND-DILATOR




- Double-ended combined device designed for gentle dilation of the cervix and depth sounding.
- 277mm length with maximum diameter of 5.5mm
- Strong but flexible medical HD Polyethylene construction
- Smooth, rounded atraumatic tip, 1 cm depth indicator marks.
- Single use, latex free, sterile easipeel packed
- UOM Each or 25/dispenser-box

Order Code

W2118SD

HEGAR CERVICAL DILATORS



- 3 sizes, double-ended
- Reinforced nylon construction
- Central thumb grip
- Single use, latex free, sterile easipeel packed
- UOM Each or 25/dispenser-box

Order Code

Product

W21131DL	3/4mm
W21132GP	5/6mm
W21133SE	7/8mm

Single Use Gynaecology Instruments

Innovative disposable devices perfect for your clinic

ENDOSAMPLER™

- Disposable curettes for sampling the endometrium. The Endosampler II™, has two holes for even easier suction. The unique design of the Endosampler allows for superior accuracy and a more reliable sample than traditional endometrial sample devices.
- Convenient for routine use, suitable for cytologic and histologic diagnosis
- Facilitates both endometrial scraping and aspiration
- Avoids endometrial sample contamination
- Semi-rigid 3mm curette minimizes possibility of perforation
- Curette conveniently bent with sharp slot on side assures better sampling from all areas of the endometrial cavity
- Backflow of collected sample is prevented by lock-spring mechanism on syringe
- Sterile easipeel packed, single use, latex free
- UOM Each or 25/dispenser-box



Order Code	Product
W22720ES	Single hole
W22722ME	Twin Hole

AMNIOHOOK AMNIOTIC MEMBRANE PERFORATOR

- Safely ruptures amniotic membrane with rounded tip hook design, minimising foetal injury risk
- Avoids contact with maternal tissue
- Single use, latex free, sterile easipeel packed



Order Code	Product
W800MP	Straight
W802CV	Curved

SOLO™ Pessaries

Highest grade surgical
silicone construction

Flexible supple pessaries
providing the perfect balance
between clinical effectiveness
and patient comfort



Full range of pessaries available
from rings to specialist styles such
as Gellhorns, Cubes and Hodge

22 different styles and a
complete range of sizes
available

Pessaries

Surgical grade silicone for New Zealand's most supple pessary range

RING WITHOUT SUPPORT

Order Code	Size
W4001	51mm
W4002	57mm
W4003	64mm
W4004	70mm
W4005	76mm
W4006	83mm
W4007	89mm
W4008	95mm
W4009	102mm



RING WITH SUPPORT

Order Code	Size
W4011	51mm
W4012	57mm
W4013	64mm
W4014	70mm
W4015	76mm
W4016	83mm
W4017	89mm
W4018	95mm
W4019	102mm



RING WITH KNOB WITHOUT SUPPORT

Order Code	Size
W4021	51mm
W4022	57mm
W4023	64mm
W4024	70mm
W4025	76mm
W4026	83mm
W4027	89mm



Pessaries

Surgical grade silicone for New Zealand's most supple pessary range

RING WITH KNOB WITH SUPPORT



Order Code	Size
W4031	51mm
W4032	57mm
W4033	64mm
W4034	70mm
W4035	76mm
W4036	83mm
W4037	89mm

DONUT



Order Code	Size
W4040	51mm
W4041	57mm
W4042	64mm
W4043	70mm
W4044	76mm
W4045	83mm
W4046	89mm
W4047	95mm

OVAL WITHOUT SUPPORT



Order Code	Size
W4051	51mm
W4052	57mm
W4053	64mm
W4054	70mm
W4055	76mm
W4056	83mm
W4057	89mm
W4058	95mm
W4059	102mm

Pessaries

Surgical grade silicone for New Zealand's most supple pessary range

OVAL WITH SUPPORT

Order Code	Size
W4061	51mm
W4062	57mm
W4063	64mm
W4064	70mm
W4065	76mm
W4066	83mm
W4067	89mm
W4068	95mm
W4069	102mm



SHAATZ

Order Code	Size
W4070	38mm
W4071	44mm
W4072	51mm
W4073	57mm
W4074	64mm
W4075	70mm
W4076	76mm
W4077	83mm
W4078	89mm



GELLHORN WITH DRAIN, LONG STEM

Order Code	Size
W4080	38mm
W4081	44mm
W4082	51mm
W4083	57mm
W4084	64mm
W4085	70mm
W4086	76mm
W4087	83mm
W4088	89mm



Pessaries

Surgical grade silicone for New Zealand's most supple pessary range

GELLHORN WITH DRAIN, SHORT STEM



Order Code	Size
W4091	44mm
W4092	51mm
W4093	57mm
W4094	64mm
W4095	70mm
W4096	76mm
W4097	83mm

CUBE WITHOUT DRAIN



Order Code	Size
W4100	25mm
W4101	29mm
W4102	33mm
W4103	37mm
W4104	41mm
W4105	45mm
W4106	50mm
W4107	56mm
W4108	63mm
W4109	70mm
W4110	75mm

CUBE WITH DRAIN



Order Code	Size
W4120	25mm
W4121	29mm
W4122	33mm
W4123	37mm
W4124	41mm
W4125	45mm
W4126	50mm
W4127	56mm
W4128	63mm
W4129	70mm
W4130	75mm

Pessaries

Surgical grade silicone for New Zealand's most supple pessary range

CERCLAGE WITHOUT DRAIN

Order Code	Size
W4141	65x32x21mm
W4142	70x32x21mm
W4143	65x35x21mm
W4144	70x35x21mm
W4145	65x35x30mm
W4146	70x35x30mm



CERCLAGE WITH DRAIN

Order Code	Size
W4151	65x32x21mm
W4152	70x32x21mm
W4153	65x35x21mm
W4154	70x35x21mm
W4155	65x35x30mm
W4156	70x35x30mm



CUP WITH SUPPORT

Order Code	Size
W4160	50mm
W4161	55mm
W4162	60mm
W4163	65mm
W4164	70mm
W4165	75mm
W4166	80mm
W4167	85mm



CUP WITHOUT SUPPORT

Order Code	Size
W4170	50mm
W4171	55mm
W4172	60mm
W4173	65mm
W4174	70mm
W4175	75mm
W4176	80mm
W4177	85mm



Pessaries

Surgical grade silicone for New Zealand's most supple pessary range

DISH WITH SUPPORT



Order Code	Size
W4180	50mm
W4181	55mm
W4182	60mm
W4183	65mm
W4184	70mm
W4185	75mm
W4186	80mm
W4187	85mm

DISH WITHOUT SUPPORT



Order Code	Size
W4190	50mm
W4191	55mm
W4192	60mm
W4193	65mm
W4194	70mm
W4195	75mm
W4196	80mm
W4197	85mm

MARLAND WITH SUPPORT



Order Code	Size
W4202	57mm
W4203	64mm
W4204	70mm
W4205	76mm
W4206	83mm
W4207	89mm
W4208	95mm

Pessaries

Surgical grade silicone for New Zealand's most supple pessary range

MARLAND WITHOUT SUPPORT

Order Code	Size
W4212	57mm
W4213	64mm
W4214	70mm
W4215	76mm
W4216	83mm
W4217	89mm
W4218	95mm




HODGE WITH SUPPORT

Order Code	Size
W4220	65mm
W4221	70mm
W4222	75mm
W4223	80mm
W4224	85mm
W4225	90mm
W4226	95mm
W4227	100mm
W4228	105mm
W4229	110mm



HODGE WITHOUT SUPPORT

Order Code	Size
W4230	65mm
W4231	70mm
W4232	75mm
W4233	80mm
W4234	85mm
W4235	90mm
W4236	95mm
W4237	100mm
W4238	105mm
W4239	110mm



General Consumables

Everyday consumables for your clinic

BASTION LATEX MEDICAL EXAMINATION GLOVES, POWDER FREE



- Natural biodegradable rubber
- Greatest material durability and flexibility
- Superior tactile sensitivity, comfort, and fit
- UOM: 100/Dispenser Box

Order Code	Product
G100161	Small
G100162	Medium
G100163	Large
G100164	Extra Large

SURGI-gel PLUS LUBRICATING GEL SACHETS



- Sterile, latex free 2.7g Sachets
- Water soluble, greaseless, transparent, non-staining, non-irritating
- High viscosity for excellent adherence
- UOM: 144/Dispenser Box

Order Code
G100120

General Consumables

Everyday consumables for your clinic

AGEL LUBRICATING GEL - 500ML

- Water soluble, greaseless, transparent, non-staining, non-irritating
- High viscosity for excellent adherence
- Latex free, clean packed
- UOM: Each

Order Code

G100100



SOLO™ EXAMINATION UNDERPADS

- Soft touch super wick absorbent
- Latex free, 600 X 400mm, 3 ply
- Individually folded, clean packed ready for use
- UOM 250/dispenser-box

Order Code Product

W301UP 60 x 40cm



Disposable Spackmann Cannula

Evenly placed perforations in the tip for equal and accurate dye distribution

Save reprocessing time and cost while eliminating all cross contamination and infection risk


Unique, single-handed operating, self-locking clamp designed for fast tenaculum adjustment

Fixed or movable acorn designs available for a perfect seal on the cervix

Guaranteed clean, guaranteed sterile

Made in the UK from premium materials

Strong, rigid stainless steel shaft for effective manipulation as required



Uterine Manipulators

Exceptional design and quality for improved manipulation and visualisation

SUREHOLD™ SPACKMANN UTERINE MANIPULATOR

- Risk of cross contamination eliminated with single use device
- For hydrotubation dye tests and uterine manipulation during laparoscopic procedures
- Unique self-locking clamp enables single handed use
- Option of fixed or movable acorn for perfect seal on cervix
- Evenly placed holes at the tip for equal dye distribution
- Non-blocking design
- 240mm length
- Single use, latex free, sterile easipeel packed
- UOM Each or 12/dispenser-box



Order Code	Product
W2121NS	Fixed Acorn
W2122MA	Moveable Acorn

The Clearview Uterine Manipulator

Single handed with the self holding rotating control handle allows the surgeon to maneuver the device independent of an assistant

Greatest range of motion of all uterine manipulators - 90° Retroversion and 120° Anteversion - and lateral movement for clear visibility

The Clearview pivots at the cervical os allowing the patient to remain in the supine position.



Single use - completely disposable

One kit contains everything you require including a sound-dilator, 12ml syringe for the tip balloon and a dye injection port

Uterine Manipulators

Exceptional design and quality for improved manipulation and visualisation

CLEARVIEW UTERINE MANIPULATOR

- The ClearView Manipulator gives the greatest range of movement of any manipulator. The rotating control handle facilitates one-handed control, enabling the surgeon to manipulate the uterus without breaking the sterile barrier.
- Available in 4 tip lengths
- Kit is complete – no additional components required
- Dye injector – luer fitting & balloon inflation valve for effective chromotubation
- Single use, latex free, sterile easipeel packed
- UOM Each or 5/dispenser-box



Order Code	Product
WUM201	6/7cm Tip
WUM202	8/9cm Tip

CLEARVIEW TOTAL UTERINE MANIPULATOR – FOR THE TLH PROCEDURE

- This model has the same great features as the ClearView above with the addition of a pre-attached occluder balloon and three sizes of colpotomy cups. The colpotomy cups are designed for improved visibility during insertion and are coloured for easy identification when making the colpotomy incision. The occluder balloon serves as a seal, preventing the loss of the pneumoperitoneum once the colpotomy incision is made.
- Includes 3 x colpotomy cups & occluder balloon for TLH, depth spacer turret, uterine sound, 60ml syringe and 1 x 10ml syringe
- Available in 4 tip lengths
- Single use, latex free, sterile easipeel packed
- UOM Each or 5/dispenser-box



Order Code	Product
W9207CC	6/7cm tip
W9209SD	8/9cm tip

About Us



Hallmark Surgical has deliberately set the bar high. The name stands for excellence and commitment to, arguably, the most responsible, and most respected profession known to man. To fall short, when a life, (or lives), are at stake – on your part, or ours – would be unthinkable. But, more than that, Hallmark Surgical backs the skill and experience of people at all levels of health care. All life is both precious and perilous. And the instruments, devices and clothing we provide are devoted to preserving that life.

Because every day, like you, we both protect and serve.

HallmarkSurgical™
When you care.

Proving
indispensable
to those
who care.