

## PeliSpec Smoke-X Single Use Diathermy Speculum

Traditional reusable insulated speculum have long proven to be a high cost, high maintenance instrument. The smoke tube, curved surfaces and small pivoting or threaded components are slow and difficult to clean, the smoke tube cannot be visually checked, and very often the insulation starts to lift or bubble after a relatively small number of processes.

The PeliSpec Smoke-X Speculum is the ideal solution. This user friendly and economical speculum range is intrinsically safe due to its special plastic construction and very effective in removing harmful surgical plume.

***A new, sterile instrument for every procedure - eliminates all cross-infection risk - saves time  
- saves cost - a consistent outcome for every procedure!***



1. Special diathermic plastic construction provides intrinsic electrical non-conductivity
2. Guaranteed sterile instrument for every procedure, eliminating any infection risk
3. Strong, wide opening blades for excellent visibility
4. The full length smoke tube with multiple ports, guarantees excellent smoke extraction.
5. Textured handles for positive grip, with fast, smooth running nut/thread mechanism
6. Three sizes are available to suit all patients: Medium, Medium-Long and Broad.
7. The high maintenance or replacement cost of reusable insulated Speculum means the PeliSpec Speculum is substantially lower cost per procedure.
8. Saves the CSD Dept significant time in difficult re-processing procedures

### Size options and order codes:

W153MM: Medium

W154ML: Medium-Long

W155BR: Broad size

***Request samples or speak with a Product Specialist!***

**AU Ph** 62 2 8355 7189

**AU Free** 1800 99 82 13

**NZ Ph** 64 9 302 0755

**NZ Free** 0800 67 37 28

**E** operate@hallmarksurgical.com

Level 6, 10 Herb Elliott Ave  
Sydney Olympic Park, NSW 2127, Australia

PO Box 8087, Symonds St  
Auckland 1150, New Zealand

**W** hallmarksurgical.com

**HallmarkSurgical™**  
When you care.