



UROLOGY



UROLOGY



ABOUT US



SINCE 20 YEARS



EXPORT IN
+70 COUNTRIES



24H/SERVICE



+ 10.000
PRODUCTS



LOCATED IN THE
WORLD CENTER
FOR MEDICAL
TECHNOLOGY

CERTIFICATED

STEMA Medizintechnik GmbH is certified in accordance with the highest standards:

- EN ISO 13485:2016 + AC:2018 + A11:2021- ISO 13485:2016
- EU 2017/745 (**MDR**)
- Council Directive 93/42/EEC Annex II excluding Section 4
- Council Directive 93/42/EEC Annex V excluding Section 3
- FDA

These standards create an essential base for our company and our partners. They also symbolize the commitment and strive to keep up the quality of our products.



PRODUCT RANGE

The current product range of STEMA Medizintechnik GmbH includes more than 10.000 different items.

We cover the following areas in surgery with high quality products:

ARTHROSCOPY
BRONCHOSCOPY
CARDIAC MIS
CONTAINER-SYSTEMS
GASTROENTEROLOGY
GYNAECOLOGY
HF-SURGERY
LAMINECTOMY
LAPAROSCOPY
MICRO MOTOR SYSTEMS
PROCTOSCOPY
SPINE-SURGERY
THORACOSCOPY
TITANIUM PLATING SYSTEM
UNITS
UROLOGY



CONTENT

UROLOGY

PRODUCT

PAGE

1. CYSTO-URETHROSCOPY

1.1 Cysto-Urethroscopy	8
1.2 Slim Cysto-Urethroscopy	25
1.3 Laser Cysto-Urethroscopy	28
1.4 Mini HD Compact Cystoscope	30
1.5 Flexible Cystoscope	31
1.6 Cystoscopy Set	33

2. RESECTOSCOPY

2.1 Monopolar Resectoscopy	37
2.2 Saline Resectoscopy	59
2.3 Laser Resectoscopy	76
2.4 Accessories	87
2.5 Resectoscopy Sets	91

3. URETHROTOMY

3.1 Urethrotome	97
3.2 Optical Urethrotomy	98
3.3 Urethrotomy Set	106

4. NEPHROSCOPY

4.1 HD Nephroscopy	109
4.2 Mini Nephroscopy	111
4.3 Standard Nephroscopy	114
4.4 Accessories	118

5. URETERO-RENSCOPY

5.1 Uretero-Renoscopes semi-rigid	123
5.2 Flexible Uretero-Renoscope	142

CONTENT

UROLOGY

PRODUCT

PAGE

6. LITHOTRIPSY

6.1 Lithotripsy	145
-------------------	-----

7. MORCELLATION

7.1 Morcellation 26 Fr.	150
7.2 Morcellation 22 Fr.	152

8. PEDIATRIC UROLOGY

8.1 Pediatric Cysto-Urethroscopy	155
8.2 Pediatric Resectoscopy	161
8.3 Pediatric Urethrotomy	169

9. ACCESSORIES

9.1 Instruments	176
9.2 Sterilization	181
9.3 Cleaning and Maintenance	183

10. INDEX

10.1 Index	186
--------------	-----

1. CYSTO-URETHROSCOPY

1.1 | CYSTO-URETHROSCOPY

CYSTOSCOPES (4.0MM)

SHEATHS

OBTURATORS

ENDOSCOPE BRIDGES

FORCEPS, SCISSORS & INSTRUMENTS

STONE EXTRACTORS

DIATHERMY LOOPS

MONOPOLAR CONNECTING CABLES

1.2 | SLIM CYSTO-URETHROSCOPY

CYSTOSCOPES (2.9MM)

SHEATHS

OBTURATORS

ENDOSCOPE BRIDGES

1.3 | LASER CYSTO-URETHROSCOPY

CYSTOSCOPES (4.0MM)

SHEATHS

OBTURATORS

1.4 | MINI HD COMPACT CYSTOSCOPES

1.5 | FLEXIBLE CYSTOSCOPES


1.6 | CYSTOSCOPY SET

1. CYSTO-URETHROSCOPY


4MM

2.9MM

1 CYSTOSCOPES



2 SHEATHS & OBTURATORS




3 ACCESSORIES



1.1 | CYSTO-URETHROSCOPY

CYSTOSCOPES AUTOCLAVABLE Ø 4.0MM



Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 0°	Cystoscope autoclavable	Ø 4mm	302mm	AA-002-001
	Cystoscope autoclavable, wide angle	Ø 4mm	302mm	AA-002-001-W
	HD Cystoscope autoclavable	Ø 4mm	302mm	AA-002-011
	HD Cystoscope autoclavable, wide angle	Ø 4mm	302mm	AA-002-011-W
	Ultra HD (4K) Cystoscope autoclavable	Ø 4mm	302mm	AA-002-021

Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 12°	Cystoscope autoclavable	Ø 4mm	302mm	AA-002-002
	HD Cystoscope autoclavable	Ø 4mm	302mm	AA-002-012
	Ultra HD (4K) Cystoscope autoclavable	Ø 4mm	302mm	AA-002-022

Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 30°	Cystoscope autoclavable	Ø 4mm	302mm	AA-002-003
	Cystoscope autoclavable, wide angle	Ø 4mm	302mm	AA-002-003-W
	HD Cystoscope autoclavable	Ø 4mm	302mm	AA-002-013
	HD Cystoscope autoclavable, wide angle	Ø 4mm	302mm	AA-002-013-W
	Ultra HD (4K) Cystoscope autoclavable	Ø 4mm	302mm	AA-002-023

Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 70°	Cystoscope autoclavable	Ø 4mm	302mm	AA-002-004
	Cystoscope autoclavable, wide angle	Ø 4mm	302mm	AA-002-004-W
	HD Cystoscope autoclavable	Ø 4mm	302mm	AA-002-014
	HD Cystoscope autoclavable, wide angle	Ø 4mm	302mm	AA-002-014-W
	Ultra HD (4K) Cystoscope autoclavable	Ø 4mm	302mm	AA-002-024

1.1 | CYSTO-URETHROSCOPY

SHEATHS WITH 2 FIXED STOPCOCK & STANDARD OBTURATOR



Ordering no.	Description	French	Colour code
AB-001-002	Cysto-Urethroscope sheath with 2 fixed stopcock and standard obturator	17 Fr.	Yellow
AB-001-003		19 Fr.	Green
AB-001-004		21 Fr.	Red
AB-001-005		23 Fr.	Blue
AB-001-006		25 Fr.	White

SHEATHS WITH CENTRAL VALVE & STANDARD OBTURATOR



Ordering no.	Description	French	Colour code
AB-002-001	Cysto-Urethroscope sheath with central valve and standard obturator	17 Fr.	Yellow
AB-002-002		19 Fr.	Green
AB-002-003		21 Fr.	Red
AB-002-004		23 Fr.	Blue
AB-002-005		25 Fr.	White

1.1 | CYSTO-URETHROSCOPY

STANDARD OBTURATORS



Ordering no.	Description	French	Colour code
AB-001-012	Standard obturator Cysto-Urethroscope sheath	17 Fr.	Yellow
AB-001-013		19 Fr.	Green
AB-001-014		21 Fr.	Red
AB-001-015		23 Fr.	Blue
AB-001-016		25 Fr.	White


VIEWING OBTURATORS




Ordering no.	Description	French	Colour code
AB-001-022	Viewing obturator Cysto-Urethroscope sheath	17 Fr.	Yellow
AB-001-023		19 Fr.	Green
AB-001-024		21 Fr.	Red
AB-001-025		23 Fr.	Blue
AB-001-026		25 Fr.	White

ADAPTER

Ordering no.	Description
AB-001-010	Cystoscope adapter for female urethroscopy




1.1 | CYSTO-URETHROSCOPY ENDOSCOPE BRIDGES




Ordering no.	Description
AB-003-002	Endoscope bridge diagnostic



Ordering no.	Description
AB-003-003	Endoscope bridge with 1 instrument channel



Ordering no.	Description
AB-003-004	Endoscope bridge with 2 instrument channels



Ordering no.	Description
AB-003-015	Endoscope bridge with 2 instrument channels, short



Ordering no.	Description
AJ-010-001	Rubber cap for instrument channel

1.1 | CYSTO-URETHROSCOPY

WORKING ELEMENTS



Ordering no.	Description
AB-003-006	Working element with Albarran deflector and 1 instrument channel
AB-003-008	Working element with Albarran deflector and 1 instrument channel, without ratchet



Ordering no.	Description
AB-003-007	Working element with Albarran deflector and 2 instrument channels
AB-003-009	Working element with Albarran deflector and 2 instrument channels, without ratchet

1.1 | CYSTO-URETHROSCOPY

OPTICAL FORCEPS AND SCISSORS



Info: Only for sheaths 19-25 Fr.



	Ordering no.	For scopes	Description
	AB-004-005	12°/30°	Biopsy forceps (down), double action
	AB-004-006	70°	Biopsy forceps (down), double action
	AB-004-001	0°/12°/30°	Biopsy forceps, double action
	AB-004-002	0°/12°/30°	Biopsy forceps, through cutting
	AB-004-004	0°/12°	Hook scissors, double action
	AB-004-010	0°/12°	Forceps, double action, for removal of stents
	AB-004-003	0°/12°	Grasping forceps, double action
	AB-004-008		Adapter for use with resectoscope sheaths

1.1 | CYSTO-URETHROSCOPY
OPTICAL FORCEPS





Info: Only for sheats 25 Fr.

	Ordering no.	For scopes	Description
	AB-004-007	0°/12°/30°	Optical stone punch
	AB-004-009		Adapter for use with resectoscope sheaths 26 - 28.5 Fr.
	AB-004-009-QL		Adapter for use with resectoscope sheaths 26 - 28.5 Fr.- STEMA Click-System

1.1 | CYSTO-URETHROSCOPY


FORCEPS FLEXIBLE | SEMI-RIGID | RIGID





	Description	Diameter	French	Flexible	Semi-rigid	Rigid
	Biopsy forceps, single action	Ø 2.3mm	7 Fr.	AC-010-004	AC-010-204	
	Grasping forceps, single action	Ø 2.3mm	7 Fr.	AC-011-004	AC-011-204	
	Biopsy forceps oval, double action	Ø 1.0mm	3 Fr.	AC-001-001	AC-001-201	AC-001-001-R
		Ø 1.2mm	4 Fr.	AC-001-008	AC-001-208	AC-001-008-R
		Ø 1.6mm	5 Fr.	AC-001-002	AC-001-202	AC-001-002-R
		Ø 2.1mm	6 Fr.	AC-001-003	AC-001-203	AC-001-003-R
		Ø 2.3mm	7 Fr.	AC-001-004	AC-001-204	AC-001-004-R
		Ø 3.0mm	9 Fr.	AC-001-005	AC-001-205	AC-001-005-R
	Biopsy forceps oval, serrated, double action	Ø 1.0mm	3 Fr.	AC-002-001	AC-002-201	AC-002-001-R
		Ø 1.6mm	5 Fr.	AC-002-002	AC-002-202	AC-002-002-R
		Ø 2.1mm	6 Fr.	AC-002-003	AC-002-203	AC-002-003-R
		Ø 2.3mm	7 Fr.	AC-002-004	AC-002-204	AC-002-004-R
		Ø 3.0mm	9 Fr.	AC-002-005	AC-002-205	AC-002-005-R
	Biopsy forceps oval, serrated, teeth interlock, double action	Ø 1.6mm	5 Fr.	AC-015-002	AC-015-202	AC-015-002-R
		Ø 2.1mm	6 Fr.		AC-015-203	AC-015-003-R
		Ø 2.3mm	7 Fr.	AC-015-004	AC-015-204	AC-015-004-R
	Grasping forceps for stone street, double action	Ø 1.2mm	4 Fr.	AC-017-001	AC-017-201	AC-017-001-R
		Ø 2.3mm	7 Fr.	AC-017-004	AC-017-204	AC-017-004-R
	Grasping forceps, alligator jaw, double action	Ø 1.0mm	3 Fr.	AC-006-001	AC-006-201	AC-006-001-R
		Ø 1.2mm	4 Fr.	AC-006-003	AC-006-203	AC-006-003-R
		Ø 1.6mm	5 Fr.	AC-006-002	AC-006-202	AC-006-002-R
		Ø 2.3mm	7 Fr.	AC-006-004	AC-006-204	AC-006-004-R
		Ø 3.0mm	9 Fr.	AC-006-005	AC-006-205	AC-006-005-R

1.1 | CYSTO-URETHROSCOPY


FORCEPS FLEXIBLE | SEMI-RIGID | RIGID




	Description	Diameter	French	Flexible	Semi-rigid	Rigid
	Grasping forceps, DJ, double action	Ø 2.3mm	7 Fr.	AC-008-004	AC-008-204	AC-008-004-R
		Ø 3.0mm	9 Fr.	AC-008-005	AC-008-205	AC-008-005-R
	Grasping forceps teeth double action	Ø 1.2mm	4 Fr.	AC-016-001	AC-016-201	AC-016-001-R
		Ø 1.6mm	5 Fr.	AC-016-002	AC-016-202	AC-016-002-R
	Grasping forceps teeth interlock, double action	Ø 1.6mm	5 Fr.	AC-020-002	AC-020-202	AC-020-002-R
		Ø 2.3mm	7 Fr.	AC-020-004	AC-020-204	AC-020-004-R
	Grasping forceps rat tooth, double action	Ø 1.0mm	3 Fr.	AC-019-001	AC-019-201	AC-019-001-R
		Ø 1.6mm	5 Fr.	AC-019-002	AC-019-202	AC-019-002-R
	Grasping forceps for stone street, double action	Ø 1.2mm	4 Fr.	AC-007-001	AC-007-201	AC-007-001-R
		Ø 1.6mm	5 Fr.	AC-007-002	AC-007-202	AC-007-002-R
	Grasping forceps for big stones, double action	Ø 1.2mm	4 Fr.	AC-009-001	AC-009-201	AC-009-001-R
		Ø 1.6mm	5 Fr.	AC-009-002	AC-009-202	AC-009-002-R
		Ø 2.3mm	7 Fr.	AC-009-004	AC-009-204	AC-009-004-R
	3-Prongs retriever	Ø 1.0mm	3 Fr.	AC-012-001		
		Ø 1.6mm	5 Fr.	AC-012-002		
		Ø 2.3mm	7 Fr.	AC-012-004		


1.1 | CYSTO-URETHROSCOPY

SCISSORS FLEXIBLE | SEMI-RIGID | RIGID



	Description	Diameter	French	Flexible	Semi-rigid	Rigid
	Hook scissors, single action	Ø 2.3mm	7 Fr.	AC-004-004	AC-004-204	AC-004-004-R
		Ø 3.0mm	9 Fr.	AC-004-005	AC-004-205	AC-004-005-R
	Scissors, single action	Ø 1.0mm	3 Fr.	AC-003-000	AC-003-200	AC-003-000-R
		Ø 1.2mm	4 Fr.	AC-003-001	AC-003-201	AC-003-001-R
		Ø 1.6mm	5 Fr.	AC-003-002	AC-003-202	AC-003-002-R
		Ø 2.1mm	6 Fr.	AC-003-003	AC-003-203	AC-003-003-R
		Ø 2.3mm	7 Fr.	AC-003-004	AC-003-204	AC-003-004-R
		Ø 3.0mm	9 Fr.	AC-003-005	AC-003-205	AC-003-005-R
	Uretero slitting scissors, single action	Ø 1.6mm	4 Fr.		AC-018-201	AC-018-001-R
		Ø 2.3mm	5 Fr.		AC-018-202	AC-018-002-R
	Scissors pointed, single action	Ø 1.6mm	5 Fr.		AC-003-302	
		Ø 2.3mm	7 Fr.		AC-003-304	
	Scissors, serrated flexible, double action	Ø 2.3mm	7 Fr.	AC-005-004	AC-005-204	AC-005-004-R

NEW VERSION OF FORCEPS WITH LUER LOCK AND ROTARY WHEEL




FEATURES

- 1 QUICK AND EASY LOCKING AND RELEASING OF THE HANDLE WITH CLICK LOCK
- 2 ROTARY WHEEL ENABLES PRECISE POSITIONING OF JAW SECTION
- 3 OVERLOAD PROTECTION AVAILABLE
- 4 SIDE-MOUNTED LUER LOCK

1.1 | CYSTO-URETHROSCOPY

INSTRUMENTS SEMI-RIGID | DETACHABLE




Ordering no.	Description
AC-100-001	HANDLE FOR DETACHABLE INSTRUMENTS

	Description	Diameter	French	Complete	Insert
	Scissors, single action	Ø 1.0mm	3 Fr.	AC-003-200-T	AC-003-200-OT
		Ø 1.2mm	4 Fr.	AC-003-201-T	AC-003-201-OT
		Ø 1.6mm	5 Fr.	AC-003-202-T	AC-003-202-OT
		Ø 2.3mm	7 Fr.	AC-003-204-T	AC-003-204-OT
		Ø 3.0mm	9 Fr.	AC-003-205-T	AC-003-205-OT
	Scissors pointed, single action	Ø 1.6mm	5 Fr.	AC-003-612-T	AC-003-612-OT
		Ø 2.3mm	7 Fr.	AC-003-614-T	AC-003-614-OT
	Hook scissors, single action	Ø 2.3mm	7 Fr.	AC-004-204-T	AC-004-204-OT
		Ø 3.0mm	9 Fr.	AC-004-205-T	AC-004-205-OT
	Biopsy forceps oval, double action	Ø 1.0mm	3 Fr.	AC-001-201-T	AC-001-201-OT
		Ø 1.2mm	4 Fr.	AC-001-206-T	AC-001-206-OT
		Ø 1.6mm	5 Fr.	AC-001-202-T	AC-001-202-OT
		Ø 2.3mm	7 Fr.	AC-001-203-T	AC-001-203-OT
		Ø 3.0mm	9 Fr.	AC-001-204-T	AC-001-204-OT
	Grasping forceps, alligator jaw, double action	Ø 1.0mm	3 Fr.	AC-006-201-T	AC-006-201-OT
		Ø 1.2mm	4 Fr.	AC-006-203-T	AC-006-203-OT
		Ø 1.6mm	5 Fr.	AC-006-202-T	AC-006-202-OT
		Ø 2.3mm	7 Fr.	AC-006-204-T	AC-006-204-OT
		Ø 3.0mm	9 Fr.	AC-006-205-T	AC-006-205-OT
	Grasping forceps for stone street, double action	Ø 1.2mm	4 Fr.	AC-007-201-T	AC-007-201-OT
		Ø 1.6mm	5 Fr.	AC-007-202-T	AC-007-202-OT
	Grasping forceps for big stones, double action	Ø 1.2mm	5 Fr.	AC-009-201-T	AC-009-201-OT
		Ø 1.6mm	5 Fr.		AC-009-202-OT

1.1 | CYSTO-URETHROSCOPY

INSTRUMENTS RIGID | DETACHABLE



Ordering no.	Description
AC-100-001	HANDLE FOR DETACHABLE INSTRUMENTS

	Description	Diameter	French	Complete	Insert
	Scissors, single action	Ø 1.0mm	3 Fr.	AC-103-200-T	AC-103-200-OT
		Ø 1.2mm	4 Fr.	AC-103-201-T	AC-103-201-OT
		Ø 1.6mm	5 Fr.	AC-103-202-T	AC-103-202-OT
		Ø 2.3mm	7 Fr.	AC-103-204-T	AC-103-204-OT
		Ø 3.0mm	9 Fr.	AC-103-205-T	AC-103-205-OT
	Scissors pointed, single action	Ø 1.6mm	5 Fr.	AC-103-612-T	AC-103-612-OT
		Ø 2.3mm	7 Fr.	AC-103-614-T	AC-103-614-OT
	Hook scissors, single action	Ø 2.3mm	7 Fr.	AC-104-204-T	AC-104-204-OT
		Ø 3.0mm	9 Fr.	AC-104-205-T	AC-104-205-OT
	Biopsy forceps oval, double action	Ø 1.0mm	3 Fr.	AC-101-201-T	AC-101-201-OT
		Ø 1.2mm	4 Fr.	AC-101-206-T	AC-101-206-OT
		Ø 1.6mm	5 Fr.	AC-101-202-T	AC-101-202-OT
		Ø 2.3mm	7 Fr.	AC-101-203-T	AC-101-203-OT
		Ø 3.0mm	9 Fr.	AC-101-204-T	AC-101-204-OT
	Grasping forceps, alligator jaw, double action	Ø 1.0mm	3 Fr.	AC-106-201-T	AC-106-201-OT
		Ø 1.2mm	4 Fr.	AC-106-203-T	AC-106-203-OT
		Ø 1.6mm	5 Fr.	AC-106-202-T	AC-106-202-OT
		Ø 2.3mm	7 Fr.	AC-106-204-T	AC-106-204-OT
		Ø 3.0mm	9 Fr.	AC-106-205-T	AC-106-205-OT
	Grasping forceps for stone street, double action	Ø 1.2mm	4 Fr.	AC-107-201-T	AC-107-201-OT
		Ø 1.6mm	5 Fr.	AC-107-202-T	AC-107-202-OT
	Grasping forceps for big stones, double action	Ø 1.2mm	5 Fr.	AC-109-201-T	AC-109-201-OT

1.1 | CYSTO-URETHROSCOPY

STONE EXTRACTORS | 4 WIRES



Ordering no.	Description
AF-001-021	Handle for dormia stone extractor 3 and 5 Fr.
AF-001-024	Handle for dormia stone extractor 7 Fr.

Ordering no.	Wires	French	Working length	Description
AF-001-001	4	3 Fr.	60cm	Dormia stone extractor, basket with tubus only
AF-001-002	4	3 Fr.	60cm	Dormia stone extractor complete with handle
AF-001-015	4	3 Fr.	90cm	Dormia stone extractor, basket with tubus only
AF-001-016	4	3 Fr.	90cm	Dormia stone extractor complete with handle
AF-001-023	4	3 Fr.	120cm	Dormia stone extractor, basket with tubus only
AF-001-022	4	3 Fr.	120cm	Dormia stone extractor complete with handle
AF-001-003	4	5 Fr.	60cm	Dormia stone extractor, basket with tubus only
AF-001-004	4	5 Fr.	60cm	Dormia stone extractor complete with handle
AF-001-017	4	5 Fr.	90cm	Dormia stone extractor, basket with tubus only
AF-001-019	4	5 Fr.	90cm	Dormia stone extractor complete with handle
AF-001-005	4	7 Fr.	60cm	Dormia stone extractor, basket with tubus only
AF-001-006	4	7 Fr.	60cm	Dormia stone extractor complete with handle

1.1 | CYSTO-URETHROSCOPY

STONE EXTRACTORS | 5-7 WIRES



Ordering no.	Description
AF-001-021	Handle for dormia stone extractor 3 and 5 Fr.
AF-001-024	Handle for dormia stone extractor 7 Fr.

Ordering no.	Wires	French	Working length	Description
AF-001-007	5	5 Fr.	60cm	Dormia stone extractor, basket with tubus only
AF-001-008	5	5 Fr.	60cm	Dormia stone extractor complete with handle
AF-001-025	5	5 Fr.	90cm	Dormia stone extractor, basket with tubus only
AF-001-018	5	5 Fr.	90cm	Dormia stone extractor complete with handle
AF-001-009	5	7 Fr.	60cm	Dormia stone extractor, basket with tubus only
AF-001-010	5	7 Fr.	60cm	Dormia stone extractor complete with handle
AF-001-011	6	5 Fr.	60cm	Dormia stone extractor, basket with tubus only
AF-001-012	6	5 Fr.	60cm	Dormia stone extractor complete with handle
AF-001-013	6	7 Fr.	60cm	Dormia stone extractor, basket with tubus only
AF-001-014	7	7 Fr.	60cm	Dormia stone extractor complete with handle

1.1 | CYSTO-URETHROSCOPY

MONOPOLAR FLEXIBLE ELECTRODES



	Description	French	Length	Ordering no.	Length	Ordering no.
	Loop electrode flexible, monopolar	5 Fr.	530mm	AD-003-003	650mm	AD-003-013
		7 Fr.	530mm	AD-003-005	650mm	AD-003-015
		9 Fr.	530mm	AD-003-007	650mm	AD-003-017
		10 Fr.	530mm	AD-003-008	650mm	AD-003-018
	Needle electrode flexible, monopolar	3 Fr.	530mm	AD-002-001	650mm	AD-002-011
		4 Fr.	530mm	AD-002-002	650mm	AD-002-012
		5 Fr.	530mm	AD-002-003	650mm	AD-002-013
		7 Fr.	530mm	AD-002-005	650mm	AD-002-015
	Button electrode flexible, monopolar	3 Fr.	530mm	AD-001-001	650mm	AD-001-011
		4 Fr.	530mm	AD-001-002	650mm	AD-001-012
		5 Fr.	530mm	AD-001-003	650mm	AD-001-013
		6 Fr.	530mm	AD-001-004	650mm	AD-001-014
		7 Fr.	530mm	AD-001-005	650mm	AD-001-015
		8 Fr.	530mm	AD-001-006	650mm	AD-001-016
		9 Fr.	530mm	AD-001-007	650mm	AD-001-017
		10 Fr.	530mm	AD-001-008	650mm	AD-001-018


1.1 | CYSTO-URETHROSCOPY

MONOPOLAR CONNECTING CABLES



CONNECTION INSTRUMENT SIDE



Unit Side	Ordering no.	Length	Diameter	Description
	AD-010-001	3m	Ø 4mm	HF Monopolar connecting cable for Berchtold, Martin, Aesculap
	AD-010-006	5m	Ø 4mm	
	AD-010-002	3m	Ø 5mm	HF Monopolar connecting cable for Erbe, ACC/ICC, Storz
	AD-010-007	5m	Ø 5mm	
	AD-010-003	3m	Ø 8mm	HF Monopolar connecting cable for STEMA HF-U Series, Valleylab/Conmed
	AD-010-008	5m	Ø 8mm	
	AD-010-024	4.5m	Ø 2mm	HF Monopolar connecting cable male for STEMA MASTER



	Ordering no.	Description
	EK-015-060	Universal adapter used for monopolar HF cables (For cables: AD-010-001 to AD-010-008)


1.2 | SLIM CYSTO-URETHROSCOPY

SCOPES



Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 30°	Cystoscope autoclavable	Ø 2.9mm	302mm	AA-003-002
	HD Cystoscope autoclavable	Ø 2.9mm	302mm	AA-003-012
	HD Cystoscope autoclavable, wide angle	Ø 2.9mm	302mm	AA-003-012-W

SHEATHS WITH 2 FIXED STOPCOCKS & OBTURATOR



Ordering no.	Description
AB-001-007	Cysto-Urethroscope sheath 15 Fr. with 2 fixed stopcocks, 5 Fr. instrument channel, incl. obturator




Ordering no.	Description
AB-001-000	Cysto-Urethroscope sheath 16 Fr. with 2 fixed stopcocks, Quick-Lock, incl. obturator


1.2 | SLIM CYSTO-URETHROSCOPY

ENDOSCOPE BRIDGES | QUICK-LOCK SYSTEM


Quick-Lock System



Ordering no.	Description
AB-003-010	Endoscope bridge diagnostic used for 16 Fr. sheath



Ordering no.	Description
AB-003-011	Endoscope bridge with 1 instrument channel used for 16 Fr. sheath



Ordering no.	Description
AB-003-012	Endoscope bridge with 2 instrument channels used for 16 Fr. sheath



Ordering no.	Description
AJ-010-001	Rubber cap for instrument channel
AJ-010-003	Rubber cap for instrument channel and URS, Ø 0,6-0,8mm

1.2 | SLIM CYSTO-URETHROSCOPY WORKING ELEMENTS



Ordering no.	Description
AB-003-013	Working element with Albarran deflector and 1 instrument channel used for 16 Fr. sheath
AB-003-014	Working element with Albarran deflector and 2 instrument channels used for 16 Fr. sheath



Ordering no.	Description
AB-003-020	Protective rod for working elements with Albarran deflector

1.3 | LASER CYSTO-URETHROSCOPY

SCOPES



Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 30°	Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-003
	HD Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-013
	ULTRA HD (4K) Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-023

OUTER SHEATH



Ordering no.	Description
AB-001-040	LASER Cysto-Urethroscope outer sheath 23 Fr. with 2 fixed stopcock used for 30° endoscope


INNER SHEATH



Ordering no.	Description
AB-001-041	LASER Cysto-Urethroscope inner sheath with 8 Fr. working channel




1.3 | LASER CYSTO-URETHROSCOPY
OBTURATORS



Ordering no.	Description
AB-001-042	Viewing obturator for LASER cysto-urethroscope sheath AB-001-040
AB-001-043	Obturator for LASER cysto-urethroscope sheath AB-001-040

1.4 | HD MINI COMPACT CYSTOSCOPE



Ordering no.	Product
AA-007-060	HD MINI COMPACT CYSTOSCOPE

Technical data

Eyepiece	30° angled ocular
Sheath	9-12 Fr.
Direction of view	0° viewing angle
Working length	220mm
Instrument channel	5 Fr.
Luer Lock	3 connections
Optic	40.000 px

Delivery contains

HD Mini compact cystoscope

1.5 | FLEXIBLE CYSTOSCOPE FLEX I

- *High luminance*
- *Exact adjustment of focus*
- *Precise up-/down angulation*



Ordering no.	Product
AA-008-002	FLEXIBLE CYSTOSCOPE FLEX I

Technical data

Cystoscope

Insertion tube	Ø 5.5mm
Working Channel	Ø 2.3mm
Working length	380mm
Field of view	100°
Pixel	18.000
Angulation up-down	210° - 120°
Accessories	Suction valve, leak tester

Delivery contains

Flexible cystoscope FLEX I, leak tester, suction valve

INCLUDED IN SCOPE OF DELIVERY:

AA-008-103

Leak tester

AA-008-105

Suction valve



AA-008-103



AA-008-105

1.5 | FLEXIBLE CYSTOSCOPE ACCESSORIES



Ordering no.	Description
AA-008-100	Bridge with Luer-Lock and 2 stop cocks for AA-008-002



Ordering no.	Description
AA-008-102	Suitcase grey for STEMA flexible scopes



Ordering no.	Description
AA-008-103	Leak tester



Ordering no.	Description
AA-008-104	Single use Luer-Lock adapter, 25 pcs./box



Ordering no.	Description
AA-008-105	Suction valve

1.6 CYSTO-URETHROSCOPY SET

1 HD HYSTEROSCOPE

Ø 4.0mm, 30° viewing angle, autoclavable, 302mm working length

1 AA-002-013



2 CYSTO-URETHROSCOPE SHEATH

21 Fr. with 2 fixed stopcock and obturator

2 AB-001-004

3 CYSTO-URETHROSCOPE SHEATH

25 Fr. with 2 fixed stopcock and obturator


3 AB-001-006



4 ENDOSCOPE BRIDGE

with 1 instrument channel

4 AB-003-003



5 WORKING ELEMENT

with Albaran deflector and 2 instrument channels

5 AB-003-007



6 BIOPSY FORCEPS

flexible, oval, double action, Ø 1.6mm, WL 40cm

6 AC-001-002



7 SCISSORS

flexible, single action, Ø 1.6mm, WL 40cm

7 AC-003-002



8 GRASPING FORCEPS ALLIGATOR

flexible, double action, Ø 2.3mm, WL 40cm

8 AC-006-004



1.6 CYSTO-URETHROSCOPY SET

9 BUTTON ELECTRODE

flexible, 5 Fr., monopolar, length 530mm



9 AD-001-003

10 NEEDLE ELECTRODE

flexible, 5 Fr., monopolar, length 530mm



10 AD-002-003



11 HF MONOPOLAR CABLE

HF monopolar connecting cable for Berchtold, Martin, Aesculap, Ø 4mm



11 AD-010-003

12 DORMIA STONE EXTRACTOR

5 wires, Ø 1.6mm, complete with handle, 60cm Working length



12 AF-001-008

13 BIOPSY FORCEPS

double action



13 AB-004-001

14 GRASPING FORCEPS

double action



14 AB-004-003

15 OPTICAL STONE PUNCH



15 AB-004-007

2. RESECTOSCOPY

2.1 | MONOPOLAR RESECTOSCOPY

24/26/28.5 FRENCH
19/22 FRENCH (SLIM)

2.2 | SALINE RESECTOSCOPY

24/26 FRENCH
19/22 FRENCH (SLIM)
17.5 FRENCH (MINI)

2.3 | LASER RESECTOSCOPY

24/26 FRENCH
19/22 FRENCH (SLIM)
17.5 FRENCH (MINI)

2.4 | ACCESSORIES

2.5 | RESECTOSCOPY SETS

2. RESECTOSCOPY

24/26 FR.
28.5 FR.

19/22 FR.

17.5 FR.


COMPLETE RESECTOSCOPE




1 CYTOSCOPIES



2 WORKING ELEMENTS



3 SHEATHS & OBTURATORS



4 ELECTRODES



2.1 | MONOPOLAR RESECTOSCOPY

CYSTOSCOPES AUTOCLAVABLE Ø 4.0MM



Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 0°	Cystoscope autoclavable	Ø 4mm	302mm	AA-002-001
	Cystoscope autoclavable, wide angle	Ø 4mm	302mm	AA-002-001-W
	HD Cystoscope autoclavable	Ø 4mm	302mm	AA-002-011
	HD Cystoscope autoclavable, wide angle	Ø 4mm	302mm	AA-002-011-W
	Ultra HD (4K) Cystoscope autoclavable	Ø 4mm	302mm	AA-002-021

Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 12°	Cystoscope autoclavable	Ø 4mm	302mm	AA-002-002
	HD Cystoscope autoclavable	Ø 4mm	302mm	AA-002-012
	Ultra HD (4K) Cystoscope autoclavable	Ø 4mm	302mm	AA-002-022

Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 30°	Cystoscope autoclavable	Ø 4mm	302mm	AA-002-003
	Cystoscope autoclavable, wide angle	Ø 4mm	302mm	AA-002-003-W
	HD Cystoscope autoclavable	Ø 4mm	302mm	AA-002-013
	HD Cystoscope autoclavable, wide angle	Ø 4mm	302mm	AA-002-013-W
	Ultra HD (4K) Cystoscope autoclavable	Ø 4mm	302mm	AA-002-023

Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 70°	Cystoscope autoclavable	Ø 4mm	302mm	AA-002-004
	Cystoscope autoclavable, wide angle	Ø 4mm	302mm	AA-002-004-W
	HD Cystoscope autoclavable	Ø 4mm	302mm	AA-002-014
	HD Cystoscope autoclavable, wide angle	Ø 4mm	302mm	AA-002-014-W
	Ultra HD (4K) Cystoscope autoclavable	Ø 4mm	302mm	AA-002-024

2.1 | MONOPOLAR RESECTOSCOPY

CONTINUOUS FLOW SHEATHS FIXED



	Ordering no.	Description	French	Colour code
	AH-001-202	Outer sheath with perforation	26 Fr.	Yellow
	AH-001-212	Outer sheath with perforation	28.5 Fr.	Black
	AH-001-204	Outer sheath with upper/lower perforation	26 Fr.	Yellow
	AH-001-203	Outer sheath with surrounding fenestrated perforation	26 Fr.	Yellow
	AH-001-205	Outer sheath with upper fenestrated perforation	26 Fr.	Yellow





Ordering no.	Description	French	Colour code
AH-001-201	Inner sheath with ceramic isolation	26 Fr.	Yellow
AH-001-211	Inner sheath with ceramic isolation	28.5 Fr.	Black

2.1 | MONOPOLAR RESECTOSCOPY

CONTINUOUS FLOW SHEATHS ROTATABLE



	Ordering no.	Description	French	Colour code
	AH-001-302	Outer sheath with perforation, rotatable	26 Fr.	Yellow
	AH-001-312	Outer sheath with perforation, rotatable	28.5 Fr.	Black
	AH-001-304	Outer sheath with upper/lower perforation, rotatable	26 Fr.	Yellow
	AH-001-306	Outer sheath with surrounding fenestrated perforation, rotatable	26 Fr.	Yellow
	AH-001-303	Outer sheath with upper fenestrated perforation, rotatable	26 Fr.	Yellow
	AH-001-305	Outer sheath with upper fenestrated perforation, rotatable	26 Fr.	Yellow



Ordering no.	Description	French	Colour code
AH-001-301	Inner sheath with ceramic isolation, rotatable	26 Fr.	Yellow
AH-001-311	Inner sheath with ceramic isolation, rotatable	28.5 Fr.	Black



Ordering no.	Description	French	Colour code
AH-001-308	Inner sheath without ceramic isolation, with laser channel up to 7 Fr., for use with AH-001-301	26 Fr.	Yellow

2.1 | MONOPOLAR RESECTOSCOPY

STANDARD SHEATHS WITH INTERMITTENT IRRIGATION



Ordering no.	Description	French	Colour code
AH-001-001	Standard resectoscope sheath with 1 fixed stopcock	24 Fr.	Yellow
AH-001-002	Standard resectoscope sheath with 1 fixed stopcock	27 Fr.	Black



Ordering no.	Description	French	Colour code
AH-001-101	Standard resectoscope sheath with central valve	24 Fr.	Yellow
AH-001-102	Standard resectoscope sheath with central valve	27 Fr.	Black

2.1 | MONOPOLAR RESECTOSCOPY OBTURATORS



Ordering no.	Description	French	For outer sheath	Colour code
AH-001-011	Standard Obturator for resectoscope sheaths	24 Fr.	26 Fr.	Yellow
AH-001-012	Standard Obturator for resectoscope sheaths	27 Fr.	28.5 Fr.	Black



Ordering no.	Description	French	For outer sheath	Colour code
AH-001-015	Deflecting Obturator for resectoscope sheaths	24 Fr.	26 Fr.	Yellow
AH-001-016	Deflecting Obturator for resectoscope sheaths	27 Fr.	28.5 Fr.	Black



Ordering no.	Description	French	For outer sheath	Colour code
AH-001-013	Viewing Obturator for resectoscope sheaths	24 Fr.	26 Fr.	Yellow




Ordering no.	Description	French	For outer sheath	Colour code
AH-001-017	Viewing Obturator for resectoscope sheaths with 5 Fr. Instrument channel	24 Fr.	26 Fr.	Yellow
AH-001-018	Viewing Obturator for resectoscope sheaths with 5 Fr. Instrument channel	27 Fr.	28.5 Fr.	Black

2.1 | MONOPOLAR RESECTOSCOPY

WORKING ELEMENTS | PASSIVE CUTTING



Ordering no.	Description
AG-001-001	Working element used for resectoscopes, passive cutting, handle closed , cutting spring assisted
AG-001-001-HF	Working element used for resectoscopes, passive cutting, handle closed , cutting spring assisted, HF connection on top




Ordering no.	Description
AG-001-002	Working element used for resectoscopes, passive cutting, handle open , cutting spring assisted

2.1 | MONOPOLAR RESECTOSCOPY

WORKING ELEMENTS | ACTIVE CUTTING



Ordering no.	Description
AG-001-003	Working element used for resectoscopes, active cutting, handle closed , cutting again spring assistance
AG-001-003-HF	Working element used for resectoscopes, active cutting, handle closed , cutting again spring assistance, HF connection on top




Ordering no.	Description
AG-001-004	Working element used for resectoscopes, active cutting, handle open , cutting spring assisted

2.1 | MONOPOLAR RESECTOSCOPY

WORKING ELEMENTS | TITAN



Ordering no.	Description
AG-001-007	Working element used for resectoscopes, passive cutting , handle closed, cutting spring assisted



Ordering no.	Description
AG-001-008	Working element used for resectoscopes, active cutting , handle closed, cutting again spring assistance

2.1 | MONOPOLAR RESECTOSCOPY

WORKING ELEMENTS | GYRUS



Ordering no.	Description
AG-001-017	GYRUS Saline Working element used for resectoscopes, active cutting, handle closed TITANIUM VERSION




Ordering no.	Description
AG-001-018	GYRUS Saline Working element used for resectoscopes, passive cutting, handle closed TITANIUM VERSION

2.1 | MONOPOLAR RESECTOSCOPY

ELECTRODES FOR MONOPOLAR RESECTOSCOPY




	Ordering no.	French	Description
	AE-001-001	24 Fr. (yellow)	Cutting loop electrode 30° angled
	AE-002-001	27 Fr. (brown)	
	AE-001-002	24 Fr. (yellow)	Cutting loop electrode angled
	AE-002-002	27 Fr. (brown)	
	AE-001-003	24 Fr. (yellow)	Cutting loop electrode straight
	AE-002-003	27 Fr. (brown)	
	AE-001-004	24 Fr. (yellow)	Rollerball electrode Ø 3mm
	AE-002-004	27 Fr. (brown)	
	AE-001-005	24 Fr. (yellow)	Rollerball electrode Ø 5mm
	AE-002-005	27 Fr. (brown)	
	AE-001-010	24 Fr. (yellow)	Roller electrode Ø 3mm
	AE-002-010	27 Fr. (brown)	
	AE-001-011	24 Fr. (yellow)	Roller electrode Ø 5mm
	AE-002-011	27 Fr. (brown)	
	AE-001-012	24 Fr. (yellow)	Spike electrode Ø 3mm
	AE-002-012	27 Fr. (brown)	
	AE-001-013	24 Fr. (yellow)	Spike electrode Ø 5mm
	AE-002-013	27 Fr. (brown)	

2.1 | MONOPOLAR RESECTOSCOPY

ELECTRODES FOR MONOPOLAR RESECTOSCOPY



	Ordering no.	French	Description
	AE-001-014	24 Fr. (yellow)	Vaporisation electrode Ø 3mm
	AE-002-014	27 Fr. (brown)	
	AE-001-015	24 Fr. (yellow)	Vaporisation electrode Ø 5mm
	AE-002-015	27 Fr. (brown)	
	AE-001-009	24 Fr. (yellow)	Conical electrode
	AE-001-016	24 Fr. (yellow)	Knife electrode
	AE-002-016	27 Fr. (brown)	
	AE-001-006	24 Fr. (yellow)	Band electrode angled
	AE-002-006	27 Fr. (brown)	
	AE-001-007	24 Fr. (yellow)	Band electrode 30° angled
	AE-002-007	27 Fr. (brown)	
	AE-001-008	24 Fr. (yellow)	Band vaporization electrode angled
	AE-001-017	24 Fr. (yellow)	Curette sharp, round (not for use with HF)



Ordering no.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives

2.1 | MONOPOLAR SLIM RESECTOSCOPY

CYSTOSCOPES AUTOCLAVABLE Ø 2.9MM



Viewing angle	Endoscope	Diameter	Working length	Ordering no.
	Cystoscope autoclavable	Ø 2.9mm	302mm	AA-003-002
	HD Cystoscope autoclavable	Ø 2.9mm	302mm	AA-003-012
30°				

2.1 | MONOPOLAR SLIM RESECTOSCOPY

CONTINUOUS FLOW SHEATHS | FIXED



Ordering no.	Description	French	Colour code
AH-002-212	Outer sheath with perforation, fixed	22 Fr.	White



Ordering no.	Description	French	Colour code
AH-002-201	Inner sheath with ceramic isolation, fixed	19 Fr.	White

2.1 | MONOPOLAR SLIM RESECTOSCOPY

CONTINUOUS FLOW SHEATHS | ROTATABLE



	Ordering no.	Description	French	Colour code
	AH-002-302	Outer sheath with perforation, rotatable	22 Fr.	White



Ordering no.	Description	French	Colour code
AH-002-301	Inner sheath with ceramic isolation, rotatable	22 Fr.	White

2.1 | MONOPOLAR SLIM RESECTOSCOPY STANDARD SHEATHS



	Ordering no.	Description	French	Colour code
	AH-002-001	Standard resectoscope sheath, 1 fixed stopcock , without obturator	19 Fr.	White

2.1 | MONOPOLAR SLIM RESECTOSCOPY OBTURATORS



Ordering no.	Description	French	Colour code
AH-002-011	Standard Obturator	19 Fr.	White




Ordering no.	Description	French	Colour code
AH-002-013	Viewing Obturator	19 Fr.	White

2.1 | MONOPOLAR SLIM RESECTOSCOPY

WORKING ELEMENTS | PASSIVE CUTTING



Ordering no.	Description
AG-002-001	Working element used for resectoscopes, passive cutting, handle closed , cutting against spring assistance




Ordering no.	Description
AG-002-002	Working element used for resectoscopes, passive cutting, handle open , cutting against spring assistance

2.1 | MONOPOLAR SLIM RESECTOSCOPY

WORKING ELEMENTS | ACTIVE CUTTING




Ordering no.	Description
AG-002-003	Working element used for slim resectoscopes, active cutting, handle closed , cutting against spring assistance



Ordering no.	Description
AG-002-004	Working element used for slim resectoscopes, active cutting, handle open , cutting against spring assistance

2.1 | MONOPOLAR SLIM RESECTOSCOPY

WORKING ELEMENTS | TITAN



Ordering no.	Description
AG-002-006	Working element used for slim resectoscopes, active cutting , handle closed, TITANIUM VERSION

2.1 | MONOPOLAR SLIM RESECTOSCOPY

ELECTRODES FOR MONOPOLAR RESECTOSCOPY




	Ordering no.	French	Description
	AE-003-001	19 Fr.	Cutting loop electrode 30° angled
	AE-003-003	19 Fr.	Cutting loop electrode straight
	AE-003-004	19 Fr.	Rollerball electrode Ø 3mm
	AE-003-010	19 Fr.	Roller electrode Ø 3mm
	AE-003-014	19 Fr.	Vaporization electrode Ø 3mm
	AE-003-016	19 Fr.	Knife electrode




Ordering no.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives


2.1 | MONOPOLAR SLIM RESECTOSCOPY

MONOPOLAR CONNECTING CABLES



CONNECTION INSTRUMENT SIDE



Unit Side	Ordering no.	Length	Diameter	Description
	AD-010-001	3m	Ø 4mm	HF Monopolar connecting cable for Berchtold, Martin, Aesculap
	AD-010-006	5m	Ø 4mm	
	AD-010-002	3m	Ø 5mm	HF Monopolar connecting cable for Erbe, ACC/ICC, Storz
	AD-010-007	5m	Ø 5mm	
	AD-010-003	3m	Ø 8mm	HF Monopolar connecting cable for STEMA HF-U Series, Valleylab/Conmed
	AD-010-008	5m	Ø 8mm	
	AD-010-024	4.5m	Ø 2mm	HF Monopolar connecting cable male for STEMA MASTER



	Ordering no.	Description
	EK-015-060	Universal adapter used for monopolar HF cables (For cables: AD-010-001 to AD-010-008)

2.1 | MONOPOLAR SLIM RESECTOSCOPY CONNECTING CABLES MONOPOLAR



CONNECTION UNIT SIDE



Instrument Side	Ordering no.	Length	Diameter	Description
	AD-010-025	4.5m	90° angled connector	HF Monopolar connecting cable for STORZ resectoscope
	AD-010-026	4.5m	Ø 3mm	HF Monopolar connecting cable for OLYMPUS resectoscope



RECOMMENDED PRODUCTS:

CI-001-030 STEMA HF-U PRO TOUCH

- For all monopolar & bipolar techniques
- with 7" touch screen display
- provided to perform the most advanced electrosurgery
- ADC / APC System

CI-001-040 DOUBLE PEDAL FOOTSWITCH

2.2 | SALINE RESECTOSCOPY

CYSTOSCOPES AUTOCLAVABLE Ø 4.0MM





Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 12°	HD Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-012

Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 30°	Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-003
	HD Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-013
	ULTRA HD Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-023

2.2 | SALINE RESECTOSCOPY

CONTINUOUS FLOW SHEATHS | ROTATABLE



	Ordering no.	Description	French	Colour code
	AH-001-302	Outer sheath with perforation, rotatable	26 Fr.	Yellow
	AH-001-312	Outer sheath with perforation, rotatable	28.5 Fr.	Black
	AH-001-304	Outer sheath with upper/lower perforation, rotatable	26 Fr.	Yellow
	AH-001-306	Outer sheath with surrounding fenestrated perforation, rotatable	26 Fr.	Yellow
	AH-001-303	Outer sheath with upper fenestrated perforation, rotatable	26 Fr.	Yellow
	AH-001-305	Outer sheath with upper fenestrated perforation, rotatable	26 Fr.	Yellow



Ordering no.	Description	French	Colour code
AH-001-301	Inner sheath with ceramic isolation, rotatable	26 Fr.	Yellow
AH-001-311	Inner sheath with ceramic isolation, rotatable	28.5 Fr.	Black

2.2 | SALINE RESECTOSCOPY OBTURATORS



Ordering no.	Description	French	For outer sheath	Colour code
AH-001-011	Standard Obturator for resectoscope sheaths	24 Fr.	26 Fr.	Yellow
AH-001-012	Standard Obturator for resectoscope sheaths	27 Fr.	28.5 Fr.	Black



Ordering no.	Description	French	For outer sheath	Colour code
AH-001-015	Deflecting Obturator for resectoscope sheaths	24 Fr.	26 Fr.	Yellow
AH-001-016	Deflecting Obturator for resectoscope sheaths	27 Fr.	28.5 Fr.	Black




Ordering no.	Description	French	For outer sheath	Colour code
AH-001-013	Viewing Obturator for resectoscope sheaths	24 Fr.	26 Fr.	Yellow



Ordering no.	Description	French	For outer sheath	Colour code
AH-001-017	Viewing Obturator for resectoscope sheaths with 5 Fr. Instrument channel	24 Fr.	26 Fr.	Yellow
AH-001-018	Viewing Obturator for resectoscope sheaths with 5 Fr. Instrument channel	27 Fr.	28.5 Fr.	Black

2.2 | SALINE RESECTOSCOPY WORKING ELEMENTS



Ordering no.	Description
AG-001-011	SALINE working element used for resectoscopes, 24 Fr., passive cutting , handle closed, cutting spring assistance




Ordering no.	Description
AG-001-012	SALINE working element used for resectoscopes, 24 Fr. , active cutting , handle closed, cutting spring assistance

2.2 | SALINE RESECTOSCOPY

ELECTRODES FOR SALINE RESECTOSCOPY




	Ordering no.	French	Description
	AE-001-030	24 Fr.	Cutting loop electrode 30° angled, strong cut
	AE-001-031	24 Fr.	Cutting loop electrode 30° angled
	AE-001-033	24 Fr.	Cutting loop electrode 90 ° angled
	AE-001-048	24 Fr.	Cutting loop electrode 90 ° angled, for 0° Scope
	AE-001-036	24 Fr.	Cutting loop electrode straight
	AE-001-032	24 Fr.	Knife electrode
	AE-001-041	24 Fr.	Band electrode, 30° angled
	AE-001-038	24 Fr.	Vaporisation spike rollerball electrode Ø 3mm
	AE-001-037	24 Fr.	Vaporisation rollerball electrode Ø 3mm
	AE-001-043	24 Fr.	Vaporisation rollerball electrode Ø 5mm
	AE-001-040	24 Fr.	Roller electrode smooth Ø 3mm
	AE-001-042	24 Fr.	Roller electrode smooth Ø 5mm


2.2 | SALINE RESECTOSCOPY

ELECTRODES FOR SALINE RESECTOSCOPY



	Ordering no.	French	Description
	AE-001-034	24 Fr.	Rollerball electrode Ø 3mm
	AE-001-035	24 Fr.	Rollerball electrode Ø 5mm
	AE-001-044	24 Fr.	Conical electrode „Mushroom“
	AE-001-045	24 Fr.	Conical electrode „Mushroom“, convex
	AE-001-046	24 Fr.	Conical electrode „Mushroom“, concav
	AE-001-047	24 Fr.	Conical electrode „Mushroom“, grooved
	AE-001-039	24 Fr.	Conical electrode

	Ordering no.	French	Description
	AE-002-031	27 Fr.	Cutting loop electrode, 30° angled



Ordering no.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives

2.2 | SALINE SLIM RESECTOSCOPY

CYSTOSCOPES AUTOCLAVABLE Ø 2.9MM



Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 12°	Cystoscope autoclavable	Ø 2.9mm	302mm	AA-003-004
	HD Cystoscope autoclavable	Ø 2.9mm	302mm	AA-003-014

Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 30°	Cystoscope autoclavable	Ø 2.9mm	302mm	AA-003-002
	HD Cystoscope autoclavable	Ø 2.9mm	302mm	AA-003-012

2.2 | SALINE SLIM RESECTOSCOPY SHEATHS



Ordering no.	Description
AH-002-302	Outer sheath 22 Fr. with perforation, rotatable



Ordering no.	Description
AH-002-301	Inner sheath 22 Fr. with ceramic isolation, rotatable

2.2 | SALINE SLIM RESECTOSCOPY OBTURATORS




Ordering no.	Description
AH-002-011	Standard obturator for resectoscope sheaths 19 Fr.




Ordering no.	Description
AH-002-013	Viewing Obturator for resectoscope sheaths 19 Fr.

2.2 | SALINE SLIM RESECTOSCOPY

WORKING ELEMENTS



Ordering no.	Description
AG-002-011	SALINE working element used for resectoscopes, 19 Fr., passive cutting , handle closed, cutting spring assistance



Ordering no.	Description
AG-002-012	SALINE working element used for resectoscopes, 19 Fr., active cutting , handle closed, cutting spring assistance

2.2 | SALINE SLIM RESECTOSCOPY

ELECTRODES FOR SALINE RESECTOSCOPY



	Ordering no.	French	Description
	AE-003-031	19 Fr.	Cutting loop electrode 30° angled
	AE-003-032	19 Fr.	Cutting loop electrode 90° angled
	AE-003-033	19 Fr.	Cutting loop electrode straight
	AE-003-034	19 Fr.	Rollerball electrode Ø 3mm
	AE-003-040	19 Fr.	Roller electrode Ø 3mm
	AE-003-036	19 Fr.	Knife electrode



Ordering no.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives




17.5 FR.

2.2 | SALINE MINI RESECTOSCOPY
CYSTOSCOPES



Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 30°	Cystoscope autoclavable	Ø 2.9mm	302mm	AA-003-002
	HD Cystoscope autoclavable	Ø 2.9mm	302mm	AA-003-012

2.2 | SALINE MINI RESECTOSCOPY WORKING ELEMENTES



Ordering no.	Description
AG-002-015	SALINE working element used for resectoscopes, 18.5 Fr., passive cutting , handle closed, cutting spring assistance



Ordering no.	Description
AG-002-013	SALINE working element used for resectoscopes, 18.5 Fr., active cutting , handle closed, cutting spring assistance

2.2 | SALINE MINI RESECTOSCOPY
SHEATHS



Ordering no.	Description
AH-002-213	Outer sheath 18.5 Fr. with perforation



Ordering no.	Description
AH-002-211	Inner sheath 17.5 Fr. with ceramic isolation, rotatable



Ordering no.	Description
AH-002-216	Inner sheath 17.5 Fr. without ceramic isolation , with laser channel up to 5 Fr., for use with AH-001-213

2.2 | SALINE MINI RESECTOSCOPY OBTURATORS



Ordering no.	Description
AH-002-015	Standard obturator 17.5 Fr.





Ordering no.	Description
AH-002-016	Viewing Obturator 17.5 Fr. with 5 fr. instrument channel

2.2 | SALINE MINI RESECTOSCOPY

ELECTRODES FOR BIPOLAR RESECTOSCOPY




	Ordering no.	French	Description
	AE-003-041	17.5 Fr.	Cutting loop electrode 30° angled
	AE-003-043	17.5 Fr.	Cutting loop electrode straight
	AE-003-044	17.5 Fr.	Rollerball electrode Ø 2mm
	AE-003-045	17.5 Fr.	Roller electrode Ø 2mm
	AE-003-046	17.5 Fr.	Knife electrode



Ordering no.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives

2.2 | SALINE RESECTOSCOPY

CONNECTING CABLES | BIPOLAR



CONNECTION INSTRUMENT SIDE



Connection Unit Side	Ordering no.	Length	Description
	AD-010-011	4.5m	Bipolar resectoscope connecting cable for Martin Maxium
	AD-010-012	4.5m	Bipolar resectoscope connecting cable for Erbe Vio
	AD-010-013	4.5m	Bipolar resectoscope connecting cable for Valleylab, Lamidey, EMC
	AD-010-015	4.5m	Bipolar resectoscope connecting cable for Erbe Vio 3

2.3 | LASER RESECTOSCOPY
CYSTOSCOPES AUTOCLAVABLE Ø 4.0MM



Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 30°	HD Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-013

2.3 | LASER RESECTOSCOPY SHEATHS | OBTURATORS



Ordering no.	Description	French	Colour code
AH-001-302	Outer sheath with perforation, rotatable	26 Fr.	Yellow



Ordering no.	Description	French	Colour code
AH-001-301	Inner sheath with ceramic isolation, rotatable	26 Fr.	Yellow



Ordering no.	Description	French	Colour code
AH-001-308	Inner sheath without ceramic isolation, with laser channel up to 7 Fr., for use with AH-001-301	26 Fr.	Yellow



Ordering no.	Description	French	For outer sheath	Colour code
AH-001-011	Standard Obturator for resectoscope sheaths	24 Fr.	26 Fr.	Yellow

2.3 | LASER RESECTOSCOPY WORKING ELEMENTS



Ordering no.	Description
AG-001-006	Working element used for resectoscopes, titanium, passive cutting, handle closed, 24/26 Fr. with probe channel 0.8mm
AG-001-009	Working element used for resectoscopes, titanium, passive cutting, handle closed, 24/26 Fr. with probe channel 1.2mm

Ordering no.	Description
AG-100-006	Adapter used for laser probes
AG-100-007	TUOHY-Borst adapter
AG-001-013	Adapator used for laser probe, working element used for LASER resectoscopes



19/22 FR.

2.3 | LASER SLIM RESECTOSCOPY
CYSTOSCOPES AUTOCLAVABLE Ø 2.9MM



Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 30°	HD Cystoscope autoclavable	Ø 2.9mm	302mm	AA-003-012

2.3 | LASER SLIM RESECTOSCOPY SHEATHS



Ordering no.	Description
AH-002-302	Outer sheath 22 Fr. with perforation, rotatable



Ordering no.	Description
AH-002-301	Inner sheath 22 Fr. with ceramic isolation, rotatable

2.3 | LASER SLIM RESECTOSCOPY OBTURATORS



Ordering no.	Description
AH-002-011	Standard obturator for resectoscope sheaths 19 Fr.



Ordering no.	Description
AH-002-013	Viewing Obturator for resectoscope sheaths 19 Fr.

2.3 | LASER SLIM RESECTOSCOPY WORKING ELEMENTS



Ordering no.	Description
AG-002-007	Working element used for resectoscopes, titanium, passive cutting, handle closed, 19/22 Fr. with probe channel 0.8mm
AG-002-008	Working element used for resectoscopes, titanium, passive cutting, handle closed, 19/22 Fr. with probe channel 1.2mm

2.3 | LASER MINI RESECTOSCOPY
CYSTOSCOPES AUTOCLAVABLE Ø 2.9MM



Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 30°	HD Cystoscope autoclavable	Ø 2.9mm	302mm	AA-003-012

2.3 | LASER MINI RESECTOSCOPY SHEATHS



Ordering no.	Description
AH-002-213	Outer sheath 18.5 Fr. with perforation



Ordering no.	Description
AH-002-211	Inner sheath 17.5 Fr. with ceramic isolation, rotatable



Ordering no.	Description
AH-002-216	Inner sheath 17.5 Fr. without ceramic isolation , with laser channel up to 5 Fr., for use with AH-001-213

2.3 | LASER SLIM RESECTOSCOPY OBTURATORS



Ordering no.	Description
AH-002-015	Standard obturator 17.5 Fr.



Ordering no.	Description
AH-002-016	Viewing Obturator 17.5 Fr. with 5 fr. instrument channel


2.3 | LASER SLIM RESECTOSCOPY WORKING ELEMENTS



Ordering no.	Description
AG-002-016	Working element used for resectoscopes, titanium, passive cutting, handle closed, 17.5 Fr. with probe channel 0.8mm

2.4 | ACCESSORIES

SUPRAPUBIC ASPIRATION



Ordering no.	Description
AB-007-001	Suprapubic aspiration set, complete set



Ordering no.	Description
AB-007-002	Trocarn sleeve 16 Fr. with pressure measuring channel



Ordering no.	Description
AB-007-003	Adjustable stabilisation disc, lateral opening



Ordering no.	Description
AB-007-008	Guide sleeve, lateral opening , for inserting a suprapubic catheter

2.4 | ACCESSORIES

SUPRAPUBIC ASPIRATION



Ordering no.	Description
AB-007-004	Trocar with cannula tip



Ordering no.	Description
AB-007-005	Trocar with pyramidal tip



Ordering no.	Description
AB-007-006	Aspiration cannula



Ordering no.	Description
AB-007-007	Overflow cannula for pressure stabilisation

2.4 | ACCESSORIES

EVACUATOR




Ordering no.	Description
AJ-002-001	Ellik evacuator complete with sheath connection and locking device
AJ-002-002	Bowl only
AJ-002-003	Suction Ball only
AJ-002-004	Cone with tube only
AJ-002-005	Sheath connection for Ellik only

BLADDER SYRINGES




Not available in the EU




Ordering no.	Description
AJ-001-000	RECORD syringe with adapter, 50ml
AJ-001-010	Replacement cylinder for RECORD syringe with adapter, 50ml
AJ-001-001	RECORD syringe with adapter, 75ml
AJ-001-011	Replacement cylinder for RECORD syringe with adapter, 75ml
AJ-001-002	RECORD syringe with adapter, 100ml
AJ-001-012	Replacement cylinder for RECORD syringe with adapter, 100ml
AJ-001-003	RECORD syringe with adapter, 150ml
AJ-001-013	Replacement cylinder for RECORD syringe with adapter, 150ml

2.4 | ACCESSORIES

TRAY INSERT



Ordering no.	Description
AJ-010-020	Tray insert for resectoscopy (without instruments, without sterilization box) , dimension 370x255x45mm (LxWxH)



RECOMMENDED PRODUCTS:

- AA-002-013** HD autoclavable Scope 30° | 4mm
- AG-001-011** SALINE working element, passive cutting, handle closed
- AH-001-302** Standard resectoscope sheath 24 Fr., with 1 fixed stopcock
- AH-001-301** Inner sheath 26 Fr. with ceramic isolation, rotatable
- AH-001-011** Standard obturator for resectoscope sheaths 24 Fr.
- AJ-010-006** Protection sleeve made of steel for sterilization and storage of electrodes and cold knives
- AE-001-031** Saline Cutting loop electrode 30° angled, 24 Fr. wire thickness 0,35mm
- AE-001-034** Saline Rollerball electrode 3mm, 24 Fr.
+ Connecting cable (depends on unit)

2.5 MONOPOLAR RESECTOSCOPY SET

1 HD CYSTOSCOPE

Ø 4.0mm, 30° viewing angle, autoclavable,
WL 302mm

1 AA-002-013



2 WORKING ELEMENT

for resectoscopes, passive cutting, handle closed

2 AG-001-001



3 OUTER SHEATH

26 Fr. with perforation, rotatable

3 AH-001-302



4 INNER SHEATH

26 Fr. with ceramic isolation, rotatable

4 AH-001-301



5 STANDARD OBTURATOR

24 Fr.

5 AH-001-011



6 CUTTING LOOP ELECTRODE

30° angled, 24 Fr.

6 AE-001-001



7 ROLLERBALL ELECTRODE 3 MM

24 Fr.

7 AE-001-004



8 KNIFE ELECTRODE

24 Fr.

8 AE-001-016



2.5 MONOPOLAR RESECTOSCOPY SET

9 HF MONOPOLAR CONNECTING CABLE
for Berchtold, Martin, Aesculap, Ø 4mm

9 AD-010-001



10 RECORD SYRINGE
100ml

10 AJ-001-002



11 ELLIK EVACUATOR
complete

11 AJ-002-001



12 BIOPSY FORCEPS
double action

12 AB-004-001




13 GRASPING FORCEPS
double action

13 AB-004-003



14 Adapter
for use with resectoscope sheaths

14 AB-004-008



15 PROTECTION SLEEVE
For sterilization and storage of electrodes

15 AJ-010-006



2.5 SALINE RESECTOSCOPY SET

1 HD CYSTOSCOPE

Ø 4.0mm, 30° viewing angle, autoclavable,
WL 302mm

1 AA-002-013



2 SALINE WORKING ELEMENT

for resectoscopes, passive cutting, handle closed

2 AG-001-011



3 OUTER SHEATH

26 Fr. with perforation, rotatable

3 AH-001-302



4 INNER SHEATH

26 Fr. with ceramic isolation, rotatable

4 AH-001-301



5 STANDARD OBTURATOR

24 Fr.

5 AH-001-011



6 CUTTING LOOP ELECTRODE

30° angled, 24 Fr.

6 AE-001-031



7 ROLLERBALL ELECTRODE 3 MM

24 Fr.

7 AE-001-034



8 KNIFE ELECTRODE

24 Fr.

8 AE-001-032



2.5 SALINE RESECTOSCOPY SET

- 9 HF BIPOLAR CONNECTING CABLE**
for Valleylab, Lamidey, EMC, 4.5mm

9 AD-010-013



- 10 RECORD SYRINGE**
100 ml

10 AJ-001-002



- 11 ELLIK EVACUATOR**
complete

11 AJ-002-001



- 12 BIOPSY FORCEPS**
double action

12 AB-004-001




- 13 GRASPING FORCEPS**
double action

13 AB-004-003



- 14 ADAPTER**
for use with resectoscope sheaths

14 AB-004-008



- 15 PROTECTION SLEEVE**
for sterilization and storage of electrodes

15 AJ-010-006



3. URETHROTOMY

3.1 | URETHROTOME

3.2 | OPTICAL URETHROTOMY

CYSTOSCOPES

SHEATHS

OBTURATORS

WORKING ELEMENTS

KNIVES

3.3 | URETHROTOMY SET

3. URETHROTOMY


4MM

2.9MM


1 CYSTOSCOPES




2 WORKING ELEMENTS



3 SHEATHS & OBTURATORS



4 KNIVES



3.1 | URETHROTOME



Ordering no.	Description
AI-002-001	OTIS urethrotome with 3 exchangeable tips inclusive 2 knives

SPARE KNIVES



Ordering no.	Description
AI-002-002	Spare knives for OTIS urethrotome


3.2 | OPTICAL URETHROTOMY

CYSTOSCOPES



Viewing angle	Endoscope	Diameter	Working length	Ordering no.
 0°	Cystoscope autoclavable	Ø 4mm	302mm	AA-002-001
	HD Cystoscope autoclavable	Ø 4mm	302mm	AA-002-011
	Ultra HD (4K) Cystoscope autoclavable	Ø 4mm	302mm	AA-002-021

SHEATH



Ordering no.	Description
AI-001-001	Urethrotome sheath 21 Fr. with instrument channel 5 Fr.
AI-001-003	Endoscope bridge with 1 instrument channel , fit with urethrotome sheath AI-001-001

3.2 | OPTICAL URETHROTOMY SUPPLEMENTARY SHEATHS



Ordering no.	Description
AI-001-002	Supplementary sheath for AI-001-001, open to side



Ordering no.	Description
AI-001-006	Supplementary sheath 25 Fr. for continuous flow

3.2 | OPTICAL URETHROTOMY

OBTURATORS



Ordering no.	Description
AI-001-004	Standard obturator for urethrotome sheath 21 Fr.



Ordering no.	Description
AI-001-005	Obturator for urethrotome sheath 21 Fr. with 5 Fr. instrument channel



Ordering no.	Description
AI-001-012	Viewing obturator for urethrotome sheath 21 Fr.

3.2 | OPTICAL URETHROTOMY

WORKING ELEMENTS | PASSIVE CUTTING



Ordering no.	Description
AG-001-001	Working element used for resectoscopes, passive cutting, handle closed , cutting spring assisted
AG-001-001-HF	Working element used for resectoscopes, passive cutting, handle closed , cutting spring assisted, HF connection on top




Ordering no.	Description
AG-001-002	Working element used for resectoscopes, passive cutting, handle open , cutting spring assisted

3.2 | OPTICAL URETHROTOMY

COLD KNIVES



	Ordering no.	Description
	AI-001-007	Cold knife straight
	AI-001-008	Cold knife straight, wave shape
	AI-001-009	Cold knife round
	AI-001-011	Cold knife hook shaped



Ordering no.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives


3.2 | OPTICAL URETHROTOMY

CYSTOSCOPES



Viewing angle	Description	Diameter	Working length	Ordering no.
 0°	Cystoscope autoclavable	Ø 2.9mm	302mm	AA-003-001
	HD Cystoscope autoclavable	Ø 2.9mm	302mm	AA-003-011

SHEATH



Ordering no.	Description
AI-005-001	Urethrotome sheath 15.5 Fr. with instrument channel 3 Fr
AI-005-003	Endoscope bridge with 1 instrument channel , fit with urethrotome sheath AI-005-001

OBTURATOR



Ordering no.	Description
AI-005-004	Standard obturator for urethrotome sheath 15.5 Fr.

3.2 | OPTICAL URETHROTOMY

WORKING ELEMENT | PASSIVE CUTTING




Ordering no.	Description
AG-001-010	Working element used for resectoscopes, passive cutting, handle closed , without HF-connection

3.2 | OPTICAL URETHROTOMY

COLD KNIVES



	Ordering no.	Description
	AI-005-007	Cold knife straight
	AI-005-008	Cold knife straight, wave shape
	AI-005-011	Cold knife hook shaped



Ordering no.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives

3.4 OPTICAL UREHTROTOMY SET

1 HD CYSTOSCOPE

Ø 4.0mm, 0° viewing angle, autoclavable,
302mm working length

1 AA-002-011



2 WORKING ELEMENT

for resectoscopes, passive cutting, handle closed

2 AG-001-001



3 URETHROTOME SHEATH

21 Fr. with instrument channel 5 Fr.

3 AI-001-001



4 SUPPLEMENTARY SHEATH

25 Fr. for continous flow

4 AI-001-002



5 STANDARD OBTURATOR

24 Fr.


5 AI-001-004



6 ENDOSCOPE BRIDGE

with 1 instrument channel

6 AI-001-003



7 COLD KNIFE

straight

7 AI-001-007



8 PROTECTION SLEEVE

for sterilization and storage of electrodes

8 AJ-010-006



4. NEPHROSCOPY

4.1 | HD NEPHROSCOPY

HD NEPHROSCOPES
HD MINI NEPHROSCOPE
SHEATHS
OBTURATORS

4.2 | MINI NEPHROSCOPY

MINI NEPHROSCOPES
SHEATHS
DILATATORS

4.3 | STANDARD NEPHROSCOPY

NEPHROSCOPES 200MM
SHEATHS
OBTURATORS

NEPHROSCOPES 250MM
SHEATHS
OBTURATORS

TELESCOPE BOUGIES
INSTRUMENTS

4. NEPHROSCOPY

1 HD NEPHROSCOPE



2 SHEATHS & OBTURATORS



3 FORCEPS, INSTRUMENTS, KNIVES



4.1 | HD NEPHROSCOPES



Ordering no.	Working length	Viewing angle	Description
AA-006-069	220mm	0°	HD Nephroscope - 45° angled ocular - Ø 18.8 Fr. - working channel Ø 12.9 Fr.
AA-006-070		15°	

SHEATHS



Ordering no.	French	Description
AB-006-170	23 Fr.	Nephroscope sheath for continous flow, rotatable stopcock (fit with nephrosopes WL 220mm / scope 15°)

OBTURATOR




Ordering no.	French	Description
AB-006-171	23 Fr.	Obturator hollow (fit with nephroscopy sheat AB-006-170 / scope 15°)

SEALING CAP


Ordering no.	Description
AA-006-081	Sealing cap, Ø 15mm/1.0mm
AA-006-083	Sealing cap, Ø 20mm/2.0mm

4.1 | MINI HD COMPACT NEPHROSCOPE



Ordering no.	Working length	Viewing angle	Description
AA-006-100	220mm	0°	Mini HD compact nephroscope with 30° angled ocular, sheath 12 Fr. , instrument channel 6.6 Fr. , 3 luerlock connections, 40.000 px

4.2 | MINI NEPHROSCOPY
MINI NEPHROSCOPES




Ordering no.	French	Working length	Viewing angle	Description
AA-006-050	Ø 19 Fr.	220mm	8°	Mini nephroscope with 45° angled ocular, instrument channel 10 Fr.
AA-006-051	Ø 13.5 Fr.	220mm	8°	Mini nephroscope with 45° angled ocular, instrument channel 6 Fr.
AA-006-052	Ø 13.5 Fr.	320mm	8°	Mini nephroscope with 45° angled ocular, instrument channel 6 Fr.

	Ordering no.	Description
	AA-006-053	Instrument bridge used for Mini-Nephroscopes (not included in Mini Nephroscope)
	AJ-010-001 AJ-010-003	Rubber cap for instrument channel Rubber cap for instrument channel and URS, Ø 0,6-0,8mm

4.2 | MINI NEPHROSCOPY

SHEATHS



Ordering no.	French	Description
AB-006-100	15/18 Fr.	Mini nephroscope sheath without perforation, continuous flow, stopcock fixed
AB-006-101	16.5/19.5 Fr.	
AB-006-102	21/24 Fr.	

DILATATORS



Ordering no.	French	Description
AB-006-110	15/18 Fr.	Dilatator with central cannulation $\varnothing$ 3 Fr.
AB-006-111	16.5/19.5 Fr.	Dilatator with 2 channels $\varnothing$ 2x3 Fr.
AB-006-112	21/24 Fr.	Dilatator with 2 channels $\varnothing$ 2x5 Fr.


APPLICATOR



Ordering no.	Description
AB-006-115	Applicator used for the mini nephroscope sheaths, incl. plunger

4.2 | MINI NEPHROSCOPY


INSTRUMENTS | FORCEPS FOR MINI NEPHROSCOPES




	Ordering no.	WL	Diameter	Description
	AB-006-242	400mm	Ø 1.8mm	Fenestrated stone grasping forceps, double action, rigid, with ring handle
	AB-006-342	500mm	Ø 1.8mm	
	AB-006-442	400mm	Ø 3.0mm	
	AB-006-243	400mm	Ø 1.8mm	Biopsy forceps, oval, single action, rigid, with ring handle
	AB-006-343	500mm	Ø 1.8mm	
	AB-006-443	400mm	Ø 3.0mm	
	AB-006-244	400mm	Ø 1.8mm	Grasping forceps alligator jaw, double action, rigid, with ring handle
	AB-006-344	500mm	Ø 1.8mm	
	AB-006-444	400mm	Ø 3.0mm	
	AB-006-245	400mm	Ø 1.8mm	Scissors, pointed, rigid, single action, with ring handle
	AB-006-345	500mm	Ø 1.8mm	
	AB-006-445	400mm	Ø 3.0mm	
	AB-006-246	400mm	Ø 1.8mm	3-Prongs retriever, rigid, with U-spring handle
	AB-006-346	500mm	Ø 1.8mm	
	AB-006-446	400mm	Ø 3.0mm	

4.3 | STANDARD NEPHROSCOPY

NEPHROSCOPES 200MM



Ordering no.	Working length	Viewing angle	Description
AA-006-001	200mm	0°	Nephroscope - 45° angled ocular - Ø 7.45mm - suction channel Ø 1.5mm - working channel Ø 5.0mm
AA-006-021		6°	
AA-006-004		15°	



Ordering no.	Working length	Viewing angle	Description
AA-006-002	200mm	0°	Nephroscope - 90° angled ocular - Ø 7.45mm - suction channel Ø 1.5mm - working channel Ø 5.0mm
AA-006-022		6°	
AA-006-005		15°	


	Ordering no.	Description
	AA-006-080	Flow valve, Ø 6.1mm, used for STEMA nephroscopes
	AA-006-081	Sealing cap, Ø 15mm/1.0mm, for AA-006-080
	AA-006-083	Sealing cap, Ø 20mm/2.0mm, for AA-006-001 - AA-006-034
	AJ-010-015	Stopcock left, large (passage 4mm), with LL-connection and tube attachment M6 x 0.35mm

4.3 | STANDARD NEPHROSCOPY

SHEATHS



Ordering no.	French	Description
AB-006-001	24 Fr.	Nephroscope sheath for continous flow, rotatable stopcock
AB-006-007	26 Fr.	Nephroscope sheath for continous flow, rotatable stopcock
AB-006-005	24 Fr.	Nephroscope sheath for continous flow, fixed stopcock
AB-006-006	26 Fr.	Nephroscope sheath for continous flow, fixed stopcock



Ordering no.	French	Description
AB-006-003	26 Fr.	Nephroscope sheath for continous flow, perforated, rotatable stopcock
AB-006-008	26 Fr.	Nephroscope sheath for continous flow, perforated, fixed stopcock


OBTURATORS




Ordering no.	French	Description
AB-006-002	24 Fr.	Obturator hollow used for nephroscopy sheaths (200mm)
AB-006-004	26 Fr.	Obturator hollow used for nephroscopy sheaths (200mm)

4.3 | STANDARD NEPHROSCOPY

NEPHROSCOPES 250MM



Ordering no.	Working length	Viewing angle	Description
AA-006-011	250mm	0°	Nephroscope - 45° angled ocular - Ø 7.45mm - suction channel Ø 1.5mm - working channel Ø 5.0mm
AA-006-031		6°	
AA-006-034		15°	



Ordering no.	Working length	Viewing angle	Description
AA-006-012	250mm	0°	Nephroscope - 90° angled ocular - Ø 7.45mm - suction channel Ø 1.5mm - working channel Ø 5.0mm
AA-006-032		6°	
AA-006-006		15°	



Ordering no.

Description

AA-006-080

Flow valve, Ø 6.1mm, used for STEMA nephroscopes

AA-006-081

Sealing cap, Ø 15mm/1.0mm, for AA-006-080

AA-006-083

Sealing cap, Ø 20mm/2.0mm, for AA-006-001 - AA-006-034

AJ-010-015


Stopcock left, large (passage 4mm), with LL-connection and tube attachment M6 x 0.35mm

4.3 | STANDARD NEPHROSCOPY

SHEATHS



Ordering no.	French	Description
AB-006-061	24 Fr.	Nephroscope sheath for continuous flow, rotatable stopcock
AB-006-067	26 Fr.	
AB-006-065	24 Fr.	Nephroscope sheath for continuous flow, fixed stopcock
AB-006-066	26 Fr.	



Ordering no.	French	Description
AB-006-063	26 Fr.	Nephroscope sheath for continuous flow, perforated, rotatable stopcock
AB-006-068	26 Fr.	Nephroscope sheath for continuous flow, perforated, fixed stopcock

OBTURATORS



Ordering no.	French	Description
AB-006-062	24 Fr.	Obturator hollow used for nephroscopy sheaths (250mm)
AB-006-064	26 Fr.	

4.4 | ACCESSORIES

TELESCOPE BOUGIES





Ordering no.	Description
AB-006-011	Set of telescope bougies consisting of: 6 dilating bougies 9, 12, 15, 18 ,21 and 24 Fr., incl. 1 rigid and 1 flexible guide probe
AB-006-012	Guide probe, rigid
AB-006-013	Guide probe, flexible

Ordering no.	French	Description
AB-006-020	9 Fr.	Dilatating sheath
AB-006-021	12 Fr.	
AB-006-022	15 Fr.	
AB-006-023	18 Fr.	
AB-006-024	21 Fr.	
AB-006-025	24 Fr.	
AB-006-026	27 Fr.	
AB-006-027	30 Fr.	

4.4 | ACCESSORIES

FORCEPS



	Ordering no.	WL	Diameter	Description
	AB-006-046	380mm	Ø 3.5mm	Stone grasping forceps, fenestrated, with U-spring handle
	AB-006-040	430mm		
	AB-006-048	380mm	Ø 3.5mm	Stone grasping forceps, fenestrated, with right handle
	AB-006-042	430mm		
	AB-006-047	380mm	Ø 3.5mm	3-Prongs retriever, larger stones, with U-spring handle
	AB-006-041	430mm		
	AB-006-049	380mm	Ø 3.5mm	Biopsy forceps, oval, with ring handle
	AB-006-043	430mm		
	AB-006-050	380mm	Ø 3.5mm	Grasping forceps alligator jaw, with ring handle
	AB-006-044	430mm		
	AB-006-051	380mm	Ø 3.5mm	Scissors, pointed, with ring handle
	AB-006-045	430mm		

4.4 | ACCESSORIES

CANNULAS

Ordering no.	Diameter	Description
AB-006-052		Puncture cannula, inner and outer cannula, 5 Pack
AB-006-053	Ø 3mm	Dilatation cannula, inner and outer cannula

5. URETERO-RENOSCOPY

5.1 | URETERO-RENOSCOPES SEMI-RIGID


HD URETERO-RENOSCOPES
URETERO-RENOSCOPES
BICHLER URETERO-RENOSCOPES
MARBERGER URETERO-RENOSCOPES
FORCEPS, SCISSORS, INSTRUMENTS
MONOPOLAR CONNECTING CABLES

5.2 | FLEXIBLE URETERO-RENOSCOPE

5. HD URETERO-RENOSCOPES

40.000 PX

17.000 PX



WORKING LENGTH



VIEWING ANGLE




FIELD OF VIEW



5.1 | URETERO-RENOSCOPES SEMI-RIGID

HD URETERO-RENOSCOPES



Ordering no.	Working length	Viewing angle	Description
AA-007-003	425mm	0°	HD Uretero-renaloscope semi-rigid - 45° angled ocular - Working Ø steps by 7.5 to 12 Fr. - Working channel 1 x 5 Fr. or 2 x 3 Fr. - 40.000 px
AA-007-016		5°	
AA-007-017		12°	

Ordering no.	Working length	Viewing angle	Description
AA-007-004	310mm	0°	HD Uretero-renaloscope semi-rigid - 45° angled ocular - Working Ø steps by 7.5 to 12 Fr. - Working channel 1 x 5 Fr. or 2 x 3 Fr. - 40.000 px
AA-007-018		5°	
AA-007-019		12°	

Ordering no.	Working length	Viewing angle	Description
AA-007-005	150mm	0°	HD Uretero-renaloscope semi-rigid - 45° angled ocular - Working Ø steps by 7.5 to 12 Fr. - Working channel 1 x 5 Fr. or 2 x 3 Fr. - 40.000 px
AA-007-036		5°	
AA-007-037		12°	

Delivery contains

HD uretero-renaloscope, brush Ø 3.0mm (length 590mm), brush Ø 3.8mm (length 350mm), sealing cap silicone (grey), Duckbill sealing

5.1 | URETERO-RENOSCOPES SEMI-RIGID

HD URETERO-RENOSCOPES



Ordering no.	Working length	Viewing angle	Description
AA-007-009	425mm	0°	HD Uretero-renaloscope semi-rigid - 45° angled ocular - Working Ø steps by 7.5 to 13.2 Fr. - Working channel 1 x 5 Fr. or 2 x 3 Fr. - 40.000 px
AA-007-028		5°	
AA-007-029		12°	

Delivery contains

HD uretero-renaloscope, brush Ø 3.0mm (length 590mm), brush Ø 3.8mm (length 350mm), sealing cap silicone (grey), Duckbill sealing

5.1 | URETERO-RENOSCOPES SEMI-RIGID

HD URETERO-RENOSCOPES



Ordering no.	Working length	Viewing angle	Description
AA-007-015	425mm	0°	HD Uretero-renaloscope semi-rigid - 45° angled ocular - Working Ø steps by 6 to 7.5 Fr. - Working channel 1 x 4 Fr. or 2 x 2.4 Fr. - 17.000 px
AA-007-034		5°	
AA-007-035		12°	

Ordering no.	Working length	Viewing angle	Description
AA-007-014	310mm	0°	HD Uretero-renaloscope semi-rigid - 45° angled ocular - Working Ø steps by 6 to 7.5 Fr. - Working channel 1 x 4 Fr. or 2 x 2.4 Fr. - 17.000 px
AA-007-032		5°	
AA-007-033		12°	

Ordering no.	Working length	Viewing angle	Description
AA-007-013	140mm	0°	HD Uretero-renaloscope semi-rigid - 45° angled ocular - Working Ø steps by 6 to 7.5 Fr. - Working channel 1 x 4 Fr. or 2 x 2.4 Fr. - 17.000 px
AA-007-030		5°	
AA-007-031		12°	

Delivery contains

HD uretero-renaloscope, brush Ø 3.0mm (length 590mm), brush Ø 3.8mm (length 350mm), sealing cap silicone (grey), Duckbill sealing

5.1 | URETERO-RENOSCOPY

HD URETERO-RENOSCOPES



Ordering no.	Working length	Viewing angle	Description
AA-007-006	425mm	0°	HD Uretero-renaloscope semi-rigid - 45° angled ocular - Working Ø steps by 4.9 Fr. to 6.5 Fr. - Working channel 1 x 3 Fr. - 17.000 px
AA-007-022		5°	
AA-007-023		12°	

Ordering no.	Working length	Viewing angle	Description
AA-007-007	310mm	0°	HD Uretero-renaloscope semi-rigid - 45° angled ocular - Working Ø steps by 4.9 Fr. to 6.5 Fr. - Working channel 1 x 3 Fr. - 17.000 px
AA-007-024		5°	
AA-007-025		12°	


Ordering no.	Working length	Viewing angle	Description
AA-007-008	150mm	0°	HD Uretero-renaloscope semi-rigid - 45° angled ocular - Working Ø steps by 4.9 Fr. to 6.5 Fr. - Working channel 1 x 3 Fr. - 17.000 px
AA-007-026		5°	
AA-007-027		12°	

Delivery contains

HD uretero-renaloscope, brush Ø 3.0mm (length 590mm), brush Ø 3.8mm (length 350mm), sealing cap silicone (grey), Duckbill sealing

5.1 | URETERO-RENOSCOPES SEMI-RIGID

URETERO-RENOSCOPES



Ordering no.	Working length	Viewing angle	Description
AA-007-001	310mm	8°	Uretero-renaloscope semi-rigid - 45° angled ocular - Working Ø steps by 7.5 Fr. to 10 Fr. - Working channel 1 x 5 Fr. or 2 x 3 Fr. - Field of view: 85°
AA-007-002	425mm		

Ordering no.	Working length	Viewing angle	Description
AA-007-011	310mm	8°	Uretero-renaloscope semi-rigid - 45° angled ocular - Working Ø steps by 6 Fr. to 7.5 Fr. - Working channel 1 x 4 Fr. or 2 x 2.4 Fr. - Field of view: 100°
AA-007-012	425mm		




Ordering no.	Description
AA-007-021	Instrument bridge used for uretero-renaloscopes (included in the delivery content)
AJ-010-001 AJ-010-003	Rubber cap for instrument channel Rubber cap for instrument channel and URS, Ø 0,6-0,8mm



5.1 | URETERO-RENSCOPES SEMI-RIGID

BICHLER URETERO-RENSCOPES



Ordering no.	Working length	Viewing angle	Description
AA-007-056	310mm	5°	BICHLER Uretero-renaloscope semi-rigid - 45° angled ocular - Working Ø steps by 4.5 Fr. to 6.5 Fr. - Working channel 1 x 3.3 Fr.
AA-007-055	430mm		


Ordering no.	Working length	Viewing angle	Description
AA-007-043	310mm	12°	BICHLER Uretero-renaloscope semi-rigid - 45° angled ocular - Working Ø steps by 6 Fr. to 7.5 Fr. - Working channel 1 x 4 Fr. or 2 x 2.4 Fr.
AA-007-044	430mm		

Ordering no.	Working length	Viewing angle	Description
AA-007-047	310mm	12°	BICHLER Uretero-renaloscope semi-rigid - 45° angled ocular - Working Ø steps by 8 Fr. to 9.8 Fr. - Working channel 1 x 5 Fr. or 2 x 3 Fr.
AA-007-048	430mm		

Ordering no.	Working length	Viewing angle	Description
AA-007-051	310mm	12°	BICHLER Uretero-renaloscope semi-rigid - 45° angled ocular - Working Ø steps by 8.5 Fr. to 11.5 Fr. - Working channel 1 x 6 Fr. or 2 x 4 Fr.
AA-007-052	430mm		

5.1 | URETERO-RENSCOPES SEMI-RIGID

MARBERGER URETERO-RENSCOPES



Ordering no.	Working length	Viewing angle	Description
AA-007-054	310mm	5°	MARBERGER Uretero-roscope semi-rigid - 45° angled ocular - Working Ø steps by 4.5 Fr. to 6.5 Fr. - Working channel 1 x 3.3 Fr.
AA-007-053	430mm		


Ordering no.	Working length	Viewing angle	Description
AA-007-041	310mm	12°	MARBERGER Uretero-roscope semi-rigid - 45° angled ocular - Working Ø steps by 6 Fr. to 7.5 Fr. - Working channel 1 x 4 Fr. or 2 x 2.4 Fr.
AA-007-042	430mm		



Ordering no.	Working length	Viewing angle	Description
AA-007-045	310mm	12°	MARBERGER Uretero-roscope semi-rigid - 45° angled ocular - Working Ø steps by 8 Fr. to 9.8 Fr. - Working channel 1 x 5 Fr. or 2 x 3 Fr.
AA-007-046	430mm		

Ordering no.	Working length	Viewing angle	Description
AA-007-049	310mm	12°	MARBERGER Uretero-roscope semi-rigid - 45° angled ocular - Working Ø steps by 8.5 Fr. to 11.5 Fr. - Working channel 1 x 6 Fr. or 2 x 4 Fr.
AA-007-050	430mm		

5.1 | URETERO-RENOSCOPY


FORCEPS FLEXIBLE | SEMI-RIGID | RIGID




	Description	Diameter	French	Flexible	Semi-rigid	Rigid
	Biopsy forceps, single action	Ø 2.3mm	7 Fr.	AC-010-014	AC-010-214	
	Grasping forceps, single action	Ø 2.3mm	7 Fr.		AC-011-214	
	Biopsy forceps oval, double action	Ø 1.0mm	3 Fr.	AC-001-011	AC-001-211	AC-001-011-R
		Ø 1.2mm	4 Fr.	AC-001-010	AC-001-210	AC-001-010-R
		Ø 1.6mm	5 Fr.	AC-001-012	AC-001-212	AC-001-012-R
		Ø 2.1mm	6 Fr.	AC-001-013	AC-001-213	AC-001-013-R
		Ø 2.3mm	7 Fr.	AC-001-014	AC-001-214	AC-001-014-R
		Ø 3.0mm	9 Fr.	AC-001-015	AC-001-215	AC-001-015-R
	Biopsy forceps oval, serrated, double action	Ø 1.0mm	3 Fr.	AC-002-011	AC-002-211	AC-002-011-R
		Ø 1.6mm	5 Fr.	AC-002-012	AC-002-212	AC-002-012-R
		Ø 2.1mm	6 Fr.	AC-002-013	AC-002-213	AC-002-013-R
		Ø 2.3mm	7 Fr.	AC-002-014	AC-002-214	AC-002-014-R
		Ø 3.0mm	9 Fr.	AC-002-015	AC-002-215	AC-002-015-R
	Biopsy forceps oval, serrated, teeth interlock, double action	Ø 1.6mm	5 Fr.	AC-015-012	AC-015-212	AC-015-012-R
		Ø 2.1mm	6 Fr.		AC-015-213	AC-015-013-R
		Ø 2.3mm	7 Fr.	AC-015-014	AC-015-214	AC-015-014-R
	Grasping forceps, double action	Ø 1.2mm	4 Fr.	AC-017-011	AC-017-211	AC-017-011-R
		Ø 2.3mm	7 Fr.	AC-017-014	AC-017-214	AC-017-014-R
	Grasping forceps, alligator jaw, double action	Ø 1.0mm	3 Fr.	AC-006-011	AC-006-211	AC-006-011-R
		Ø 1.2mm	4 Fr.	AC-006-013	AC-006-213	AC-006-013-R
		Ø 1.6mm	5 Fr.	AC-006-012	AC-006-212	AC-006-012-R
		Ø 2.3mm	7 Fr.	AC-006-014	AC-006-214	AC-006-014-R
		Ø 3.0mm	9 Fr.	AC-006-015	AC-006-215	AC-006-015-R

5.1 | URETERO-RENOSCOPY


FORCEPS FLEXIBLE | SEMI-RIGID | RIGID





	Description	Diameter	French	Flexible	Semi-rigid	Rigid
	Grasping forceps, DJ, double action	Ø 2.3mm	7 Fr.	AC-008-014	AC-008-214	AC-008-014-R
		Ø 3.0mm	9 Fr.	AC-008-015	AC-008-215	AC-008-015-R
	Grasping forceps teeth double action	Ø 1.2mm	4 Fr.	AC-016-011	AC-016-211	AC-016-011-R
		Ø 1.6mm	5 Fr.	AC-016-012	AC-016-212	AC-016-012-R
	Grasping forceps teeth interlock, double action	Ø 1.6mm	5 Fr.	AC-020-012	AC-020-212	AC-020-012-R
		Ø 2.3mm	7 Fr.	AC-020-014	AC-020-214	AC-020-014-R
	Grasping forceps rat tooth, double action	Ø 1.0mm	3 Fr.	AC-019-011	AC-019-211	AC-019-011-R
		Ø 1.6mm	5 Fr.	AC-019-012	AC-019-212	AC-019-012-R
	Grasping forceps for stone street, double action	Ø 1.2mm	4 Fr.	AC-007-011	AC-007-211	AC-007-011-R
		Ø 1.6mm	5 Fr.	AC-007-012	AC-007-212	AC-007-012-R
	Grasping forceps for big stones, double action	Ø 1.2mm	4 Fr.	AC-009-011	AC-009-211	AC-009-011-R
		Ø 1.6mm	5 Fr.	AC-009-012	AC-009-212	AC-009-012-R
		Ø 2.3mm	7 Fr.	AC-009-014	AC-009-214	AC-009-014-R
	3-Prongs retriever	Ø 1.0mm	3 Fr.	AC-012-011		
		Ø 1.6mm	5 Fr.	AC-012-012		
		Ø 2.3mm	7 Fr.	AC-012-014		


5.1 | URETERO-RENOSCOPY

SCISSORS FLEXIBLE | SEMI-RIGID | RIGID



	Description	Diameter	French	Flexible	Semi-rigid	Rigid
	Scissors, single action	Ø 1.0mm	3 Fr.	AC-003-011	AC-003-211	AC-003-011-R
		Ø 1.2mm	4 Fr.	AC-003-010	AC-003-210	AC-003-010-R
		Ø 1.6mm	5 Fr.	AC-003-012	AC-003-212	AC-003-012-R
		Ø 2.1mm	6 Fr.	AC-003-013	AC-003-213	AC-003-013-R
		Ø 2.3mm	7 Fr.	AC-003-014	AC-003-214	AC-003-014-R
		Ø 3.0mm	9 Fr.	AC-003-015	AC-003-215	AC-003-015-R
	Uretero slitting scissors, single action	Ø 1.6mm	4 Fr.		AC-018-211	AC-018-011-R
		Ø 2.3mm	5 Fr.		AC-018-212	AC-018-012-R
	Hook scissors, single action	Ø 2.3mm	7 Fr.	AC-004-014	AC-004-214	AC-004-014-R
		Ø 3.0mm	9 Fr.	AC-004-015	AC-004-215	AC-004-015-R
	Scissors, serrated flexible, double action					
		Ø 2.3mm	7 Fr.	AC-005-014	AC-005-214	AC-005-014-R

NEW VERSION OF FORCEPS WITH LUER LOCK AND ROTARY WHEEL




FEATURES

- 1 QUICK AND EASY LOCKING AND RELEASING OF THE HANDLE WITH CLICK LOCK
- 2 ROTARY WHEEL ENABLES PRECISE POSITIONING OF JAW SECTION
- 3 OVERLOAD PROTECTION AVAILABLE
- 4 SIDE-MOUNTED LUER LOCK

5.1 | URETERO-RENSCOPY

INSTRUMENTS SEMI-RIGID | DETACHABLE




Ordering no.	Description
AC-100-001	HANDLE FOR DETACHABLE INSTRUMENTS

	Description	Diameter	French	Complete	Insert
	Scissors, single action	Ø 1.0mm	3 Fr.	AC-003-210-T	AC-003-210-OT
		Ø 1.2mm	4 Fr.	AC-003-211-T	AC-003-211-OT
		Ø 1.6mm	5 Fr.	AC-003-212-T	AC-003-212-OT
		Ø 2.3mm	7 Fr.	AC-003-214-T	AC-003-214-OT
		Ø 3.0mm	9 Fr.	AC-003-215-T	AC-003-215-OT
	Scissors pointed, single action	Ø 1.6mm	5 Fr.	AC-004-622-T	AC-004-622-OT
		Ø 2.3mm	7 Fr.	AC-004-624-T	AC-004-624-OT
	Hook scissors, single action	Ø 2.3mm	7 Fr.	AC-004-214-T	AC-004-214-OT
		Ø 3.0mm	9 Fr.	AC-004-215-T	AC-004-215-OT

5.1 | URETERO-RENOSCOPY

INSTRUMENTS SEMI-RIGID | DETACHABLE




Ordering no.	Description
AC-100-001	HANDLE FOR DETACHABLE INSTRUMENTS

	Description	Diameter	French	Complete	Insert
	Biopsy forceps oval, double action	Ø 1.0mm	3 Fr.	AC-001-211-T	AC-001-211-OT
		Ø 1.2mm	4 Fr.	AC-001-210-T	AC-001-210-OT
		Ø 1.6mm	5 Fr.	AC-001-212-T	AC-001-212-OT
		Ø 2.3mm	7 Fr.	AC-001-214-T	AC-001-214-OT
		Ø 3.0mm	9 Fr.	AC-001-215-T	AC-001-215-OT
	Biopsy forceps oval, serrated, double action	Ø 1.2mm	4 Fr.	AC-016-211-T	AC-016-211-OT
		Ø 1.6mm	5 Fr.	AC-016-212-T	AC-016-212-OT
	Grasping forceps, alligator jaw, double action	Ø 1.0mm	3 Fr.	AC-006-211-T	AC-006-211-OT
		Ø 1.2mm	4 Fr.	AC-006-213-T	AC-006-213-OT
		Ø 1.6mm	5 Fr.	AC-006-212-T	AC-006-212-OT
		Ø 2.3mm	7 Fr.	AC-006-214-T	AC-006-214-OT
		Ø 3.0mm	9 Fr.	AC-006-215-T	AC-006-215-OT
	Grasping forceps for stone street, double action	Ø 1.2mm	4 Fr.	AC-007-211-T	AC-007-211-OT
		Ø 1.6mm	5 Fr.	AC-007-212-T	AC-007-212-OT
	Grasping forceps for big stones, double action	Ø 1.2mm	5 Fr.	AC-009-211-T	AC-009-211-OT
		Ø 1.6mm	5 Fr.	AC-009-212-T	AC-009-212-OT

5.1 | URETERO-RENOSCOPY

INSTRUMENTS RIGID | DETACHABLE




Ordering no.	Description
AC-100-001	HANDLE FOR DETACHABLE INSTRUMENTS

	Description	Diameter	French	Complete	Insert
	Scissors, single action	Ø 1.0mm	3 Fr.	AC-103-210-T	AC-103-210-OT
		Ø 1.2mm	4 Fr.	AC-103-211-T	AC-103-211-OT
		Ø 1.6mm	5 Fr.	AC-103-212-T	AC-103-212-OT
		Ø 2.3mm	7 Fr.	AC-103-214-T	AC-103-214-OT
		Ø 3.0mm	9 Fr.	AC-103-215-T	AC-103-215-OT
	Scissors pointed, single action	Ø 1.6mm	5 Fr.	AC-104-622-T	AC-104-622-OT
		Ø 2.3mm	7 Fr.	AC-104-624-T	AC-104-624-OT
	Hook scissors, single action	Ø 2.3mm	7 Fr.	AC-104-214-T	AC-104-214-OT
		Ø 3.0mm	9 Fr.	AC-104-215-T	AC-104-215-OT
	Biopsy punch forceps, single action	Ø 1.6mm	5 Fr.	AC-114-212-T	AC-114-212-OT

5.1 | URETERO-RENSCOPY

INSTRUMENTS RIGID | DETACHABLE



Ordering no.	Description
AC-100-001	HANDLE FOR DETACHABLE INSTRUMENTS

	Description	Diameter	French	Complete	Insert
	Biopsy forceps oval, double action	3 Fr. Ø 1.0mm	3 Fr.	AC-101-211-T	AC-101-211-OT
		Ø 1.2mm	4 Fr.	AC-101-210-T	AC-101-210-OT
		Ø 1.6mm	5 Fr.	AC-101-212-T	AC-101-212-OT
		Ø 2.3mm	7 Fr.	AC-101-214-T	AC-101-214-OT
		Ø 3.0mm	9 Fr.	AC-101-215-T	AC-101-215-OT
	Biopsy forceps oval, double action	Ø 1.2mm	4 Fr.	AC-116-211-T	AC-116-211-OT
		Ø 1.6mm	5 Fr.	AC-116-212-T	AC-116-212-OT
	Grasping forceps, alligator jaw, double action	Ø 1.0mm	3 Fr.	AC-106-211-T	AC-106-211-OT
		Ø 1.2mm	4 Fr.	AC-106-213-T	AC-106-213-OT
		Ø 1.6mm	5 Fr.	AC-106-212-T	AC-106-212-OT
		Ø 2.3mm	7 Fr.	AC-106-214-T	AC-106-214-OT
		Ø 3.0mm	9 Fr.	AC-106-215-T	AC-106-215-OT
	Grasping forceps for stone street, double action		4 Fr.	AC-107-211-T	AC-107-211-OT
			5 Fr.	AC-107-212-T	AC-107-212-OT
	Grasping forceps for big stones, double action	4 Fr. Ø 1.2mm	5 Fr.	AC-109-211-T	AC-109-211-OT
		Ø 1.6mm	5 Fr.	AC-109-212-T	AC-109-212-OT

5.1 | URETERO-RENSCOPY

STONE EXTRACTORS | 4 WIRES



Ordering no.	Description
AF-001-021	Handle for dormia stone extractor 3 and 5 Fr.
AF-001-024	Handle for dormia stone extractor 7 Fr.

Ordering no.	Wires	French	Working length	Description
AF-001-001	4	3 Fr.	60cm	Dormia stone extractor, basket with tubus only
AF-001-002	4	3 Fr.	60cm	Dormia stone extractor complete with handle
AF-001-015	4	3 Fr.	90cm	Dormia stone extractor, basket with tubus only
AF-001-016	4	3 Fr.	90cm	Dormia stone extractor complete with handle
AF-001-023	4	3 Fr.	120cm	Dormia stone extractor, basket with tubus only
AF-001-022	4	3 Fr.	120cm	Dormia stone extractor complete with handle
AF-001-003	4	5 Fr.	60cm	Dormia stone extractor, basket with tubus only
AF-001-004	4	5 Fr.	60cm	Dormia stone extractor complete with handle
AF-001-017	4	5 Fr.	90cm	Dormia stone extractor, basket with tubus only
AF-001-019	4	5 Fr.	90cm	Dormia stone extractor complete with handle
AF-001-005	4	7 Fr.	60cm	Dormia stone extractor, basket with tubus only
AF-001-006	4	7 Fr.	60cm	Dormia stone extractor complete with handle

5.1 | URETERO-RENSCOPY

STONE EXTRACTORS | 5-7 WIRES



Ordering no.	Description
AF-001-021	Handle for dormia stone extractor 3 and 5 Fr.
AF-001-024	Handle for dormia stone extractor 7 Fr.

Ordering no.	Wires	French	Working length	Description
AF-001-007	5	5 Fr.	60cm	Dormia stone extractor, basket with tubus only
AF-001-008	5	5 Fr.	60cm	Dormia stone extractor complete with handle
AF-001-025	5	5 Fr.	90cm	Dormia stone extractor, basket with tubus only
AF-001-018	5	5 Fr.	90cm	Dormia stone extractor complete with handle
AF-001-009	5	7 Fr.	60cm	Dormia stone extractor, basket with tubus only
AF-001-010	5	7 Fr.	60cm	Dormia stone extractor complete with handle
AF-001-011	6	5 Fr.	60cm	Dormia stone extractor, basket with tubus only
AF-001-012	6	5 Fr.	60cm	Dormia stone extractor complete with handle
AF-001-013	6	7 Fr.	60cm	Dormia stone extractor, basket with tubus only
AF-001-014	7	7 Fr.	60cm	Dormia stone extractor complete with handle

5.1 | URETERO-RENOSCOPY

MONOPOLAR FLEXIBLE ELECTRODES




	Description	French	Length	Ordering no.	Length	Ordering no.
	Loop electrode flexible, monopolar	5 Fr.	530mm	AD-003-003	650mm	AD-003-013
		7 Fr.	530mm	AD-003-005	650mm	AD-003-015
		9 Fr.	530mm	AD-003-007	650mm	AD-003-017
		10 Fr.	530mm	AD-003-008	650mm	AD-003-018
	Needle electrode flexible, monopolar	3 Fr.	530mm	AD-002-001	650mm	AD-002-011
		4 Fr.	530mm	AD-002-002	650mm	AD-002-012
		5 Fr.	530mm	AD-002-003	650mm	AD-002-013
		7 Fr.	530mm	AD-002-005	650mm	AD-002-015
	Button electrode flexible, monopolar	3 Fr.	530mm	AD-001-001	650mm	AD-001-011
		4 Fr.	530mm	AD-001-002	650mm	AD-001-012
		5 Fr.	530mm	AD-001-003	650mm	AD-001-013
		6 Fr.	530mm	AD-001-004	650mm	AD-001-014
		7 Fr.	530mm	AD-001-005	650mm	AD-001-015
		8 Fr.	530mm	AD-001-006	650mm	AD-001-016
		9 Fr.	530mm	AD-001-007	650mm	AD-001-017
		10 Fr.	530mm	AD-001-008	650mm	AD-001-018

5.1 | URETERO-RENOSCOPY


MONOPOLAR CONNECTING CABLES



CONNECTION INSTRUMENT SIDE



Unit Side	Ordering no.	Length	Diameter	Description
	AD-010-001	3m	Ø 4mm	HF Monopolar connecting cable for Berchtold, Martin, Aesculap
	AD-010-006	5m	Ø 4mm	
	AD-010-002	3m	Ø 5mm	HF Monopolar connecting cable for Erbe, ACC/ICC, Storz
	AD-010-007	5m	Ø 5mm	
	AD-010-003	3m	Ø 8mm	HF Monopolar connecting cable for STEMA HF-U Series, Valleylab/Conmed
	AD-010-008	5m	Ø 8mm	
	AD-010-024	4.5m	Ø 2mm	HF Monopolar connecting cable male for STEMA MASTER



Ordering no.	Description
 EK-015-060	Universal adapter used for monopolar HF cables (For cables: AD-010-001 to AD-010-008)

5.2 | FLEXIBLE URETERO-RENSCOPE FLEX II

- *High luminance*
- *Exact adjustment of focus*
- *Precise up-/down angulation*



Ordering no.	Product
AA-008-003	FLEXIBLE URETERO-RENSCOPE FLEX II

Technical data

Cystoscope

Insertion tube	Ø 2.8mm
Working Channel	Ø 1.2mm
Working length	680mm
Field of view	88°
Pixel	18.000
Angulation up-down	275°
Accessories	Suction valve, leak tester

Delivery contains

Flexible uretero-remoscope FLEX I, leak tester, suction valve

INCLUDED IN SCOPE OF DELIVERY:

AA-008-103

Leak tester

AA-008-105

Suction valve



AA-008-103



AA-008-105

5.2 | FLEXIBLE URETERO-RENSCOPE ACCESSORIES



Ordering no.	Description
AA-008-100	Bridge with Luer-Lock and 2 stop cocks for AA-008-002

Ordering no.	Description
AA-008-101	Single use cleaning brushes, 1.0mm, 120cm length, package of 250 pcs



Ordering no.	Description
AA-008-102	Suitcase grey for STEMA flexible scopes



Ordering no.	Description
AA-008-103	Leak tester



Ordering no.	Description
AA-008-104	Single use Luer-Lock adapter, 25 pcs./box



Ordering no.	Description
AA-008-105	Suction valve

6. LITHOTRIPSY

6.1 | LITHOTRIPSY
CYSTOSCOPES
SHEATHS
OBTURATORS
WORKING ELEMENTS

6.1 | LITHOTRIPSY

CYSTOSCOPES



Viewing angle	Endoscope	Diameter	Working length	Ordering no.
0°	Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-001
	HD Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-011

Viewing angle	Endoscope	Diameter	Working length	Ordering no.
30°	Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-003
	HD Cystoscope autoclavable	Ø 4.0mm	302mm	AA-002-013

SHEATHS



Ordering no.	Description
AB-002-010	Lithotripsy sheath with central valve used for AB-004-011




Ordering no.	Description
AB-002-013	Lithotripsy sheath used for AB-004-011

6.1 | LITHOTRIPSY

SHEATHS | OBTURATORS



Ordering no.	Description
AB-002-011	Diagnostic sheath with working channel for instruments 7 Fr.



Ordering no.	Description
AB-002-014	Adapter for use with Ellik



Ordering no.	Description
AB-002-012	Obturator for AB-002-010 + AB-002-013

6.1 | LITHOTRIPSY

WORKING ELEMENT



Ordering no.	Description
AB-004-011	Working element used for punch lithotripsy

LITHOTRIPTOR



Ordering no.	Description
AB-004-012	Mechanical lithotripter with handle and turning screw, Ø 24 Fr., 2 stopcocks with luer lock adapter, for use with scope AA-002-004

7. MORCELLATION

7.1 | MORCELLATION 26 FR.
NEPHROSCOPE
SHEATHS
OBTURATORS

7.2 | MORCELLATION 22 FR.
MORCESCOPE
SHEATHS
OBTURATORS


7. MORCELLATION

26 FR.

22 FR.


1 HD NEPHROSCOPE

2 SHEATHS
& OBTURATORS



3 MORCESCOPE

4 SHEATHS
& OBTURATORS






26 FR.

7.1 | MORCELLATION
NEPHROSCOPE



Ordering no.	Working length	Viewing angle	Description
AA-006-012	250mm	0°	Nephroscope with 90° angled ocular - Ø 7.45mm - Suction channel Ø 1.5mm - Working channel Ø 5.0mm

	Ordering no.	Description
	AB-004-010-QL	Adapter for use with nephroscope and resectoscope sheaths 26 Fr. STEMA Click-System
	AA-006-080	Flow valve, Ø 6.1mm, used for STEMA nephroscopes
	AA-006-081	Sealing cap, Ø 15mm/1.0mm, for AA-006-080
	AA-006-083	Sealing cap, Ø 20mm/2.0mm, for AA-006-001 - AA-006-034

7.1 | MORCELLATION

SHEATHS | OBTURATORS



Ordering no.	Description	French	Colour code
AH-001-302	Outer sheath with perforation, rotatable	26 Fr.	Yellow



Ordering no.	Description	French	Colour code
AH-001-301	Inner sheath with ceramic isolation, rotatable	26 Fr.	Yellow



Ordering no.	Description	French	For outer sheath	Colour code
AH-001-011	Standard Obturator for resectoscope sheaths	24 Fr.	26 Fr.	Yellow

22 FR.

7.2 | MORCELLATION
MORCESCOPE



Ordering no.	Working length	Direction of view	Description
AA-006-072	180mm	0°	Morcescope with 45° angled ocular - 22 Fr. Outer tube - Working channel 15 Fr. - Field of view 85° - 40.000 px

Ordering no.	Description
AA-006-084	Sealing cap white, Ø 13.5mm/10 pieces

7.2 | MORCELLATION

SHEATHS | OBTURATORS



Ordering no.	Description	French	Colour code
AH-002-302	Outer sheath with perforation, rotatable	22 Fr.	White



Ordering no.	Description	French	Colour code
AH-002-301	Inner sheath with ceramic isolation, rotatable	22 Fr.	White



Ordering no.	Description	French	For outer sheath	Colour code
AH-002-011	Standard Obturator for resectoscope sheaths	19 Fr.	22 Fr.	White



Ordering no.	Description	French	For outer sheath	Colour code
AH-002-013	Viewing Obturator	19 Fr.	22 Fr.	White

8. PEDIATRIC UROLOGY

8.1 | PEDIATRIC CYSTO-URETHROSCOPY

CYSTOSCOPES
SHEATHS
OBTURATORS
WORKING ELEMENTS
ENDOSCOPE BRIDGES
FORCEPS
MONOPOLAR ELECTRODES
CONNECTING CABLES

8.2 | PEDIATRIC RESECTOSCOPY

CYSTOSCOPES
SHEATHS
OBTURATORS
WORKING ELEMENTS
ENDOSCOPE BRIDGES
FORCEPS
MONOPOLAR ELECTRODES
CONNECTING CABLES

8.3 | PEDIATRIC URETHROTOMY

CYSTOSCOPES
SHEATHS
OBTURATORS
WORKING ELEMENTS
ENDOSCOPE BRIDGES
COLD KNIVES AND ACCESSORIES
CONNECTING CABLES

8.1 | PEDIATRIC CYSTO-URETHROSCOPY

CYSTOSCOPES | SHEATHS | OBTURATORS

CYSTOSCOPES



Ordering no.	Endoscope	Viewing angle	Diameter	Working length
AA-005-001	Pediatric Cystoscope autoclavable	0°	Ø 2.7mm	187.5mm
AA-005-003	Pediatric Cystoscope autoclavable	30°	Ø 2.7mm	187.5mm
AA-005-004	Pediatric Cystoscope autoclavable	70°	Ø 2.7mm	187.5mm

SHEATHS



Ordering no.	Description	French	Colour code
AB-001-001	Cysto-Urethroscope sheath with 2 fixed stopcocks , 1 x 3 Fr. instrument channel , with obturator	13 Fr.	Red
AB-001-001-13	Cysto-Urethroscope sheath with 2 fixed stopcocks , 1 x 5 Fr. or 2 x 4 Fr. instrument channel , with obturator	13 Fr.	Red

STANDARD OBTURATORS



Ordering no.	Description	For outer sheath	Colour code
AB-001-011	Standard obturator for Cysto-urethroscope sheaths	13 Fr.	Red

8.1 | PEDIATRIC CYSTO-URETHROSCOPY

CYSTOSCOPES | SHEATHS | OBTURATORS

WORKING ELEMENT



Ordering no.	Description	For outer sheath	Colour code
AB-003-005	Working element with Albaran deflector	13 Fr.	Red

ENDOSCOPE BRIDGE



Ordering no.	Description	For outer sheath	Colour code
AB-003-001	Endoscope bridge	13 Fr.	Red

8.1 | PEDIATRIC CYSTO-URETHROSCOPY

CYSTOSCOPES | SHEATHS | OBTURATORS

CYSTOSCOPES



Ordering no.	Endoscope	Viewing angle	Diameter	Working length
AA-004-011	Pediatric Cystoscope autoclavable	0°	Ø 2.0mm	187mm
AA-004-012	Pediatric Cystoscope autoclavable	30°	Ø 2.0mm	187mm

SHEATHS



Ordering no.	Description	French	Colour code
AB-001-001-9	Cysto-Urethroscope sheath with 2 fixed stopcocks , 1x3 FR. instrument channel , with obturator	9.5 Fr.	Blue
AB-001-001-11	Cysto-Urethroscope sheath with 2 fixed stopcocks , 1x3 FR. instrument channel , with obturator	11 Fr.	Green

STANDARD OBTURATORS



Ordering no.	Description	For outer sheath	Colour code
AB-001-011-9	Standard obturator for Cysto-urethroscope sheaths	9.5 Fr.	Blue
AB-001-011-11	Standard obturator for Cysto-urethroscope sheaths	11 Fr.	Green

8.1 | PEDIATRIC CYSTO-URETHROSCOPY


CYSTOSCOPES | SHEATHS | OBTURATORS

CYSTOSCOPES



Ordering no.	Endoscope	Viewing angle	Diameter	Working length
AA-004-013	Pediatric Cystoscope autoclavable	0°	Ø 1.9mm	145mm
AA-004-014	Pediatric Cystoscope autoclavable	30°	Ø 1.9mm	145mm

SHEATHS



Ordering no.	Description	French	Colour code
AB-001-001-7	Cysto-Urethroscope sheath with 2 fixed stopcocks	7.5 Fr.	Yellow


STANDARD OBTURATORS



Ordering no.	Description	For outer sheath	Colour code
AB-001-011-7	Standard obturator for Cysto-urethroscope sheaths	7.5 Fr.	Yellow

8.1 | PEDIATRIC CYSTO-URETHROSCOPY

FORCEPS | ELECTRODES



	Description	WL	Diameter	French	Flexible	Semi-rigid
	Biopsy forceps oval, double action	28cm	Ø1.0mm	3 Fr.	AC-001-021	AC-001-221
	Biopsy forceps oval, serrated, teeth interlock, double action	28cm	Ø1.0mm	3 Fr.	AC-015-021	AC-015-221
	Grasping forceps, alligator jaw, double action	28cm	Ø1.0mm	3 Fr.	AC-006-021	AC-006-221
	Grasping forceps, rat tooth, double action	28cm	Ø1.0mm	3 Fr.	AC-019-021	
	Scissors, Single action	28cm	Ø1.0mm	3 Fr.	AC-003-021	AC-003-221

ELECTRODES



	Ordering no.	French	Diameter	Description
	AD-002-021	3 Fr.	280mm	Needle electrode flexible, monopolar
	AD-001-021	3 Fr.	280mm	Button electrode flexible, monopolar


8.1 | PEDIATRIC CYSTO-URETHROSCOPY


CONNECTING CABLES | MONOPOLAR



CONNECTION INSTRUMENT SIDE



Unit Side	Ordering no.	Length	Diameter	Description
	AD-010-001	3m	Ø 4mm	HF Monopolar connecting cable for Berchtold, Martin, Aesculap
	AD-010-006	5m	Ø 4mm	
	AD-010-002	3m	Ø 5mm	HF Monopolar connecting cable for Erbe, ACC/ICC, Storz
	AD-010-007	5m	Ø 5mm	
	AD-010-003	3m	Ø 8mm	HF Monopolar connecting cable for STEMA HF-U Series, Valleylab/Conmed
	AD-010-008	5m	Ø 8mm	
	AD-010-024	4.5m	Ø 2mm	HF Monopolar connecting cable male for STEMA MASTER



	Ordering no.	Description
	EK-015-060	Universal adapter used for monopolar HF cables (For cables: AD-010-001 to AD-010-008)

8.2 | PEDIATRIC RESECTOSCOPY

RESECTOSCOPES | SHEATHS | OBTURATORS

CYSTOSCOPES



Ordering no.	Endoscope	Viewing angle	Diameter	Working length
AA-005-001	Pediatric Cystoscope autoclavable	0°	Ø 2.7mm	187.5mm

SHEATHS



Ordering no.	Description	French	Colour code
AH-003-001	Standard resectoscope sheath with 1 rotatable stopcock , without obturator	13 Fr.	Red

STANDARD OBTURATORS



Ordering no.	Description	For outer sheath	Colour code
AH-003-011	Standard Obturator for Cysto-urethroscope sheaths	13 Fr.	Red

8.2 | PEDIATRIC RESECTOSCOPY

WORKING ELEMENT | ENDOSCOPE BRIDGE



Ordering no.	Description	For outer sheath
AG-003-001	Working element used for pediatric resectoscopes, passive cutting , handle closed	13 Fr.




Ordering no.	Description	For outer sheath
AI-003-003	Endoscope bridge with 3 Fr. instrument channel	13 Fr.

8.2 | PEDIATRIC RESECTOSCOPY ELECTRODES



	Ordering no.	French	Color code	Description
	AE-004-002	13 Fr.	Red	Cutting loop electrode angled
	AE-004-018	13 Fr.	Red	Cutting loop electrode small angled
	AE-004-004	13 Fr.	Red	Rollerball electrode
	AE-004-016	13 Fr.	Red	Knife electrode angled
	AE-004-019	13 Fr.	Red	Hook electrode without ball



Ordering no.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives

8.2 | PEDIATRIC RESECTOSCOPY

RESECTOSCOPES | SHEATHS | OBTURATORS

CYSTOSCOPES



Ordering no.	Endoscope	Viewing angle	Diameter	Working length
AA-004-011	Pediatric Cystoscope autoclavable	0°	Ø 2.0mm	187.5mm

SHEATHS



Ordering no.	Description	French	Colour code
AH-005-001	Standard resectoscope sheath with 1 rotatable stopcock , without obturator	9 Fr.	Blue
AH-004-001	Standard resectoscope sheath with 1 rotatable stopcock , without obturator	11 Fr.	Green

STANDARD OBTURATORS



Ordering no.	Description	For outer sheath	Colour code
AH-005-011	Standard obturator for resectoscope sheaths	9 Fr.	Blue
AH-004-011	Standard obturator for resectoscope sheaths	11 Fr.	Green

8.2 | PEDIATRIC RESECTOSCOPY

WORKING ELEMENTS



Ordering no.	Description	For outer sheath
AG-005-001	Working element used for pediatric resectoscopes, passive cutting , handle closed	9 Fr.
AG-004-001	Working element used for pediatric resectoscopes, passive cutting , handle closed	11 Fr.

8.2 | PEDIATRIC RESECTOSCOPY ELECTRODES




	Ordering no.	French	Color code	Description
	AE-006-002	9 Fr.	Blue	Cutting loop electrode angled
	AE-005-002	11 Fr.	Green	
	AE-006-018	9 Fr.	Blue	Cutting loop electrode small angled
	AE-005-018	11 Fr.	Green	
	AE-006-004	9 Fr.	Blue	Rollerball electrode
	AE-005-004	11 Fr.	Green	
	AE-006-016	9 Fr.	Blue	Knife electrode angled
	AE-005-016	11 Fr.	Green	
	AE-006-019	9 Fr.	Blue	Hook electrode without ball
	AE-005-019	11 Fr.	Green	



Ordering no.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives


8.1 | PEDIATRIC RESECTOSCOPY

CONNECTING CABLES | MONOPOLAR



CONNECTION INSTRUMENT SIDE



Unit Side	Ordering no.	Length	Diameter	Description
	AD-010-001	3m	Ø 4mm	HF Monopolar connecting cable for Berchtold, Martin, Aesculap
	AD-010-006	5m	Ø 4mm	
	AD-010-002	3m	Ø 5mm	HF Monopolar connecting cable for Erbe, ACC/ICC, Storz
	AD-010-007	5m	Ø 5mm	
	AD-010-003	3m	Ø 8mm	HF Monopolar connecting cable for STEMA HF-U Series, Valleylab/Conmed
	AD-010-008	5m	Ø 8mm	
	AD-010-024	4.5m	Ø 2mm	HF Monopolar connecting cable male for STEMA MASTER



	Ordering no.	Description
	EK-015-060	Universal adapter used for monopolar HF cables (For cables: AD-010-001 to AD-010-008)

8.1 | PEDIATRIC RESECTOSCOPY

CONNECTING CABLES | MONOPOLAR



CONNECTION UNIT SIDE



Instrument Side	Ordering no.	Length	Diameter	Description
	AD-010-025	4.5m	90° angled connector	HF Monopolar connecting cable for STORZ resectoscope
	AD-010-026	4.5m	Ø 3mm	HF Monopolar connecting cable for OLYMPUS resectoscope



RECOMMENDED PRODUCTS:

CI-001-030 STEMA HF-U PRO TOUCH

- For all monopolar & bipolar techniques
- with 7" touch screen display
- provided to perform the most advanced electrosurgery
- ADC / APC System

CI-001-040 DOUBLE PEDAL FOOTSWITCH

8.3 | PEDIATRIC URETHROTOMY

CYSTOSCOPES | SHEATHS | OBTURATORS

CYSTOSCOPES



Ordering no.	Endoscope	Viewing angle	Diameter	Working length
AA-005-001	Pediatric Cystoscope autoclavable	0°	Ø2.7mm	187.5mm

SHEATHS



Ordering no.	Description	French	Colour code
AI-003-001	Standard urethrotome sheath with 2 fixed stopcocks , with obturator	13 Fr.	Red

STANDARD OBTURATORS



Ordering no.	Description	For outer sheath	Colour code
AI-003-004	Standard Obturator for urethrotome sheath	13 Fr.	Red

8.3 | PEDIATRIC URETHROTOMY

WORKING ELEMENT | ENDOSCOPE BRIDGE



Ordering no.	Description	For outer sheath
AG-003-001	Working element used for pediatric urethrosopes, passive cutting , handle open	13 Fr.




Ordering no.	Description	For outer sheath
AI-003-003	Endoscope bridge with 3 Fr. instrument channel	13 Fr.

8.3 | PEDIATRIC URETHROTOMY

COLD KNIVES



	Ordering no.	French	Color code	Description
	AI-003-007	13 Fr.	Red	Cold knife straight
	AI-003-009	13 Fr.	Red	Cold knife round
	AI-003-010	13 Fr.	Red	Cold knife ring shaped
	AI-003-011	13 Fr.	Red	Cold knife hook shaped



Ordering no.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives

8.3 | PEDIATRIC URETHROTOMY

CYSTOSCOPES | SHEATHS | OBTURATORS

CYSTOSCOPES



Ordering no.	Endoscope	Viewing angle	Diameter	Working length
AA-004-011	Pediatric Cystoscope autoclavable	0°	Ø2.0mm	187.5mm

SHEATHS



Ordering no.	Description	French	Colour code
AI-004-001	Standard urethrotome sheath with 2 fixed stopcocks , with obturator	11 Fr.	Green

STANDARD OBTURATORS



Ordering no.	Description	For outer sheath	Colour code
AI-004-004	Standard Obturator for urethrotome sheath	11 Fr.	Green

8.3 | PEDIATRIC URETHROTOMY

WORKING ELEMENT | ENDOSCOPE BRIDGE




Ordering no.	Description	For outer sheath
AG-004-001	Working element used for pediatric urethrosopes, passive cutting , handle open	11 Fr.

8.3 | PEDIATRIC URETHROTOMY

COLD KNIVES



	Ordering no.	French	Color code	Description
	AI-004-007	11 Fr.	Green	Cold knife straight
	AI-004-009	11 Fr.	Green	Cold knife round
	AI-004-010	11 Fr.	Green	Cold knife ring shaped
	AI-004-011	11 Fr.	Green	Cold knife hook shaped



Ordering no.	Description
AJ-010-006	Protection sleeve made of steel for sterilization and storage of electrodes and cold knives

9. ACCESSORIES

9.1 | INSTRUMENTS

DILATATING BOUGIES
METAL CATHETERS
ACCESSORIES CATHETERS

9.2 | STERILIZATION

9.3 | CLEANING AND MAINTENANCE

9.1 | INSTRUMENTS | DILATATING BOUGIES



VAN BUREN

DITTEL,
STRAIGHTDITTEL,
CURVED

French	Ordering no. VAN BUREN	Ordering no. DITTEL straight	Ordering no. DITTEL curved
8 Fr.	AJ-003-001	AJ-003-021	AJ-003-050
9 Fr.		AJ-003-022	AJ-003-051
10 Fr.	AJ-003-002	AJ-003-023	AJ-003-052
11 Fr.		AJ-003-024	AJ-003-053
12 Fr.	AJ-003-003	AJ-003-025	AJ-003-054
13 Fr.		AJ-003-026	AJ-003-055
14 Fr.	AJ-003-004	AJ-003-027	AJ-003-056
15 Fr.		AJ-003-028	AJ-003-057
16 Fr.	AJ-003-005	AJ-003-029	AJ-003-058
17 Fr.		AJ-003-030	AJ-003-059
18 Fr.	AJ-003-006	AJ-003-031	AJ-003-060
19 Fr.		AJ-003-032	AJ-003-061
20 Fr.	AJ-003-007	AJ-003-033	AJ-003-062
21 Fr.		AJ-003-034	AJ-003-063
22 Fr.	AJ-003-008	AJ-003-035	AJ-003-064
23 Fr.		AJ-003-036	AJ-003-065
24 Fr.	AJ-003-009	AJ-003-037	AJ-003-066
25 Fr.		AJ-003-038	AJ-003-067
26 Fr.	AJ-003-010	AJ-003-039	AJ-003-068
27 Fr.		AJ-003-040	AJ-003-069
28 Fr.	AJ-003-011	AJ-003-041	AJ-003-070
29 Fr.		AJ-003-042	AJ-003-071
30 Fr.	AJ-003-012	AJ-003-043	AJ-003-072
31 Fr.			
32 Fr.	AJ-003-013	AJ-003-044	AJ-003-074
33 Fr.			
34 Fr.	AJ-003-014	AJ-003-045	
35 Fr.			
36 Fr.	AJ-003-015	AJ-003-046	
37 Fr.			
38 Fr.	AJ-003-016		

9.1 | INSTRUMENTS | DILATATOR BOUGIES




GUYON

GUYON			
French	Ordering no. GUYON	French	Ordering no. GUYON
8 Fr.	AJ-003-170	25 Fr.	AJ-003-187
9 Fr.	AJ-003-171	26 Fr.	AJ-003-188
10 Fr.	AJ-003-172	27 Fr.	AJ-003-189
11 Fr.	AJ-003-173	28 Fr.	AJ-003-190
12 Fr.	AJ-003-174	29 Fr.	AJ-003-191
13 Fr.	AJ-003-175	30 Fr.	AJ-003-192
14 Fr.	AJ-003-176	31 Fr.	AJ-003-193
15 Fr.	AJ-003-177	32 Fr.	AJ-003-194
16 Fr.	AJ-003-178	33 Fr.	AJ-003-195
17 Fr.	AJ-003-179	34 Fr.	AJ-003-196
18 Fr.	AJ-003-180	35 Fr.	AJ-003-197
19 Fr.	AJ-003-181	36 Fr.	AJ-003-198
20 Fr.	AJ-003-182	37 Fr.	AJ-003-199
21 Fr.	AJ-003-183	38 Fr.	AJ-003-200
22 Fr.	AJ-003-184	39 Fr.	AJ-003-201
23 Fr.	AJ-003-185	40 Fr.	AJ-003-202
24 Fr.	AJ-003-186	50 Fr.	AJ-003-203

HEGAR

HEGAR	
French	Ordering no. HEGAR
3/4 Fr.	AJ-003-151
5/6 Fr.	AJ-003-152
7/8 Fr.	AJ-003-153
9/10 Fr.	AJ-003-154
11/12 Fr.	AJ-003-155
13/14 Fr.	AJ-003-156
15/16 Fr.	AJ-003-157
17/18 Fr.	AJ-003-158

9.1 | INSTRUMENTS | METAL CATHETERS



FOR MEN

FOR WOMEN

COXETER
PROSTATA

French	Ordering no.	Ordering no.	Ordering no.
6 Fr.		AJ-003-100	
7 Fr.	AJ-003-081		
8 Fr.	AJ-003-082	AJ-003-102	AJ-003-120
9 Fr.			
10 Fr.	AJ-003-084	AJ-003-104	AJ-003-121
11 Fr.			
12 Fr.	AJ-003-086	AJ-003-106	AJ-003-122
13 Fr.			
14 Fr.	AJ-003-088	AJ-003-108	AJ-003-123
15 Fr.			
16 Fr.	AJ-003-090	AJ-003-110	AJ-003-124
17 Fr.			
18 Fr.	AJ-003-092	AJ-003-112	AJ-003-125
19 Fr.			
20 Fr.	AJ-003-094	AJ-003-114	AJ-003-126
21 Fr.			
22 Fr.	AJ-003-096	AJ-003-116	AJ-003-127
23 Fr.			
24 Fr.	AJ-003-098	AJ-003-118	AJ-003-128

9.1 | INSTRUMENTS

CATHETER INTRODUCER



	Ordering no.	Length	Description
	AJ-003-142	39cm	Catheter introducer Guyon straight
	AJ-003-143	39cm	Catheter introducer Guyon curved

MEATUS CLAMPS

	Ordering no.	Length	Description
	AJ-003-140	8cm	Meatus clamp, Stockmann
	AJ-003-141	11cm	Meatus clamp, Strauss

SECRETION SPOON




	Ordering no.	Length	Description
	AJ-006-012	26.5cm	Secretion spoon Guttman


9.1 | INSTRUMENTS | TUBING CONNECTORS

	Ordering no.	Description
	AJ-010-011	Luer-Lock tubing connector, male Luer-Lock inside Ø 3.0mm, tubing connector Ø 6.0mm
	AJ-010-012	Luer-Lock tubing connector with stopcock
	AJ-010-013	Luer-Lock tubing connector, male Luer-Lock inside , tubing connector Ø 9.0mm
	AJ-010-018	Stopcock left with Luer-Lock, male, tubing connector Ø 9mm

9.2 | STERILIZATION | WIRE BASKET




Ordering no.	Dimension (l x w x h)	Description
FA-010-001	290x80x52mm	Wire basket for 1 endoscope with lid
FA-010-002	460x80x52mm	
FA-010-003	670x80x52mm	
FA-010-004	290x80x52mm	Wire basket for 2 endoscopes with lid
FA-010-005	460x80x52mm	



Ordering no.	Dimension (l x w x h)	Description
FA-010-006	640x150x77mm	Wire basket for 1 endoscope and 1 cold light cable with lid

9.2 | STERILIZATION | ENDOSCOPE CONTAINER



Ordering no.	Colour	Dimension (l x w x h)	Description
FA-001-001	Silver	450x70x70mm	Endoscope container, lid and bottom perforated (without filter)
FA-001-001-R	Red		
FA-001-001-G	Green		
FA-002-001	Silver	250x60x50mm	Endoscope container, lid and bottom perforated (without filter)
FA-002-001-R	Red		
FA-002-001-G	Green		



Ordering no.	Description
FA-001-000	Protection tube for endoscopes, plastic, 10 pcs, 19cm

9.3 | CLEANING AND MAINTENANCE

CLEANING PISTOL SET



Ordering no.	Description
EJ-030-010	Cleaning pistol complete inclusive 8 different tips

CLEANING BRUSH



Ordering no.	Diameter	Description
EJ-030-001	Ø 3mm	Cleaning brush, length 420mm , package of 10 pieces
EJ-030-002	Ø 5mm	Cleaning brush, length 330mm , package of 10 pieces
EJ-030-003	Ø 10mm	

GREASE & PASTE



Ordering no.	Weight	Description
EJ-030-031	3g	Special grease for stopcocks, valves and instruments
EJ-030-032	3g	Polishing paste for scopes, 50 pieces/ pack

ARE YOU INTERESTED
IN OUR PRODUCTS?
WE WILL GLADLY
SUBMIT YOU
AN OFFER.

(if available)

Company Name:

Do you have a question about our products?

10. INDEX

AA-002-001	8 , 37 , 98 , 145	AA-007-004	123	AB-001-001-7	158
AA-002-001-W	8 , 37	AA-007-005	123	AB-001-001-9	157
AA-002-002	8 , 37	AA-007-006	126	AB-001-001-11	157
AA-002-003	8 , 28 , 37 , 59 , 145	AA-007-007	126	AB-001-001-13	155
AA-002-003-W	8 , 37	AA-007-008	126	AB-001-002	9
AA-002-004	8 , 37	AA-007-009	124	AB-001-003	9
AA-002-004-W	8 , 37	AA-007-011	127	AB-001-004	9
AA-002-011	8 , 37 , 98 , 145	AA-007-012	127	AB-001-004	33
AA-002-011	106	AA-007-013	125	AB-001-005	9
AA-002-011-W	8 , 37	AA-007-014	125	AB-001-006	9
AA-002-012	8 , 37 , 59	AA-007-015	125	AB-001-006	33
AA-002-013	8 , 28 , 37 , 59 , 76 , 145	AA-007-016	123	AB-001-007	25
AA-002-013	33 , 90 , 91 , 93	AA-007-017	123	AB-001-010	10
AA-002-013-W	8 , 37	AA-007-018	123	AB-001-011	155
AA-002-014	8 , 37	AA-007-019	123	AB-001-011-7	158
AA-002-014-W	8 , 37	AA-007-021	127	AB-001-011-9	157
AA-002-021	8 , 37 , 98	AA-007-022	126	AB-001-011-11	157
AA-002-022	8 , 37	AA-007-023	126	AB-001-012	10
AA-002-023	8 , 28 , 37 , 59	AA-007-024	126	AB-001-013	10
AA-002-024	8 , 37	AA-007-025	126	AB-001-014	10
AA-003-001	103	AA-007-026	126	AB-001-015	10
AA-003-002	25 , 48 , 65 , 70	AA-007-027	126	AB-001-016	10
AA-003-004	65	AA-007-028	124	AB-001-022	10
AA-003-011	103	AA-007-029	124	AB-001-023	10
AA-003-012	25 , 48 , 65 , 70 , 79 , 83	AA-007-030	125	AB-001-024	10
AA-003-012-W	25	AA-007-031	125	AB-001-025	10
AA-003-014	65	AA-007-032	125	AB-001-026	10
AA-004-011	157 , 164 , 172	AA-007-033	125	AB-001-040	28
AA-004-012	157	AA-007-034	125	AB-001-041	28
AA-004-013	158	AA-007-035	125	AB-001-042	29
AA-004-014	158	AA-007-036	123	AB-001-043	29
AA-005-001	155 , 161 , 169	AA-007-037	123	AB-002-001	9
AA-005-003	155	AA-007-041	129	AB-002-002	9
AA-005-004	155	AA-007-042	129	AB-002-003	9
AA-006-001	114	AA-007-043	128	AB-002-004	9
AA-006-002	114	AA-007-044	128	AB-002-005	9
AA-006-004	114	AA-007-045	129	AB-002-010	145
AA-006-005	114	AA-007-046	129	AB-002-011	146
AA-006-006	116	AA-007-047	128	AB-002-012	146
AA-006-011	116	AA-007-048	128	AB-002-013	145
AA-006-012	116 , 150	AA-007-049	129	AB-002-014	146
AA-006-021	114	AA-007-050	129	AB-003-001	156
AA-006-022	114	AA-007-051	128	AB-003-002	11
AA-006-031	116	AA-007-052	128	AB-003-003	11 , 33
AA-006-032	116	AA-007-053	129	AB-003-004	11
AA-006-034	116	AA-007-054	129	AB-003-005	156
AA-006-050	111	AA-007-055	128	AB-003-006	12
AA-006-051	111	AA-007-056	128	AB-003-007	12 , 33
AA-006-052	111	AA-007-060	30	AB-003-008	12
AA-006-053	111	AA-008-002	31	AB-003-009	12
AA-006-069	109	AA-008-003	142	AB-003-010	26
AA-006-070	109	AA-008-100	32 , 143	AB-003-011	26
AA-006-072	152	AA-008-101	143	AB-003-012	26
AA-006-080	114 , 116 , 150	AA-008-102	32 , 143	AB-003-013	27
AA-006-081	109 , 114 , 116 , 150	AA-008-103	32 , 143	AB-003-014	27
AA-006-083	109 , 114 , 116 , 150	AA-008-103	31 , 142	AB-003-015	11
AA-006-084	152	AA-008-104	32 , 143	AB-003-020	27
AA-006-100	110	AA-008-105	32 , 143	AB-004-001	13 , 34 , 92 , 94
AA-007-001	127	AA-008-105	31 , 142	AB-004-002	13
AA-007-002	127	AB-001-000	25	AB-004-003	13 , 34 , 92 , 94
AA-007-003	123	AB-001-001	155	AB-004-004	13

10. INDEX

AB-004-005	13	AB-006-243	113	AC-001-206-T	19
AB-004-006	13	AB-006-244	113	AC-001-208	15
AB-004-007	14, 34	AB-006-245	113	AC-001-210	130
AB-004-008	13, 92, 94	AB-006-246	113	AC-001-210-OT	135
AB-004-009	14	AB-006-342	113	AC-001-210-T	135
AB-004-009-QL	14	AB-006-343	113	AC-001-211	130
AB-004-010	13	AB-006-344	113	AC-001-211-OT	135
AB-004-010-QL	150	AB-006-345	113	AC-001-211-T	135
AB-004-011	147	AB-006-346	113	AC-001-212	130
AB-004-012	147	AB-006-442	113	AC-001-212-OT	135
AB-006-001	115	AB-006-443	113	AC-001-212-T	135
AB-006-002	115	AB-006-444	113	AC-001-213	130
AB-006-003	115	AB-006-445	113	AC-001-214	130
AB-006-004	115	AB-006-446	113	AC-001-214-OT	135
AB-006-005	115	AB-007-001	87	AC-001-214-T	135
AB-006-006	115	AB-007-002	87	AC-001-215	130
AB-006-007	115	AB-007-003	87	AC-001-215-OT	135
AB-006-008	115	AB-007-004	88	AC-001-215-T	135
AB-006-011	118	AB-007-005	88	AC-001-221	159
AB-006-012	118	AB-007-006	88	AC-002-001	15
AB-006-013	118	AB-007-007	88	AC-002-001-R	15
AB-006-020	118	AB-007-008	87	AC-002-002	15
AB-006-021	118	AC-001-001	15	AC-002-002-R	15
AB-006-022	118	AC-001-001-R	15	AC-002-003	15
AB-006-023	118	AC-001-002	15, 33	AC-002-003-R	15
AB-006-024	118	AC-001-002-R	15	AC-002-004	15
AB-006-025	118	AC-001-003	15	AC-002-004-R	15
AB-006-026	118	AC-001-003-R	15	AC-002-005	15
AB-006-027	118	AC-001-004	15	AC-002-005-R	15
AB-006-040	119	AC-001-004-R	15	AC-002-011	130
AB-006-041	119	AC-001-005	15	AC-002-011-R	130
AB-006-042	119	AC-001-005-R	15	AC-002-012	130
AB-006-043	119	AC-001-008	15	AC-002-012-R	130
AB-006-044	119	AC-001-008-R	15	AC-002-013	130
AB-006-045	119	AC-001-010	130	AC-002-013-R	130
AB-006-046	119	AC-001-010-R	130	AC-002-014	130
AB-006-047	119	AC-001-011	130	AC-002-014-R	130
AB-006-048	119	AC-001-011-R	130	AC-002-015	130
AB-006-049	119	AC-001-012	130	AC-002-015-R	130
AB-006-050	119	AC-001-012-R	130	AC-002-201	15
AB-006-051	119	AC-001-013	130	AC-002-202	15
AB-006-052	120	AC-001-013-R	130	AC-002-203	15
AB-006-053	120	AC-001-014	130	AC-002-204	15
AB-006-061	117	AC-001-014-R	130	AC-002-205	15
AB-006-062	117	AC-001-015	130	AC-002-211	130
AB-006-063	117	AC-001-015-R	130	AC-002-212	130
AB-006-064	117	AC-001-021	159	AC-002-213	130
AB-006-065	117	AC-001-201	15	AC-002-214	130
AB-006-066	117	AC-001-201-OT	19	AC-002-215	130
AB-006-067	117	AC-001-201-T	19	AC-003-000	17
AB-006-068	117	AC-001-202	15	AC-003-000-R	17
AB-006-100	112	AC-001-202-OT	19	AC-003-001	17
AB-006-101	112	AC-001-202-T	19	AC-003-001-R	17
AB-006-102	112	AC-001-203	15	AC-003-002	17, 33
AB-006-110	112	AC-001-203-OT	19	AC-003-002-R	17
AB-006-111	112	AC-001-203-T	19	AC-003-003	17
AB-006-112	112	AC-001-204	15	AC-003-003-R	17
AB-006-115	112	AC-001-204-OT	19	AC-003-004	17
AB-006-170	109	AC-001-204-T	19	AC-003-004-R	17
AB-006-171	109	AC-001-205	15	AC-003-005	17
AB-006-242	113	AC-001-206-OT	19	AC-003-005-R	17

10. INDEX

AC-003-010	132	AC-004-204-OT	19	AC-006-212-OT	135
AC-003-010-R	132	AC-004-204-T	19	AC-006-212-T	135
AC-003-011	132	AC-004-205	17	AC-006-213	130
AC-003-011-R	132	AC-004-205-OT	19	AC-006-213-OT	135
AC-003-012	132	AC-004-205-T	19	AC-006-213-T	135
AC-003-012-R	132	AC-004-214	132	AC-006-214	130
AC-003-013	132	AC-004-214-OT	134	AC-006-214-OT	135
AC-003-013-R	132	AC-004-214-T	134	AC-006-214-T	135
AC-003-014	132	AC-004-215	132	AC-006-215	130
AC-003-014-R	132	AC-004-215-OT	134	AC-006-215-OT	135
AC-003-015	132	AC-004-215-T	134	AC-006-215-T	135
AC-003-015-R	132	AC-004-622-OT	134	AC-006-221	159
AC-003-021	159	AC-004-622-T	134	AC-007-001	16
AC-003-200	17	AC-004-624-OT	134	AC-007-001-R	16
AC-003-200-OT	19	AC-004-624-T	134	AC-007-002	16
AC-003-200-T	19	AC-005-004	17	AC-007-002-R	16
AC-003-201	17	AC-005-004-R	17	AC-007-011	131
AC-003-201-OT	19	AC-005-014	132	AC-007-011-R	131
AC-003-201-T	19	AC-005-014-R	132	AC-007-012	131
AC-003-202	17	AC-005-204	17	AC-007-012-R	131
AC-003-202-OT	19	AC-005-214	132	AC-007-201	16
AC-003-202-T	19	AC-006-001	15	AC-007-201-OT	19
AC-003-203	17	AC-006-001-R	15	AC-007-201-T	19
AC-003-204	17	AC-006-002	15	AC-007-202	16
AC-003-204-OT	19	AC-006-002-R	15	AC-007-202-OT	19
AC-003-204-T	19	AC-006-003	15	AC-007-202-T	19
AC-003-205	17	AC-006-003-R	15	AC-007-211	131
AC-003-205-OT	19	AC-006-004	15, 33	AC-007-211-OT	135
AC-003-205-T	19	AC-006-004-R	15	AC-007-211-T	135
AC-003-210	132	AC-006-005	15	AC-007-212	131
AC-003-210-OT	134	AC-006-005-R	15	AC-007-212-OT	135
AC-003-210-T	134	AC-006-011	130	AC-007-212-T	135
AC-003-211	132	AC-006-011-R	130	AC-008-004	16
AC-003-211-OT	134	AC-006-012	130	AC-008-004-R	16
AC-003-211-T	134	AC-006-012-R	130	AC-008-005	16
AC-003-212	132	AC-006-013	130	AC-008-005-R	16
AC-003-212-OT	134	AC-006-013-R	130	AC-008-014	131
AC-003-212-T	134	AC-006-014	130	AC-008-014-R	131
AC-003-213	132	AC-006-014-R	130	AC-008-015	131
AC-003-214	132	AC-006-015	130	AC-008-015-R	131
AC-003-214-OT	134	AC-006-015-R	130	AC-008-204	16
AC-003-214-T	134	AC-006-021	159	AC-008-205	16
AC-003-215	132	AC-006-201	15	AC-008-214	131
AC-003-215-OT	134	AC-006-201-OT	19	AC-008-215	131
AC-003-215-T	134	AC-006-201-T	19	AC-009-001	16
AC-003-221	159	AC-006-202	15	AC-009-001-R	16
AC-003-302	17	AC-006-202-OT	19	AC-009-002	16
AC-003-304	17	AC-006-202-T	19	AC-009-002-R	16
AC-003-612-OT	19	AC-006-203	15	AC-009-004	16
AC-003-612-T	19	AC-006-203-OT	19	AC-009-004-R	16
AC-003-614-OT	19	AC-006-203-T	19	AC-009-011	131
AC-003-614-T	19	AC-006-204	15	AC-009-011-R	131
AC-004-004	17	AC-006-204-OT	19	AC-009-012	131
AC-004-004-R	17	AC-006-204-T	19	AC-009-012-R	131
AC-004-005	17	AC-006-205	15	AC-009-014	131
AC-004-005-R	17	AC-006-205-OT	19	AC-009-014-R	131
AC-004-014	132	AC-006-205-T	19	AC-009-201	16
AC-004-014-R	132	AC-006-211	130	AC-009-201-OT	19
AC-004-015	132	AC-006-211-OT	135	AC-009-201-T	19
AC-004-015-R	132	AC-006-211-T	135	AC-009-202	16
AC-004-204	17	AC-006-212	130	AC-009-202-OT	19

10. INDEX

AC-009-204	16	AC-017-014	130	AC-103-201-OT	20
AC-009-211	131	AC-017-014-R	130	AC-103-201-T	20
AC-009-211-OT	135	AC-017-201	15	AC-103-202-OT	20
AC-009-211-T	135	AC-017-204	15	AC-103-202-T	20
AC-009-212	131	AC-017-211	130	AC-103-204-OT	20
AC-009-212-OT	135	AC-017-214	130	AC-103-204-T	20
AC-009-212-T	135	AC-018-001-R	17	AC-103-205-OT	20
AC-009-214	131	AC-018-002-R	17	AC-103-205-T	20
AC-010-004	15	AC-018-011-R	132	AC-103-210-OT	136
AC-010-014	130	AC-018-012-R	132	AC-103-210-T	136
AC-010-204	15	AC-018-201	17	AC-103-211-OT	136
AC-010-214	130	AC-018-202	17	AC-103-211-T	136
AC-011-004	15	AC-018-211	132	AC-103-212-OT	136
AC-011-204	15	AC-018-212	132	AC-103-212-T	136
AC-011-214	130	AC-019-001	16	AC-103-214-OT	136
AC-012-001	16	AC-019-001-R	16	AC-103-214-T	136
AC-012-002	16	AC-019-002	16	AC-103-215-OT	136
AC-012-004	16	AC-019-002-R	16	AC-103-215-T	136
AC-012-011	131	AC-019-011	131	AC-103-612-OT	20
AC-012-012	131	AC-019-011-R	131	AC-103-612-T	20
AC-012-014	131	AC-019-012	131	AC-103-614-OT	20
AC-015-002	15	AC-019-012-R	131	AC-103-614-T	20
AC-015-002-R	15	AC-019-021	159	AC-104-204-OT	20
AC-015-003-R	15	AC-019-201	16	AC-104-204-T	20
AC-015-004	15	AC-019-202	16	AC-104-205-OT	20
AC-015-004-R	15	AC-019-211	131	AC-104-205-T	20
AC-015-012	130	AC-019-212	131	AC-104-214-OT	136
AC-015-012-R	130	AC-020-002	16	AC-104-214-T	136
AC-015-013-R	130	AC-020-002-R	16	AC-104-215-OT	136
AC-015-014	130	AC-020-004	16	AC-104-215-T	136
AC-015-014-R	130	AC-020-004-R	16	AC-104-622-OT	136
AC-015-021	159	AC-020-012	131	AC-104-622-T	136
AC-015-202	15	AC-020-012-R	131	AC-104-624-OT	136
AC-015-203	15	AC-020-014	131	AC-104-624-T	136
AC-015-204	15	AC-020-014-R	131	AC-106-201-OT	20
AC-015-212	130	AC-020-202	16	AC-106-201-T	20
AC-015-213	130	AC-020-204	16	AC-106-202-OT	20
AC-015-214	130	AC-020-212	131	AC-106-202-T	20
AC-015-221	159	AC-020-214	131	AC-106-203-OT	20
AC-016-001	16	AC-101-201-OT	20	AC-106-203-T	20
AC-016-001-R	16	AC-101-201-T	20	AC-106-204-OT	20
AC-016-002	16	AC-101-202-OT	20	AC-106-204-T	20
AC-016-002-R	16	AC-101-202-T	20	AC-106-205-OT	20
AC-016-011	131	AC-101-203-OT	20	AC-106-205-T	20
AC-016-011-R	131	AC-101-203-T	20	AC-106-211-OT	137
AC-016-012	131	AC-101-204-OT	20	AC-106-211-T	137
AC-016-012-R	131	AC-101-204-T	20	AC-106-212-OT	137
AC-016-201	16	AC-101-206-OT	20	AC-106-212-T	137
AC-016-202	16	AC-101-206-T	20	AC-106-213-OT	137
AC-016-211	131	AC-101-210-OT	137	AC-106-213-T	137
AC-016-211-OT	135	AC-101-210-T	137	AC-106-214-OT	137
AC-016-211-T	135	AC-101-211-OT	137	AC-106-214-T	137
AC-016-212	131	AC-101-211-T	137	AC-106-215-OT	137
AC-016-212-OT	135	AC-101-212-OT	137	AC-106-215-T	137
AC-016-212-T	135	AC-101-212-T	137	AC-107-201-OT	20
AC-017-001	15	AC-101-214-OT	137	AC-107-201-T	20
AC-017-001-R	15	AC-101-214-T	137	AC-107-202-OT	20
AC-017-004	15	AC-101-215-OT	137	AC-107-202-T	20
AC-017-004-R	15	AC-101-215-T	137	AC-107-211-OT	137
AC-017-011	130	AC-103-200-OT	20	AC-107-211-T	137
AC-017-011-R	130	AC-103-200-T	20	AC-107-212-OT	137

10. INDEX

AC-107-212-T	137	AE-001-002	46	AE-003-040	69
AC-109-201-OT	20	AE-001-003	46	AE-003-041	74
AC-109-201-T	20	AE-001-004	46, 91	AE-003-043	74
AC-109-211-OT	137	AE-001-005	46	AE-003-044	74
AC-109-211-T	137	AE-001-006	47	AE-003-045	74
AC-109-212-OT	137	AE-001-007	47	AE-003-046	74
AC-109-212-T	137	AE-001-008	47	AE-004-002	163
AC-114-212-OT	136	AE-001-009	47	AE-004-004	163
AC-114-212-T	136	AE-001-010	46	AE-004-016	163
AC-116-211-OT	137	AE-001-011	46	AE-004-018	163
AC-116-211-T	137	AE-001-012	46	AE-004-019	163
AC-116-212-OT	137	AE-001-013	46	AE-005-002	166
AC-116-212-T	137	AE-001-014	47	AE-005-004	166
AD-001-001	23, 140	AE-001-015	47	AE-005-016	166
AD-001-002	23, 140	AE-001-016	47, 91	AE-005-018	166
AD-001-003	23, 34, 140	AE-001-017	47	AE-005-019	166
AD-001-004	23, 140	AE-001-030	63	AE-006-002	166
AD-001-005	23, 140	AE-001-031	63, 90, 93	AE-006-004	166
AD-001-006	23, 140	AE-001-032	63, 93	AE-006-016	166
AD-001-007	23, 140	AE-001-033	63	AE-006-018	166
AD-001-008	23, 140	AE-001-034	64, 90, 93	AE-006-019	166
AD-001-011	23, 140	AE-001-035	64	AF-001-001	21, 138
AD-001-012	23, 140	AE-001-036	63	AF-001-002	21, 138
AD-001-013	23, 140	AE-001-037	63	AF-001-003	21, 138
AD-001-014	23, 140	AE-001-038	63	AF-001-004	21, 138
AD-001-015	23, 140	AE-001-039	64	AF-001-005	21, 138
AD-001-016	23, 140	AE-001-040	63	AF-001-006	21, 138
AD-001-017	23, 140	AE-001-041	63	AF-001-007	22, 139
AD-001-018	23, 140	AE-001-042	63	AF-001-008	22, 139
AD-001-021	159	AE-001-043	63	AF-001-008	34
AD-002-001	23, 140	AE-001-044	64	AF-001-009	22, 139
AD-002-002	23, 140	AE-001-045	64	AF-001-010	22, 139
AD-002-003	23, 34, 140	AE-001-046	64	AF-001-011	22, 139
AD-002-005	23, 140	AE-001-047	64	AF-001-012	22, 139
AD-002-011	23, 140	AE-001-048	63	AF-001-013	22, 139
AD-002-012	23, 140	AE-002-001	46	AF-001-014	22, 139
AD-002-013	23, 140	AE-002-002	46	AF-001-015	21, 138
AD-002-015	23, 140	AE-002-003	46	AF-001-016	21, 138
AD-002-021	159	AE-002-004	46	AF-001-017	21, 138
AD-003-003	23, 140	AE-002-005	46	AF-001-018	22, 139
AD-003-005	23, 140	AE-002-006	47	AF-001-019	21, 138
AD-003-007	23, 140	AE-002-007	47	AF-001-021	21, 22, 138, 139
AD-003-008	23, 140	AE-002-010	46	AF-001-022	21, 138
AD-003-013	23, 140	AE-002-011	46	AF-001-023	21, 138
AD-003-015	23, 140	AE-002-012	46	AF-001-024	21, 22, 138, 139
AD-003-017	23, 140	AE-002-013	46	AF-001-025	22, 139
AD-003-018	23, 140	AE-002-014	47	AG-001-001	42, 101
AD-010-001	24, 57, 92, 141, 160, 167	AE-002-015	47	AG-001-001	91, 106
AD-010-002	24, 57, 141, 160, 167	AE-002-016	47	AG-001-001-HF	42, 101
AD-010-003	24, 34, 57, 141, 160, 167	AE-002-031	64	AG-001-002	42, 101
AD-010-006	24, 57, 141, 160, 167	AE-003-001	56	AG-001-003	43
AD-010-007	24, 57, 141, 160, 167	AE-003-003	56	AG-001-003-HF	43
AD-010-008	24, 57, 141, 160, 167	AE-003-004	56	AG-001-004	43
AD-010-011	75	AE-003-010	56	AG-001-006	78
AD-010-012	75	AE-003-014	56	AG-001-007	44
AD-010-013	75, 94	AE-003-016	56	AG-001-008	44
AD-010-015	75	AE-003-031	69	AG-001-009	78
AD-010-024	24, 57, 141, 160, 167	AE-003-032	69	AG-001-010	104
AD-010-025	58, 168	AE-003-033	69	AG-001-011	62
AD-010-026	58, 168	AE-003-034	69	AG-001-011	90, 93
AE-001-001	46, 91	AE-003-036	69	AG-001-012	62

10. INDEX

AG-001-013	78	AH-002-213	72, 84	AJ-003-002	176
AG-001-017	45	AH-002-216	72, 84	AJ-003-003	176
AG-001-018	45	AH-002-301	66, 80	AJ-003-004	176
AG-002-001	53	AH-002-301	50, 153	AJ-003-005	176
AG-002-002	53	AH-002-302	66, 80	AJ-003-006	176
AG-002-003	54	AH-002-302	50, 153	AJ-003-007	176
AG-002-004	54	AH-003-001	161	AJ-003-008	176
AG-002-006	55	AH-003-011	161	AJ-003-009	176
AG-002-007	82	AH-004-001	164	AJ-003-010	176
AG-002-008	82	AH-004-011	164	AJ-003-011	176
AG-002-011	68	AH-005-001	164	AJ-003-012	176
AG-002-012	68	AH-005-011	164	AJ-003-013	176
AG-002-013	71	AI-001-001	98, 106	AJ-003-014	176
AG-002-015	71	AI-001-002	99, 106	AJ-003-015	176
AG-002-016	86	AI-001-003	98, 106	AJ-003-016	176
AG-003-001	162, 170	AI-001-004	100, 106	AJ-003-021	176
AG-004-001	165, 173	AI-001-005	100	AJ-003-022	176
AG-005-001	165	AI-001-006	99	AJ-003-023	176
AG-100-006	78	AI-001-007	102	AJ-003-024	176
AG-100-007	78	AI-001-007	106	AJ-003-025	176
AH-001-001	40	AI-001-008	102	AJ-003-026	176
AH-001-002	40	AI-001-009	102	AJ-003-027	176
AH-001-011	41, 61, 77, 90, 91, 93, 151	AI-001-011	102	AJ-003-028	176
		AI-001-012	100	AJ-003-029	176
AH-001-012	41, 61	AI-002-001	97	AJ-003-030	176
AH-001-013	41, 61	AI-002-002	97	AJ-003-031	176
AH-001-015	41, 61	AI-003-001	169	AJ-003-032	176
AH-001-016	41, 61	AI-003-003	162, 170	AJ-003-033	176
AH-001-017	41, 61	AI-003-004	169	AJ-003-034	176
AH-001-018	41, 61	AI-003-007	171	AJ-003-035	176
AH-001-101	40	AI-003-009	171	AJ-003-036	176
AH-001-102	40	AI-003-010	171	AJ-003-037	176
AH-001-201	38	AI-003-011	171	AJ-003-038	176
AH-001-202	38	AI-004-001	172	AJ-003-039	176
AH-001-203	38	AI-004-004	172	AJ-003-040	176
AH-001-204	38	AI-004-007	174	AJ-003-041	176
AH-001-205	38	AI-004-009	174	AJ-003-042	176
AH-001-211	38	AI-004-010	174	AJ-003-043	176
AH-001-212	38	AI-004-011	174	AJ-003-044	176
AH-001-301	39, 60, 77, 90, 91, 93, 151	AI-005-001	103	AJ-003-045	176
		AI-005-003	103	AJ-003-046	176
AH-001-302	39, 60, 77, 151	AI-005-004	103	AJ-003-050	176
AH-001-302	90, 91, 93	AI-005-007	105	AJ-003-051	176
AH-001-303	39, 60	AI-005-008	105	AJ-003-052	176
AH-001-304	39, 60	AI-005-011	105	AJ-003-053	176
AH-001-305	39, 60	AJ-001-000	89	AJ-003-054	176
AH-001-306	39, 60	AJ-001-001	89	AJ-003-055	176
AH-001-308	39, 77	AJ-001-002	89	AJ-003-056	176
AH-001-311	39, 60	AJ-001-002	92, 94	AJ-003-057	176
AH-001-312	39, 60	AJ-001-003	89	AJ-003-058	176
AH-002-001	51	AJ-001-010	89	AJ-003-059	176
AH-002-011	67, 81	AJ-001-011	89	AJ-003-060	176
AH-002-011	52	AJ-001-012	89	AJ-003-061	176
AH-002-011	153	AJ-001-013	89	AJ-003-062	176
AH-002-013	67, 81	AJ-002-001	89	AJ-003-063	176
AH-002-013	52, 153	AJ-002-001	92, 94	AJ-003-064	176
AH-002-015	73, 85	AJ-002-002	89	AJ-003-065	176
AH-002-016	73, 85	AJ-002-003	89	AJ-003-066	176
AH-002-201	49	AJ-002-004	89	AJ-003-067	176
AH-002-211	72, 84	AJ-002-005	89	AJ-003-068	176
AH-002-212	49	AJ-003-001	176	AJ-003-069	176

10. INDEX

AJ-003-070	176	AJ-003-128	178	AJ-010-003	26 , 111 , 127
AJ-003-071	176	AJ-003-140	179	AJ-010-006	47 , 56 , 64 , 69 , 74 , 90 , 92 , 94 , 102 , 105 , 106 , 163 , 166 , 171 , 174
AJ-003-072	176	AJ-003-141	179		
AJ-003-074	176	AJ-003-142	179	AJ-010-011	180
AJ-003-081	178	AJ-003-143	178 , 179	AJ-010-012	180
AJ-003-082	178	AJ-003-151	177	AJ-010-013	180
AJ-003-084	178	AJ-003-152	177	AJ-010-015	114 , 116
AJ-003-086	178	AJ-003-153	177	AJ-010-018	180
AJ-003-088	178	AJ-003-154	177	AJ-010-020	90
AJ-003-090	178	AJ-003-155	177	CI-001-030	58 , 168
AJ-003-092	178	AJ-003-156	177	CI-001-040	58 , 168
AJ-003-094	178	AJ-003-157	177	EJ-030-001	183
AJ-003-096	178	AJ-003-158	177	EJ-030-002	183
AJ-003-098	178	AJ-003-187	177	EJ-030-003	183
AJ-003-100	178	AJ-003-188	177	EJ-030-010	183
AJ-003-102	178	AJ-003-189	177	EJ-030-031	183
AJ-003-104	178	AJ-003-190	177	EJ-030-032	183
AJ-003-106	178	AJ-003-191	177	EK-015-060	24 , 57 , 141 , 160 , 167
AJ-003-108	178	AJ-003-192	177	FA-001-000	182
AJ-003-110	178	AJ-003-193	177	FA-002-001	182
AJ-003-112	178	AJ-003-194	177	FA-002-001-G	182
AJ-003-114	178	AJ-003-195	177	FA-002-001-R	182
AJ-003-116	178	AJ-003-196	177	FA-010-001	181
AJ-003-118	178	AJ-003-197	177	FA-010-002	181
AJ-003-120	178	AJ-003-198	177	FA-010-003	181
AJ-003-121	178	AJ-003-199	177	FA-010-004	181
AJ-003-122	178	AJ-003-200	177	FA-010-005	181
AJ-003-123	178	AJ-003-201	177	FA-010-006	181
AJ-003-124	178	AJ-003-202	177		
AJ-003-125	178	AJ-003-203	177		
AJ-003-126	178	AJ-006-012	179		
AJ-003-127	178	AJ-010-001	11 , 26 , 111 , 127		

WE COVER THE FOLLOWING AREAS IN SURGERY WITH HIGH QUALITY PRODUCTS

ARTHROSCOPY

BRONCHOSCOPY

GASTROENTEROLOGY

GYNAECOLOGY

ENDOSCOPY

ENT

LAMINECTOMY

LAPAROSCOPY

MICROMOTOR SYSTEM

OMS

PROCTOSCOPY

THORACOSCOPY

UROLOGY

SPINE

CONTAINER SYSTEMS

STERILIZATION



STEMA MEDIZINTECHNIK GMBH | TAKE-OFF GEWERBEPARK 74 | 78579 NEUHAUSEN OB ECK
INFO@STEMA-MEDIZINTECHNIK.DE | WWW.STEMA-MEDIZINTECHNIK.DE | +49 7461 910 55 09

10/24 | PRINTED IN GERMANY | COPYRIGHT BY STEMA MEDIZINTECHNIK GMBH |
WE RESERVE THE RIGHT TO MAKE ALTERATIONS